**Stryker ESCS BV Dubai Branch**

Baniyas Road, Twin Towers 1101-1102

P.O.Box 41446, Deira – Dubai

United Arab Emirates

P +971 4 205 4222

F +971 4 205 4299

www.stryker.com

**27/03/2016**

**URGENT: Field Safety Notice**

**FSCA identifier:** Product Field Action RA2016-036  
**Type of Action:** Field Safety Corrective Action: Advisory Notice  
**Description:** GDC 360  
**Catalog #:** M00334811300 & M0033471020SR0  
**Lot Code:** 17074659, 17723781, 17735100

Dear customer:

Stryker Neurovascular has initiated this advisory notice for the devices identified above. The intent of this letter is to list all known potential hazards associated with the below noted issue and list the risk mitigation factors associated with the use of the product.

**Issue:**

Stryker Neurovascular has become aware that the incorrect DFU was included with the above listed product. The general GDC DFU was added to the product instead of the GDC 360 specific DFU. Both DFU's are equivalent

**Potential Risk:**

No adverse health consequences are anticipated.  
There is no risk to patients associated with previous use of affected devices.

**Risk Mitigation:**

Follow the instructions below under Immediate Actions Required.

Our records indicate that you have received at least one of the subject devices and you are therefore affected by this action.

We request that you read this notice carefully and complete the following actions:

1. Immediately check your internal inventory.
2. The affected items may continue to be used.
3. Please follow the instructions listed in the Directions For Use (DFU) accompanying this notice.
4. Circulate this Field Safety Notice internally to all interested/affected parties.
5. Maintain awareness of this notice internally until all required actions have been completed within your facility
6. Inform Stryker if any of the subject devices have been distributed to other organizations.
  - a) *Please provide contact details so that Stryker can inform the recipients appropriately*
7. Please inform Stryker of any adverse events concerning the use of the subject devices.

Please comply with any local regulations concerning the notification of adverse events to your National or local Competent Authorities.
8. Complete the attached customer response form. It may be that you no longer have any physical inventory on site. Completing this form will allow us to update our records and will also negate the need for us to send any further unnecessary communications on this matter. Therefore please complete even if you no longer have any of the subject devices in your physical inventory.
9. Return the completed form to your nominated Stryker Representative (indicated below) for this PFA.

*We request that you respond to this notice within 7 calendar days from the date of receipt. The target date for completion of this action is 01 May 2016 and your timely response will enable us to ensure that we meet this target.*

*Your designated contact person for this action is given below. Should you have any queries concerning this matter please do not hesitate to contact them directly*

*Name: Azza Ibrahim  
Position: RAQA Manager MEA  
Email: [azza.ibrahim@stryker.com](mailto:azza.ibrahim@stryker.com)  
Telephone: 00971 4205 4204*

On behalf of Stryker we thank you sincerely for your help and support in completing this action within the target date and regret any inconvenience that may be cause. We would like to reassure you that Stryker is committed to ensuring that only conforming devices, meeting our high internal quality standards, remaining on the market.

Yours Sincerely,

**Azza Ibrahim**  
RA/QA Manager MEA

**Stryker ESCS BV - Dubai Branch**  
Baniyas Road, Twin Towers, 1101&1102  
P.O. Box 41446  
Dubai - UAE

Tel: +971 4 205 4222 (Ext: 4204)  
Direct: +971 4 205 4204  
Fax: +971 4 205 4299  
Cell: +971 50 100 7354  
E-mail: [azza.ibrahim@stryker.com](mailto:azza.ibrahim@stryker.com)

**Stryker ESCS BV Dubai Branch**  
Baniyas Road, Twin Towers 1101-1102  
P.O.Box 41446, Deira – Dubai  
United Arab Emirates  
P +971 4 205 4222  
F +971 4 205 4299  
www.stryker.com

**STRYKER® NEUROVASCULAR**  
**FIELD SAFETY CORRECTIVE ACTION ACKNOWLEDGMENT FORM**

**FSCA identified:** Product Field Action RA2016-036  
**Type of Action:** Field Safety Corrective Action: Advisory Notice  
**Description:** GDC 360  
**Catalog #:** M00334811300 & M0033471020SR0  
**Lot Code:** 17074659, 17723781, 17735100

I have received the notification from Stryker stating that they have initiated a product field action for the above referenced product.

I acknowledge receipt of the of this **URGENT FIELD SAFETY NOTICE: ADVISORY NOTICE**

**Comments:** \_\_\_\_\_

<b>Form completed by:</b>			
<b>Contact Name</b>		<b>Facility</b>	
<b>Contact address</b>		<b>Signature</b>	
		<b>Phone</b>	
<b>Date</b>		<b>Email</b>	

---

Please send this signed and dated form to: [azza.ibrahim@stryker.com](mailto:azza.ibrahim@stryker.com)

# GDC<sup>®</sup> 360<sup>o</sup>

## Detachable Coil

<b>Directions for Use</b>	<b>2</b>
<b>Instrucciones de uso</b>	<b>6</b>
<b>Mode d'emploi</b>	<b>10</b>
<b>Gebrauchsanweisung</b>	<b>14</b>
<b>Istruzioni per l'uso</b>	<b>18</b>
<b>Gebruiksaanwijzing</b>	<b>22</b>
<b>Instruções de Utilização</b>	<b>26</b>

# GDC® 360°

## Detachable Coil

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Stryker Neurovascular representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

#### DEVICE DESCRIPTION

Stryker Neurovascular's GDC 360° Detachable Coils are comprised of three coil types: **GDC-10 360° Coil SR**, **GDC-10 360° Soft Coil SR** and **GDC-18 360° Coil**. GDC-10 360° Coil SR and GDC-10 360° Soft Coil SR are stretch resistant coils (refer to the Stretch Resistant section below for a description of this feature). GDC-18 360° Coils are non-stretch resistant coils.

Each of the coil types is manufactured from a platinum/tungsten alloy wire. The distal end of the main coil is formed such that there is a smaller distal loop at the end of the main coil to facilitate placement of the coil. The diameter of the distal loop is 75% that of the rest of the main coil loops.

Stretch Resistant GDC Detachable Coils: In addition to the standard GDC 360° device (GDC-18 360° Coil), two stretch resistant devices (GDC-10 360° Coil and GDC-10 360° Soft Coil) are available which incorporate a length of multi-strand material through the center of the coil to help resist stretching.

#### INDICATIONS FOR USE/INTENDED USE

GDC 360° Detachable Coils are intended for embolization of those intracranial aneurysms that – because of their morphology, their location, or the patient's general medical condition – are considered by the treating neurosurgical team to be a) very high risk for management by traditional operative techniques, or b) inoperable. GDC 360° Detachable Coils are also intended for embolization of other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae.

GDC 360° Detachable Coils are also intended for arterial and venous embolizations in the peripheral vasculature.

This device should only be used by physicians who have received appropriate training in interventional neuroradiology or interventional radiology and preclinical training on the use of this device as established by Stryker Neurovascular.

#### CONTRAINDICATIONS

None known.

#### ADVERSE EVENTS

Potential complications include, but are not limited to: post-embolization syndrome, hematoma, hemorrhage, vessel perforation, emboli (foreign, thromboembolic), ischemia, vasospasm, revascularization, inadequate occlusion, and neurological deficits including stroke and possibly death.

#### HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

The GDC 360° Detachable Coil with delivery wire and introducer sheath is packaged inside a protective, plastic dispenser coil and held in place with a plastic wire retainer. The assembly is packaged inside a protective pouch which will remain sterile unless opened or damaged.

**Warning:** Do not use if sterile packaging has been breached or damaged. Do not use if past "use by" date on package.

#### Handling & Storage

Store in a cool, dry, dark, place.

#### REQUIRED ADDITIONAL ITEMS

**Note:** When using a Power Supply (M00345100740 or M00345100840), a second, backup Power Supply is required for all procedures.

When using the InZone® Detachment System (M00345100940 or M00345100950), there must be enough detachment systems in inventory to finish the case, plus one additional system for backup.

In addition to the GDC 360° Detachable Coils and Stryker Neurovascular power sources, the following items are required:

- If using a Power Supply (M00345100740 or M00345100840):
  - One set of Detachable Coil Connecting Cables, part number M00345110240
  - Two fresh 9-Volt batteries for the Power Supply for each patient procedure
- If using the InZone Detachment System (M00345100940 or M00345100950):
  - One IZDS® Connecting Cable, part number M00345110250
- 6F (2 mm) non-tapered guiding catheter to facilitate Stryker Neurovascular 2-tip infusion catheter access to the vessel
- Continuous flush setup with two rotating hemostatic valves (RHVs), three bags of appropriate flush, one 3-way stopcock, and one 1-way stopcock
- Sterile 20 or 22 gauge (0.9 or 0.7 mm) stainless steel hypodermic needle with no Teflon coating (to provide ground for the current delivery system)
- Alcohol-dampened gauze to clean the proximal end of the delivery wire prior to attaching the InZone Detachment System or, if detaching with the Power Supply, prior to attaching the red cable of the Detachable Coil Connecting Cables (M00345110240)
- Stryker Neurovascular "–10" or "–18" infusion catheter with 2 tip markers

For GDC-10 360° Detachable Coils and GDC-10 360° Soft Detachable Coils compatible catheters are:

- FasTracker®-10, 2 tip marker Microcatheter
- Excelsior® 1018®, 2 tip marker Microcatheter
- Excelsior SL-10®, 2 tip marker Microcatheter
- Tracker® Excel™-14, 2 tip marker Microcatheter

For GDC-18 360° Coils compatible catheters are:

- Excelsior 1018, 2 tip marker Microcatheter
- FasTracker-18, 2 tip marker Microcatheter
- Stryker Neurovascular 0.25 mm (0.010 in), 0.36 mm (0.014 in) or 0.41 mm (0.016 in) steerable guidewires

#### PREPARATIONS FOR USE

##### Coil Size & Catheter Selection

Correct coil selection increases GDC 360° Detachable Coil effectiveness and patient safety. Occlusive efficiency is, in part, a function of compactness and overall mass. In order to choose the optimum coil for any given lesion, examine preembolization angiograms. The appropriate coil size should be chosen based upon angiographic assessment of the diameter of the vessel, aneurysm dome and/or ostium.

**Caution:** Exercise caution during placement of small coils at the base of any aneurysm. When accessing aneurysms, the diameter of the first or second coils should never be less than the width of the ostium, given the potential for these coils to migrate.

**Caution:** GDC-18 360° Detachable Coils should only be delivered through Stryker Neurovascular "–18" 2-tip infusion catheter. GDC-10 360° Detachable Coils and GDC-10 360° Soft Detachable Coils should only be delivered through Stryker Neurovascular "–10" infusion catheters. If advanced through a Stryker Neurovascular "–18" catheter, the GDC-10 Detachable Coil could fold back upon itself resulting in jamming, stretching or breaking.

**Caution:** The GDC-18 360° Detachable Coil is designed to be delivered through a Stryker Neurovascular 2-tip infusion catheter. The compatibility of the GDC Detachable Coil system with other catheters and with other coil delivery devices has not been established.

#### Continuous Flush Setup

**In order to achieve optimal performance of the GDC 360° Detachable Coils, and to reduce the risk of thromboembolic complications, it is critical that a continuous infusion of appropriate flush solution be maintained** between a) the femoral sheath and guiding catheter, b) the 2-tip infusion and guiding catheters, and c) the 2-tip infusion catheter and Stryker Neurovascular guidewires and GDC 360° delivery wire. Continuous flush also reduces the potential for thrombus formation on, and crystallization of infusate around, the coil detachment zone.

1. Attach an RHV to the hub of the guiding catheter. Attach a 3-way stopcock to the side arm of the RHV, then connect a line for continuous infusion of appropriate solution.

- Attach a second RHV to the hub of the 2-tip infusion catheter. Attach a 1-way stopcock to the side arm of the RHV, then connect a line for continuous flushing of appropriate solution. One drop from the pressure bag every 3-5 seconds is suggested.
- Check that all fittings are secure so that air is not introduced into guiding or 2-tip infusion catheters during continuous flush.


Figure 1. Example of Continuous Flush Setup

#### INZONE® DETACHMENT SYSTEM PREPARATION FOR USE (FOR M00345100940 OR M00345100950)

If using the InZone Detachment System to detach GDC 360° Coils, refer to the InZone Detachment System Directions for Use packaged with each device for “Required Additional Items” and “Removal from Packaging” instructions.

#### DIRECTIONS FOR USE

Carefully catheterize the lesion to be treated. The access system should include a guiding catheter of sufficient inner diameter (ID) to accept a 2-tip infusion catheter and to permit adequate contrast infusion around the 2-tip infusion catheter for fluoroscopic road mapping. Measure the size of the aneurysm to be treated and select an appropriately sized GDC 360° Detachable Coil.

- Slowly and simultaneously remove the GDC 360° Coil and introducer sheath from the dispenser coil. Inspect proximal section of delivery wire for irregularities. If irregularities exist, replace with a new GDC 360° Detachable Coil.

**Warning:** Damaged delivery wires may cause vessel injury or unpredictable distal tip response during coil deployment.

- To unlock the Twist-Lock mechanism on the introducer sheath, gently pinch the introducer sheath on both sides of the Twist-Lock mechanism and rotate proximal side counter-clockwise (Figure 2). This releases the delivery wire/coil assembly to move freely within the introducer sheath.


Figure 2. Unlock via Counter-Clockwise Rotation


Figure 3. Lock via Clockwise Rotation

- Slowly advance the coil out of the introducer sheath and into the palm of a gloved hand and inspect for:
  - Irregularities at the coil detachment zone, e.g., delivery wire kinking
  - Loss of circular memory (shape)
 If any of these irregularities exist, replace with a new coil.

**Warning:** Utilization of damaged coils may affect coil delivery to, and stability inside, the vessel or aneurysm, possibly resulting in coil migration and/or stretching.

- Gently immerse the coil, including the coil junction, in heparinized saline. Take care not to stretch the coil during maneuvers in order to preserve the circular memory. While still immersed in the heparinized saline, gently retract the distal tip of the coil into the introducer sheath approximately 1-2 cm. Verify that the coil distal tip is at the tapered end of introducer sheath.
- Lock the Twist-Lock mechanism with a clockwise rotation (see Figure 3) and insert the tapered distal end of the introducer sheath through the RHV and into the hub of the 2-tip infusion microcatheter until the sheath is firmly seated. A slight buckling of the introducer sheath indicates proper positioning. Tighten the RHV around the introducer sheath to prevent the back flow of blood, but not so tight as to damage the coil during its introduction into the catheter.
- Unlock the Twist-Lock mechanism and transfer the GDC 360° Detachable Coil into the catheter by advancing the delivery wire in a smooth, continuous manner. This procedure is best accomplished by two people: one to maintain the introducer sheath inside the RHV and 2-tip infusion catheter hub and one to hold the sheath straight and advance the coil.

Once the flexible portion of the delivery wire has entered the catheter shaft, loosen the RHV and remove the introducer sheath over the delivery wire's proximal end. Once completed, tighten the RHV around the delivery wire. Leaving the introducer sheath in place will interrupt normal infusion of flush solution and allow back flow of blood into the catheter. Do not discard the introducer sheath until after the GDC 360° Detachable Coil has been positioned and detached.

- Visually verify that the flush solution is infusing normally. Once confirmed, loosen the RHV enough to advance the delivery wire, but not so much as to compromise the continuous infusion.
- Advance the coil under fluoroscopy and position carefully at the desired site. If coil placement is unsatisfactory, slowly withdraw by pulling on the delivery wire, then slowly advance again to reposition the coil. If the coil size is inappropriate, remove and replace with an appropriately sized coil. (See Coil Selection section.)

#### Precautions

Advance and retract the coil slowly and smoothly, especially in tortuous anatomy. Remove the coil if unusual friction or “scratching” is noted. If friction is noted in a second coil, carefully examine both coil and the catheter for possible damage such as catheter shaft buckling or kinking, or improperly fused joint.

Do not advance the coil with force **if the coil becomes lodged** within or outside the 2-tip infusion catheter. Determine the cause of resistance and remove the system when necessary.

If resistance is encountered when withdrawing the delivery wire, draw back on the infusion catheter simultaneously until the delivery wire can be removed without resistance.

If resistance is noted during coil delivery, remove the catheter/coil system and check for possible damage to the catheter.

**Warning:** Do not rotate delivery wire during or after delivery of the coil into the aneurysm. **Rotating the GDC 360° Detachable Coil delivery wire may result in a stretched coil or premature detachment of the coil from the delivery wire, which could result in coil migration.**

- Continue to advance the GDC 360° Detachable Coil until the radiopaque proximal marker on the delivery wire is exactly distal to the proximal marker on the 2-tip infusion catheter (see Figure 4 below). Tighten O-ring on RHV to prevent movement of the delivery wire and proceed with detachment per instructions below.


Figure 4. Optimal Alignment of Radiopaque Markers

**Warning:** Verify repeatedly that the distal shaft of the catheter is not under stress before detaching the GDC® 360° Detachable coil.

Axial compression or tension forces could be stored in the 2-tip infusion catheter causing the tip to move during coil delivery. Catheter tip movement could cause the aneurysm or vessel to rupture.

**Warning:** Advancing the delivery wire beyond the catheter tip once the coil has been detached involves risk of aneurysm or vessel perforation.

**Warning:** If undesirable movement of the coil can be seen under fluoroscopy following coil placement and prior to detachment, remove the coil and replace with another more appropriately sized coil. Movement of the coil may indicate the coil could migrate once it is detached. Angiographic controls should also be performed prior to detachment to ensure that the coil mass is not protruding into the parent vessel.

#### DETACHMENT PROCEDURE

**Caution:** Use only the InZone® Detachment System (M00345100940 or M00345100950) or a Power Supply (M00345100740 or M00345100840) to detach GDC 360° Coils. Do not use other power sources.

#### InZone Detachment System

If using the InZone Detachment System (M00345100940 or M00345100950) to detach GDC 360° Coils, refer to the Directions for Use packaged with each InZone Detachment System and the section "Using the InZone Detachment System with GDC Detachable Coils or Matrix® Detachable Coils".

#### Power Supply

If using a Power Supply (M00345100740 or M00345100840) to detach GDC 360° Coils, refer to the following instructions:

#### WARNING

Increased detachment times may occur when:

- Delivery wire/catheter markers are not properly aligned.
- There is improper setup of continuous flush.
- Other embolic agents are present.
- Tungsten and/or other non-platinum coils are already present within the aneurysm. It is therefore recommended that GDC 360° Detachable Coils not be used in the presence of non-platinum coils.
- Detachable Coil Connecting Cables (M00345110240) have been resterilized. Detachable Coil Connecting Cables are supplied for one use only and should be discarded after each procedure.

**Caution:** Do not sterilize the Power Supply.

Power Supply (M00345100740 or M00345100840) is powered by two 9-Volt alkaline batteries which shall be replaced before each procedure. The batteries shall also be replaced if the Battery Check display illuminates during power supply operation. The batteries have an expected life of 4-6 hours.

Before proceeding with electrolytic detachment, make sure the Power Supply has been mounted to a 19 mm (0.75 in) diameter IV stand that has been securely fastened to the angio table. To mount power supply, pull back spring clip, align unit's concave section with pole, gently release clip. Note: If IV pole is too narrow, attachment adapter is packaged with power supply.

Refer to the Power Supply Directions for Use packaged with each device for Preventative Maintenance instructions as well as error and warning messages.

Ensure the Power Supply is mounted securely and contains fresh batteries.

When the GDC 360° Detachable Coil has been placed as desired, proceed with detachment per the following:

- Insert sterile uncoated hypodermic needle size 20 or 22 gauge (0.9 mm or 0.7 mm) at the groin.

**Caution:** Needle must not be coated.

- Remove the Detachable Coil Connecting Cables (M00345110240) from their pouch, and insert the plug of the sterile black connecting cable into the black terminal of the Power Supply.
- Clip the connector end of the same sterile black connecting cable onto the sterilized hypodermic needle (the patient return electrode). Wipe clean and dry the exposed proximal section of the delivery wire.
- Insert the plug of the sterile red connecting cable into the red terminal of the Power Supply.
- Clip the connector of the sterile red connecting cable onto the proximal end of the GDC 360° Detachable Coil delivery wire no more than 5 cm from the proximal end of the delivery wire. Assure the delivery wire is resting on a dry, clean surface and that the delivery wire remains clean and dry.


Figure 5. Detachable Coil Power Supply

**Caution:** Verify that red Power Supply connecting cable is connected to the GDC Detachable Coil's delivery wire, and the black Power Supply connecting cable is connected to the sterilized hypodermic needle. If these connections are reversed, GDC Coils will not detach.

- Confirm again under fluoroscopy that the proximal marker of the delivery wire is exactly distal to the proximal marker of the 2-tip infusion catheter (per Figure 4). (Position may have changed during attachment of the Detachable Coil Connecting Cables.)
- Press On/Off button located top right of the Power Supply (Figure 5). All displays will light and the system will perform a self check.
- If, during coil detachment, the patient experiences discomfort at the site of the patient return electrode, current flow can be reduced by pressing the CURRENT (touch-sensitive) icon. The default current setting is 1 mA. To select a different current setting, press the current button.

**Note:** Switching current ranges during a procedure:

- Has no effect on the functionality of the unit but may increase detachment time
- Will not reset timer
- WILL slightly affect the voltage required for that setting

- Again, visually verify that the continuous flush solution is infusing normally.
- Monitor Power Supply displays:
  - Refer to the Directions for Use packaged with each Power Supply for descriptions of the power supply indicators and displays.

**Caution:** In some ECG equipment, perturbations may be observed immediately before illumination of detachment indicator light on the Power Supply.

- Power supply circuitry is designed to monitor, detect, and signal the separation of the coil from the delivery wire. Once detachment is detected, the following should take place:
  - Current flow stops
  - All displays freeze, indicating the latest voltage, current levels and time immediately prior to separation
  - System beeps five times (1/2 sec on and 1/2 sec off)

At this point, the Power Supply has paused current flow to allow time for the clinician, utilizing fluoroscopy, to confirm coil detachment.

**Caution:** DO NOT PRESS THE ON/OFF BUTTON until after verification under fluoroscopy that the coil has detached.

- i. Once the detachment of the coil has been signaled, **VERIFY UNDER FLUOROSCOPY THAT THE COIL HAS DETACHED:** Slowly pull back on the delivery wire while watching fluoroscopy to make sure the coil does not move. In the unlikely event the coil moves, perform the following:
  - Advance the delivery wire to reestablish the wire and microcatheter marker alignment.
  - To resume current flow, press the CURRENT button on the Power Supply.
  - Verify coil separation under fluoroscopy as above. Repeat as necessary.
- m. Once coil detachment has been detected and fluoroscopically confirmed, disconnect the red connector from the delivery wire, and slowly withdraw the wire from the catheter.
- n. Turn system off by depressing On/Off button.
- o. Repeat the above steps if additional coil placement is required.

**Caution:** If the patient experiences pain at the site of the patient return electrode, or if detachment times are increasing, replace the return electrode with a new needle at a new insertion site.

- p. Once procedure is over, discard the Detachable Coil Connecting Cables and store Power Supply in a clean, dark, dry and secure location.

**Warning:** Detachable Coil Connecting Cables are intended for one use only. Do not resterilize and/or reuse. Resterilization could corrode the Detachable Coil Connecting Cables, resulting in increased detachment times.

**Warning:** Dispose of used Detachable Coil Connecting Cables and coil devices in accordance with hospital, administrative and/or local government policy.

#### MAGNETIC RESONANCE IMAGING (MRI) COMPATIBILITY (NEUROVASCULAR USE)

Non-clinical testing and analysis have demonstrated that when used in the neurovasculature, GDC® 360° Detachable Coils are MR Conditional. A patient with GDC 360° Detachable Coils can be safely scanned immediately after placement of the coils, under the following conditions:

- static magnetic field of 1.5 and 3.0 Tesla
- spatial gradient field up to 2500 Gauss/cm (25 Tesla/m)
- normal operational mode for gradients and SAR (maximum whole body averaged specific absorption rate (SAR) of lower than 2.0 W/kg and maximum head SAR of lower than 3.2 W/kg) for a total active MR scan time (with RF exposure) of 15 minutes or less per scan sequence.

In an analysis based on the temperature rises in non-clinical testing of spherical aneurysm models, the GDC 360° Detachable Coils were determined to produce an in-vivo temperature rise of 4°C or less for 15 minutes of MR scanning in normal mode operation in 1.5 T and 3 T MR systems. The GDC 360° Detachable Coils should not migrate in this MRI environment.

Temperature testing was not conducted in arteriovenous malformations or fistulae models.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

#### MAGNETIC RESONANCE IMAGING (MRI) COMPATIBILITY (PERIPHERAL USE)

Non-clinical testing and analysis have demonstrated that when used in peripheral vasculature, GDC 360° Detachable Coils are MR Conditional. A patient can be safely scanned immediately after placement of the coils, under the following conditions:

- static magnetic field of 1.5 and 3.0 Tesla
- spatial gradient field up to 2500 Gauss/cm (25 Tesla/m)
- normal operational mode for gradients and SAR (maximum whole body averaged specific absorption rate (SAR) of lower than 2.0 W/kg and maximum head SAR of lower than 3.2 W/kg) for a total active MR scan time (with RF exposure) of 15 minutes or less per scan sequence.

In an analysis based on the temperature rises in non-clinical testing of a bundle of 30 test coils with a length of 109 mm and the calculated SAR in the patient during an MR scan, the test coils were determined to produce an in-vivo temperature rise of 7.4°C or lower for 15 minutes of MR scanning in normal operational mode in 1.5 T and 3 T MR systems. The GDC 360° Detachable Coils should not migrate in this MRI environment.

The SAR limit of 2 W/kg applies only for coil placement in the torso. A reduction of SAR limits may be appropriate for coil placement in the legs or arms. MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

#### BATTERY CHANGE & DISPOSAL

Remove the battery access door. Replace spent batteries with two new 9-Volt alkaline batteries in the compartment as shown, noting CORRECT POLARITY per Figure 6. Snap battery door securely into place. Dispose of spent batteries in accordance with hospital, administrative and/or local government policy. Remove batteries when unit not in use.

**Caution:** Incorrect placement of batteries could result in unit malfunction.


Figure 6. Battery Compartment for Power Supply

#### PRECAUTIONS

High quality, digital subtraction **fluoroscopic road mapping is mandatory** to achieve safe catheterization of the aneurysm or vessel and correct placement of the first coil. *With smaller aneurysms this is particularly important.*

If resistance is encountered while withdrawing a coil which is at an acute angle relative to the catheter tip, it is possible to avoid coil stretching or breaking by carefully repositioning the distal tip of the catheter at the ostium of the aneurysm, or just slightly inside the parent artery. By doing so, the aneurysm and artery act to “funnel” the coil back into the catheter.

Take care not to puncture gloves or sterile drape while handling the delivery wire.

Multiple embolization procedures may be required to achieve the desired occlusion of some aneurysms or vessels.

#### WARNINGS

If coil repositioning is necessary, take special care to retract coil under fluoroscopy in a one-to-one motion with the delivery wire. If the coil does not move with a one-to-one motion, or repositioning is difficult, the coil has been stretched and could possibly break. Gently remove and discard both the catheter and coil.

Due to the delicate nature of the coils, the tortuous vascular pathways which lead to certain aneurysms and vessels, and the varying morphologies of intracranial aneurysms, a coil may occasionally stretch while being maneuvered. Stretching is a precursor to potential malfunctions, such as coil breakage and migration.

The long term effect of this product on extravascular tissues has not been established so care should be taken to retain this device in the intravascular space.

#### WARRANTY

Stryker Neurovascular warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Stryker Neurovascular’s control directly affect the instrument and the results obtained from its use. Stryker Neurovascular’s obligation under this warranty is limited to the repair or replacement of this instrument and Stryker Neurovascular shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Stryker Neurovascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **Stryker Neurovascular assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

**Caution:** The safety and performance characteristics of the GDC 360° System (GDC 360° Detachable Coils, Power Supply, InZone® Detachment System, delivery systems, and accessories), when used with another manufacturer’s devices (whether coils, coil delivery devices, catheters, guidewires, and/or other accessories) has NOT been established. Due to the potential incompatibility of non-Stryker Neurovascular components with the GDC 360° System, the use of another manufacturer’s device(s) with the GDC 360° System is not recommended.

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Copyright © 2013 Stryker
Stryker France S.A.S.
Stryker Australia Pty Ltd