

FDA Home³ Medical Devices⁴ Databases⁵

Class 2 Device Recall CSA Medical truFreeze Spray Kit




[6 510\(K\)](#) | [DeNovo⁶](#) | [Registration & Listing⁷](#) | [Adverse Events¹⁰](#) | [Recalls¹¹](#) | [PMA¹²](#) | [HDE¹³](#) | [Classification¹⁴](#) | [Standards¹⁵](#)
[CFR Title 21¹⁶](#) | [Radiation-Emitting Products¹⁷](#) | [X-Ray Assembler¹⁸](#) | [Medsun Reports¹⁹](#) | [CLIA²⁰](#) | [TPLC²¹](#) | [Inspections²²](#)

New Search

[Back to Search Results](#)

**Class 2 Recall
CSA Medical truFreeze Spray Kit**



Date Posted	April 23, 2015
Recall Status¹	Open
Recall Number	Z-1513-2015
Recall Event ID	70851²⁴
Premarket Notification 510(K) Number	K143625²⁵
Product Classification	Carrier, Ligature²⁶ - Product Code GEJ²⁷
Product	CSA Medical truFreeze Console- Cryosurgical Unit Cryogenic Surgical Device Model: CC3-01
Code Information	Serial Numbers: 01-00106 through 01-00201
Recalling Firm/Manufacturer	CSA Medical 91 Hartwell Ave Lexington, Massachusetts 02421-3137
Manufacturer Reason for Recall	TruFreeze Console caused a higher rate of liquid nitrogen (cryogen) to be delivered and may cause: stricture, scarring, bradycardia, or pneumothorax
FDA Determined Cause²	DESIGN: Software Design
Action	CSA Medical issued letter dated 3/25/15 advising users of the problem. Users provided with: The mitigations available to the active venting procedures coupled with the extremely unlikely probability of injury, rare risk, thus allowing the physician to continue with active venting procedures. In regards a passive venting procedure may not have sufficient mitigation to allow the physician to identify the potential hazard with sufficient time to preclude a potential for injury. Therefore, passive users will be instructed not to use the system until a software improvement is put into place. Additionally, no catheters or consoles will be shipped to passive venting users of truFreeze. A response form to be signed and returned confirming receipt of the notification. Questions contact: Stephen Mascioli, MD 781-538-4755 smascioli@csamedical.com.
Quantity in Commerce	82 units
Distribution	Nationwide
Total Product Life Cycle	TPLC Device Report²⁸

¹ For details about termination of a recall see [Code of Federal Regulations \(CFR\) Title 21 §7.55²⁹](#)

² Per FDA policy, recall cause determinations are subject to modification up to the point of termination of the recall.

510(K) Database [510\(K\)s with Product Code = GEJ and Original Applicant = CSA MEDICAL, INC.³⁰](#)

Links on this page:

1. <http://www.addthis.com/bookmark.php?u508=true&v=152&username=fdamain>