

Teleflex Medical IDA Business & Technology Park Dublin Road, Athlone Co. Westmeath, Ireland

24th August 2015

URGENT - FIELD SAFETY NOTICE

Type of Action:	Recall				
Teleflex Reference:	005-2015				
Flexi-Slip™ Endotracheal Tube Stylet					
Commercial Name	Product Code	Lot Numbers			
FLEXI-SLIP ENDOTRACHEAL TUBE STYLET WITH SOFT DISTAL TIP	502501	Refer to Appendix			
FLEXISLIP STYLET, STERILE PACK	503700-000060				
	503700-06	2			

Dear Customer,

Teleflex have initiated a voluntary Field Safety Corrective Action for the above listed products.

Description of the problem

Teleflex are recalling the products referenced above following receipt of reports of the plastic coating of the stylet splitting and/or breaking off of the stylet. This may result in a piece of plastic totally or partially occluding the patient's airway and impairing ventilation, or necessitating invasive removal procedures in order to prevent complications such as atelectasis or pneumonia. No patient injuries have been reported related to this issue.

FIELD SAFETY CORRECTIVE ACTION INSTRUCTIONS

ADVICE ON ACTION TO BE TAKEN BY MEDICAL STAFF

- We request that you check your inventory for product within the scope of this field action.
 Users should cease use and distribution of stock of the affected product batch and quarantine immediately.
- 2. If you do not have stock in scope of this field action as referred to in above table then mark the according checkbox on the Acknowledgement form (Appendix 1) and return the form to the fax number or e-Mail-address mentioned below.
- 3. If you have stock from the affected product as referred to in above table, mark the according checkbox on the Acknowledgement form (Appendix 1). Contact customer service by calling the phone number mentioned in Section 6 who will issue you with a return number. Write this return number into the respective field in the Acknowledgement form.
- 4. Complete 'Appendix 1' for all products in your possession and under control. Return this form immediately to Customer Service.
- 5. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

INSTRUCTION FOR DISTRIBUTORS OF AFFECTED PRODUCT

1. If you are a distributor, provide this field safety notice to all of your customers who have received product in scope of this Field Action. Your customer is then required to complete the acknowledgement form and return this to you.



- As a Distributor you are required to confirm to Teleflex that you have completed the field activity outlined above. Upon completion of your actions, please forward the completed Acknowledgement Form to Customer Service.
- Please be aware that all European Economic Area/Switzerland (EEA/CH) and Turkey
 Member State Competent Authorities in which Teleflex distribute directly will be notified by
 Teleflex.
- 4. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/TK area, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

Teleflex

Teleflex informs all customers, employees of Teleflex and distributors on this Field Action.

Transmission of this Field Safety Notice

This notice should be passed on to all persons who need to be aware within your organization or to any organization where the potentially affected devices have been transferred. Please consider end users, clinicians, risk managers, supply chain/distribution centres etc. in the circulation of this notice.

Maintain awareness of this notice until all required actions have been completed in your organisation

Contact reference person

Should you require any further information or support concerning this issue, please contact:

Customer Service

 Contact:
 Nicole Morawiec
 Telephone:
 +41 (0) 31 818 40 90

 FAX:
 +41 (0) 31 818 40 93
 E-mail: nicole.morawiec@teleflex.com

Please be advised that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities to which Teleflex distribute directly will be notified by Teleflex.

Teleflex is committed to providing high quality, safe and effective products. We sincerely apologize for any inconvenience this action may cause your operations. If you have any other questions, feel free to contact your local sales representative or Customer Service.

For and on behalf of Teleflex.

Karen Boylan

Karen Boylan VP Global, RAQA



Appendix 1

FIELD SAFETY CORRECTIVE ACTION

Teleflex Ref. 005-2015

Acknowledgement Form							
URGENT ATTENTION REQUIRED							
Return completed form immediately to:							
FAX:	+41 (0) 31 81	8 40 93	E-mail: nicole.n	norawiec@teleflex.com			
Please check applicate	ole box:						
☐ We confirm receipt of this FSN and completed the required actions contained therein. We confirm that our inventory does NOT include products affected by this Field Action.		☐ We confirm receipt of this FSN and completed the required actions contained therein. We confirm our inventory DOES include products affected by this Field Action. The use and further distribution of the affected products has been stopped. All products are on hold and the quantity stated below will be returned. Return Authorisation No					
Please CLEARLY prin	t the below returr	information	on:				
Name of Affected Products: FLEXI-SLIP ENDOTRACHEAL TUBE STYLET WITH SOFT DISTAL TI FLEXISLIP STYLET, STERILE PACK				DFT DISTAL TIP			
Product Number (Size)		Lot Number		Quantity (Returning)			
Return Instructions: • Please label produce • Include a copy of the Returns excluding ALL	is form (including F	RAN Numbe	er) with product returns.				
Institution Name - (Ho	spital, Health Car	e Organisa	tion, etc.)				
		Т					
Institution Address:			Email Address:				
			Di vi				
Form completed by:			Phone Number:				
Print Name :			Institution Stamp:				
Signature :							
Date:							



Appendix 2									
Product Code	Lot	Product Code	Lot	Product Code	Lot	Product Code	Lot	Product Code	Lot
	12EE20		13LG20	502501	15BG05		13GG07	503700-000060	14HE32
	12FE26		14AG25		15CG24		13GG33		14IE36
	12GE27		14BG27 14CG09 14CG19 14DG24		15DE18		13HG05		14JE43
	12GE30				15DG38		13IG21		14JG05
	12IE36				15EG08		13IG23		14KG11
	12IE37				15FG21		13JG06		14LG06
	12IE39		14FG03		15GG21		13KG04		15AG28
	12JE41		14FG21		12GE27		13LG06		15AG34
	12KE48		14GG01		12HE33		13LG20		15BG05
502501	12LG25	502501	14GG03 14HE32	503700-000060	12IE37	503700-	13LG27		15CG24
	12LG29	302301			12JE40	000060	14AG04		15DE18
	13AG21		14HE34		12KE48		14BG24		15DG38
	13BG26		14IE36		12LE50		14CG09		15EG12
	13BG36		14JE43		13AT24		14DG01		15FG21
	13DG06		14JG05		13AT43		14DG24		15GG21
	13DG24		14KG11		13CG05		14FE25	503700-06	12GE27
	13GG06		14KG27		13DG09		14FG03		12HE33
	13GG33		14LG06		13EG19		14FG06		12IE37
	13KG04		15AG28		13FG19		14FG21		13AT43
	13LG06		15AG34		13GG06		14GG01		