## Information Urgente de sécurité

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
# UPDATING CLINICIAN PROGRAMMER SOFTWARE

### STEP 1: PREPARATION

**Power:** Ensure the clinician programmer is plugged in or has more than 25 percent battery power prior to installing any updates.

**Wi-Fi:** Connect to Wi-Fi.

**Close Apps:** Close any open applications.

**Time:** Depending on the network bandwidth, this process could take up to 30 minutes.

### STEP 2: UPDATING AGENT TO HUB

Medtronic applications can be updated using the Hub app , which replaced Agent. If you see the Agent app, please follow these instructions to update to Hub.


1. Navigate to Google Play Store app.
2. Tap menu icon.
3. Select **My Apps & Games**.
4. Update **Intelligent Hub** app.

### STEP 3: APPLICATION UPDATE PLANNING

1. Click on **Hub** app. Then select **App Catalog** to view apps available for installation or update.

2. Only update therapy apps that you are educated on how to use.
3. Determine which therapy apps need updating.
4. Maintain up-to-date supporting applications. Check for updates for:
  - PDSApplication.** This is the Patient Data Services Application that supports reporting.
  - CommManager.** This software supports the communicator. Updates to the communicator require use of the communicator cable. See further instructions on second page.

*Note: More than one application can update at the same time.*

### STEP 4: APPLICATION UPDATE PROCESS

1. While in **App Catalog**, click on **Update** (if an update is available) to update the apps you use.
2. Click on the **Confirm Installation** pop-up.

The screen will display **Processing** while the application is updating and will not change when installation is complete.

3. To check for installation completion, go to **Managed Apps** in Hub.
4. Once installation completes, confirm the correct software version in the Managed Apps or in the app's "About" section.

**Medtronic**  
Further, Together

## UPDATING COMMUNICATOR AND COMM MANAGER APP

After updating the communicator software, it is **crucial** to pair the communicator to the clinician programmer because the original coupling is lost.

### STEPS FOR UPDATING THE COMMUNICATOR

- 1 Remove tablet cover and any accessories preventing USB connection of the communicator cable.
- 2 After updating **Comm Manager** in the **App Catalog**, open the **Communicator Updater** application.
- 3 Power the communicator on and connect it to the tablet with the pairing cord.
- 4 Select **Update (or Recover) Communicator**.

- While the communicator is updating, **DO NOT DISCONNECT THE COMMUNICATOR**.
  - The update will take several minutes to complete. Please wait until the **Update Successful** screen appears.
- 1 Once you have reached the **Update Successful** screen, select **Done** and proceed to communicator pairing.

### PAIRING THE COMMUNICATOR

Ensure the communicator is powered on and connected to the tablet. An implantable device is not necessary for pairing the communicator.

- 1 Launch the therapy application and attempt to connect to an implantable device.
- 2 If the screen says, "Searching for Device" and not "Searching for Communicator," the communicator is now successfully paired.
  - The LED between the communicator and tablet icons will be solid green once connected.
  - If you encounter any issues, contact Digital Connectivity.

## SPECIALIZED IT SUPPORT

Medtronic Digital Connectivity handles non-therapy IT-related questions such as:

- Password resets
- Lost/Stolen reports
- Connectivity issues
- Cellular/Wi-Fi
- Bluetooth
- Printing issues (not related to Model 8840)
- Software updates
- Unresponsive Samsung device programmer issues (e.g., not powering up)
- Unable to download therapy app

### HOW TO CONTACT DIGITAL CONNECTIVITY SUPPORT:

- Technical Services will provide live support during the hours of 9 AM to 5 PM CET.
- On-call support will be provided 24/7 (English only).
- Contacts details: [rs.tsneuro@Medtronic.com](mailto:rs.tsneuro@Medtronic.com) or +31 (0)45 566 88 44

### USER ID CARD

Use the user ID card that came with the clinician programmer when calling Digital Connectivity.

See the device manual for detailed information regarding the instructions for indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [medtronic.eu](http://medtronic.eu)

For applicable products, consult instructions for use on [www.medtronic.com/manuals](http://www.medtronic.com/manuals). Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat Reader® with the browser.

## Medtronic

Europe  
Medtronic International Trading Sàrl  
Route du Molliou 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

United Kingdom/Ireland  
Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 (0)1923 212213  
Fax: +44 (0)1923 241004

UC201912239EE © Medtronic 2019.  
All rights reserved. Printed in Europe.

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S.A.S au capital de 6 564 721,44 € - R.C.S. Nanterre 722 008 232 - A.P.E. 4652 Z - T.V.A. : FR 61 722 008 232