

Ardon, le 8 juin 2021  
Courriel avec accusé de réception

**Destinataires :**

Correspondant Local de Matéiovigilance.

↳ Diffusion à réaliser pour information auprès :  
du Directeur de l'Etablissement de Santé,  
de tous les utilisateurs.

**Dispositifs médicaux concernés :**

**VKMO 11000** - Oxygénéateur QUADROX-i Neonatal HMO 11000  
avec Réservoir de cardiectomie veineux à coque rigide VHK 11000  
de référence\* :

- **701049279**,
- **701070444**.

**Objet :**

Intégrité de l'emballage stérile des bouchons étanches  
potentiellement compromise.

*Division ACT - Acute Care Therapies*



- Exemple de VHK -



- Exemple de VKMO -

Madame, Monsieur,

Par la présente nous souhaitons vous informer d'une action initiée par le fabricant Maquet Cardiopulmonary GmbH, Allemagne, concernant un certain nombre de dispositifs VHK Neonatal (Réservoir de cardiectomie veineux pédiatrique à coque rigide) et VKMO Neonatal (QUADROX-i en combinaison avec un Réservoir de cardiectomie veineux à coque rigide).

En effet d'après nos bases de traçabilité, nous avons déterminé que votre établissement avait réceptionné un ou plusieurs de ces dispositifs.

Nous vous transmettons par conséquent cette notification de sécurité (traduction en français) : vous pourrez ainsi prendre connaissance des informations relatives aux mesures à prendre par votre établissement.

Par ailleurs, nous vous saurions gré de bien vouloir **compléter le formulaire d'Accusé de réception Client** ci-joint, et **nous le retourner dans les meilleurs délais** par courrier électronique (qrc.fr@getinge.com) ; même si vous ne possédez plus de produits concernés par cette notification.

Pour le cas où certains dispositifs aient été transmis à une autre organisation, veuillez lui faire suivre ce courrier et en informer *Getinge France* par le biais de ce même formulaire.

**A titre d'information, les établissements de santé concernés par la présente notification de sécurité recevront ce courrier de *Getinge France* par courriel, avec accusé de réception.**

Cette notification a également fait l'objet d'une information auprès de l'ANSM.

L'ensemble de l'équipe *Getinge France - Division Acute Care Therapies* reste à votre entière disposition pour tout renseignement complémentaire relatif au contenu de cette lettre, et pour toute difficulté rencontrée lors de la mise en œuvre de cette action corrective.

Dans cette attente, nous vous prions de croire, Madame, Monsieur, en l'expression de nos sincères salutations.

---

Bénédicte Parisot  
Directrice QRC France  
Getinge France

\* Les références indiquées en objet sont les seules concernées en France.

Pièces jointes :

- Notification de sécurité (FSCA) - Maquet Cardiopulmonary GmbH (FSCA-2021-05-07 - traduction).
- Formulaire Accusé de réception Client.
- Annexe I - Liste des produits concernés (traduction).

08-06-2021

**URGENT - NOTIFICATION DE SÉCURITÉ**

**Objet :** FSCA-2021-05-07 Problème d'intégrité de la barrière stérile des accessoires "bouchons étanches" du VHK et du VKMO Neonatal

**Produits concernés :**

- Combinaison de Réservoir de cardiectomie veineux Neonatal à coque rigide et d'Oxygénateur QUADROX-i Neonatal avec ou sans revêtement
- Réservoir de cardiectomie veineux Neonatal à coque rigide

Identification	Référence	Désignation produit
VHK 11000	701048596	Réservoir de cardiectomie veineux Neonatal à coque rigide
BO-VHK 11000	701051430	Réservoir de cardiectomie veineux Neonatal à coque rigide, revêtement SOFTLINE
BO-VHK 11000-J	701054173	Réservoir de cardiectomie veineux Neonatal à coque rigide, revêtement SOFTLINE, Japon
VKMO 10000	701050109	Oxygénateur QUADROX-i Neonatal HMO 10000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000
BO-VKMO 10000	701053443	Oxygénateur QUADROX-i Neonatal HMO 10000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000, SOFTLINE
VKMO 11000	701049279	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000
BO-VKMO 11000	701053444	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000, SOFTLINE
BO-VKMO 11000 -J	701063847	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000, SOFTLINE, Japon
VKMO 10000	701070440	Oxygénateur QUADROX-i Neonatal HMO 10000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000
VKMO 11000	701070444	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000
BO-VKMO 11000	701071077	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000, SOFTLINE

**Numéros de lot concernés :** Voir Annexe I ci-jointe - Liste des produits concernés

Très chers Clients,

Le Réservoir de cardiectomie veineux Neonatal à coque rigide (VHK), tout comme sa combinaison avec l'Oxygénateur QUADROX-i Neonatal (VKMO), ont été développés pour être utilisés lors des interventions de chirurgie cardiaque. Ils sont utilisés dans le circuit de circulation extracorporelle pour collecter, stocker et filtrer le sang aspiré du champ opératoire ainsi que le sang veineux drainé du patient.

Lors des tests de vérification de la conception des VKMO et VHK Neonatal, une altération potentielle de l'emballage stérile de l'accessoire "bouchons étanches" a été identifiée. La Figure 1 montre la poche concernée, contenant les bouchons étanches.



Figure 1 : Poche stérile contenant les bouchons étanches.

Les Figures 2 et 3 illustrent le VKMO et les accessoires placés dans l'insert en mousse de polystyrène à l'intérieur de la boîte en carton.



Figure 2 : VKMO emballé dans un insert en mousse de polystyrène



Figure 3 : La poche enroulée qui contient les bouchons étanches dans l'insert en mousse de polystyrène

Les bouchons étanches concernés sont emballés dans un emballage stérile individuel placé en dehors du VKMO/VHK. Les bouchons étanches sont utilisés pour remplacer les bouchons aérés pendant les procédures de drainage veineux assisté par le vide (VAVD) et le drainage assisté par le vide des plaies post-opératoire. Ces bouchons doivent être stériles dans les deux applications, car ils sont en contact avec le flux sanguin.

L'exposition du flux sanguin ou du champ stérile (opératoire) à des bouchons étanches non stériles pourrait entraîner les conséquences immédiates et/ou à long terme (blessures ou maladies) suivantes :

- Exposition du champ (opératoire) stérile à un produit non stérile susceptible de compromettre l'ensemble du champ opératoire.
- Exposition du système vasculaire d'un patient à un produit non stérile entraînant une inflammation et/ou une infection.
- Septicémie due à la propagation non contrôlée et inconnue d'une infection dans d'autres régions du corps.

Maquet Cardiopulmonary n'a reçu aucune réclamation ni aucun rapport d'événements indésirables dus à des dommages du système de barrière stérile des bouchons étanches des produits susmentionnés.

**Action corrective :**

- Veuillez isoler et renvoyer immédiatement tous les produits concernés de votre stock à votre représentant Getinge local.

**Étant donné qu'aucun produit de remplacement n'est disponible, vous voudrez bien en outre tenir compte des éléments suivants :**

- Si l'emballage contenant les bouchons étanches n'est pas ouvert et mis au rebut et que les bouchons étanches ne sont pas nécessaires, par exemple si le drainage veineux par gravité est préféré au drainage veineux assisté par le vide (VAVD) ou si aucun drainage assisté par le vide de plaie n'est requis sur la base de l'expertise et de l'appréciation clinique, les produits peuvent être utilisés en toute sécurité.

**Recommandations sur les actions à entreprendre par l'utilisateur :**

- Selon nos données de surveillance, votre stock actuel peut inclure des produits concernés par cette action. Veuillez vous référer à l'Annexe I - Liste des produits concernés.
- **Attention** : Veuillez comparer les numéros de lot des étiquettes de votre stock avec les lots concernés conformément à l'Annexe I ci-jointe - Liste des produits concernés.
- Veuillez compléter et signer l'Accusé de réception Client ci-joint, puis le retourner à votre représentant Getinge local.
- Veuillez contacter votre représentant Getinge local afin d'obtenir l'avis correspondant.

**Documents référencés/ pièces jointes :**

- Annexe I Liste des produits concernés
- Accusé de réception Client

**Transmission de cette notification de sécurité :**

- Cette notification doit être transmise à l'ensemble des personnes impliquées au sein de votre établissement et à toute organisation où les dispositifs potentiellement affectés ont été transférés.
- Gardez toujours cet avis à l'esprit, tout comme les mesures qui en découlent, et ce pendant un certain temps afin de s'assurer de l'efficacité de l'action corrective.

Nous tenons à nous excuser pour la gêne susceptible d'être occasionnée, et ferons tout notre possible pour que cette action soit effectuée le plus rapidement possible.

Tel que requis, nous transmettrons cette notification aux Autorités compétentes pertinentes.

**Si vous avez des questions supplémentaires ou souhaitez obtenir de plus amples informations, n'hésitez pas à contacter votre représentant Getinge local.**

Sincères salutations,

**Managing Director**

**Safety Officer**

Maquet Cardiopulmonary GmbH  
Kehler Str. 31  
76437 Rastatt  
ALLEMAGNE

**Référence FA :** FSCA-2021-05-07  
**Objet de la FA :** FSCA-2021-05-07 Problème d'intégrité de la barrière stérile des accessoires "bouchons étanches" du VHK et du VKMO Neonatal  
**Produits concernés\* :**

- Combinaison de Réservoir de cardiectomie veineux Neonatal à coque rigide et d'Oxygénateur QUADROX-i Neonatal avec ou sans revêtement
- Réservoir de cardiectomie veineux Neonatal à coque rigide

Identification	Référence	Désignation produit
VKMO 11000	701049279	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000
VKMO 11000	701070444	Oxygénateur QUADROX-i Neonatal HMO 11000 avec Réservoir de cardiectomie veineux Neonatal à coque rigide VHK 11000

\* Les références indiquées sont les seules concernées en France.

Vous trouverez les lots concernés à l'Annexe I - Liste des produits concernés ci-joint.

Par la présente, je soussigné accuse réception de la notification ci-dessus mentionnée, et atteste avoir compris les actions requises dans ce cadre. Je confirme également avoir distribué cette Notification de sécurité au personnel concerné.

- Nous n'avons pas de produit concerné en **stock**.
- Oui, nous avons un/des produits concernés en **stock** et nous le(s) retournons à Getinge.
- Oui, nous avons un/des produits concernés en **stock** et nous l(es) utiliserons sans capuchons étanches. Le(Les) produit(s) concerné(s) n'a(n'ont) pas été en contact avec les capuchons étanches ou avec l'emballage de ces derniers.

Si oui :

Référence	Numéro de lot	Désignation	Quantité

FRANCE

Pays

Hôpital/clinique (adresse complète)

Date

Nom (fonction)

Signature

Document à renvoyer dûment complété :

Commentaires :

Par courriel : [qrc.fr@getinge.com](mailto:qrc.fr@getinge.com)

Adresse : [Getinge France \(Département QRC\)](#)  
 Parc de Limère - Av. de la Pomme de Pin  
 CS 10008 ARDON  
 45074 ORLEANS CEDEX 2



**La présente Annexe est une pièce jointe relative à la notification de sécurité de référence FSCA-2021-05-07**

Le tableau ci-dessous liste l'ensemble des numéros de lots de dispositifs distribués, impactés par la Notification de sécurité FSCA-2021-05-07.

Pays	Référence	Désignation Produit	Lot
Albania	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Argentina	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
Argentina	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309191
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261636
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261638
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261641
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261789
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268216
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268335
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268342
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271866
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281430
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92283354
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286829
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92287154
Argentina	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292949
Austria	701048596	VHK 11000 Neonatal Venous Hardshell Cardiotomy Reservoir	92295088



Pays	Référence	Désignation Produit	Lot
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268330
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271860
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271862
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92275269
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281428
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92287651
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292933
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292934
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309150
Brazil	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309191
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271860
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669

Pays	Référence	Désignation Produit	Lot
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309150
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309191
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309283
Canada	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309356
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261635
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261788
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268218
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268334
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268341
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268342
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281430
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292945
Canada	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292949
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92272926
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92273502
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92279802
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92279803
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92282088
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92286354
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290769
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290770

Pays	Référence	Désignation Produit	Lot
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290771
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92309112
Canada	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	3000144590
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261627
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261632
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268213
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268214
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268330
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268331
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271861
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271862
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92275269
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281428
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92287651
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92290475
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292933
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669

Pays	Référence	Désignation Produit	Lot
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305671
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92308546
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309150
China	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309356
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261641
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261788
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268188
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268218
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268345
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271863
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271864
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281430
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281433
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286827
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286828
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92287154
China	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000w	92295091
Colombia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
Ecuador	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281428
Estonia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
Estonia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669

Pays	Référence	Désignation Produit	Lot
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261630
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261632
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268213
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268329
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271860
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271862
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92275269
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286825
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92290475
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92291433
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
France	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
France	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305401
France	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403

Pays	Référence	Désignation Produit	Lot
France	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305404
France	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261630
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268329
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271860
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271861
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271862
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281428
Germany	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92291433
Germany	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268345
Germany	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92255923
Germany	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92285488
Germany	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290770
Germany	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290771
Germany	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305401



Pays	Référence	Désignation Produit	Lot
Germany	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
Germany	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92301531
Germany	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92301532
Germany	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92310473
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261635
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261636
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261640
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261789
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268188
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268215
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268217
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268218
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268219
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268333
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268334
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268335
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271868
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92273127
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92283354
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286826
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286829



Pays	Référence	Désignation Produit	Lot
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286830
India	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92287154
Ireland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
Ireland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261630
Ireland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
Ireland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
Ireland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286825
Ireland	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403
Ireland	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
Italy	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Italy	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
Italy	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92290475
Italy	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
Italy	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92272926
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92279803
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92279804
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92282088
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92286354
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92288038
Italy	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290769

Pays	Référence	Désignation Produit	Lot
Italy	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305401
Italy	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403
Italy	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305404
Italy	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
Italy	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92301532
Italy	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92309199
Italy	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92309200
Italy	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92310473
Japan	701051430	BO-VHK 11000 Neonatal Venous Hardshell Cardiotomy Reservoir, SOFTLINE coated	92298798
Japan	701051430	BO-VHK 11000 Neonatal Venous Hardshell Cardiotomy Reservoir, SOFTLINE coated	92308927
Japan	701051430	BO-VHK 11000 Neonatal Venous Hardshell Cardiotomy Reservoir, SOFTLINE coated	92309785
Japan	701051430	BO-VHK 11000 Neonatal Venous Hardshell Cardiotomy Reservoir, SOFTLINE coated	92311338
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92286365
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92288895
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92288896
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92292624
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92292625
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92292626
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92297016
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92297017

Pays	Référence	Désignation Produit	Lot
Japan	701053443	BO-VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92297018
Japan	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92255923
Japan	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92275882
Japan	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92285488
Japan	701054173	BO-VHK 11000-J Neonatal Venous Hardshell Cardiotomy Reservoir, SOFTLINE coated, Japan	92284145
Japan	701063847	BO-VKMO 11000 -J QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE, Japan	92286236
Lebanon	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
Lebanon	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92283354
Malaysia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Malaysia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286829
Mexico	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
Mexico	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261627
Mexico	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
Mexico	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
Netherlands	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
Netherlands	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92287651
Netherlands	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268341
Netherlands	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92255923
Netherlands	701071077	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92301531
Paraguay	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011

Pays	Référence	Désignation Produit	Lot
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261636
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261637
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261638
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268188
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271866
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92285801
Poland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286206
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261632
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268213
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286825
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92291433
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
Portugal	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669

Pays	Référence	Désignation Produit	Lot
Portugal	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403
Portugal	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305404
Portugal	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261638
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261641
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268215
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268216
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268218
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268335
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268342
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271864
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271868
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281433
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292949
South Korea	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292950
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268329
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332

Pays	Référence	Désignation Produit	Lot
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271860
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271861
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92275269
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286825
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92290475
Russia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
Russia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261641
Russia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268216
Russia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268334
Russia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268341
Russia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268342
Saudi Arabia	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
Saudi Arabia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261635
Saudi Arabia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
Saudi Arabia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92283354
Saudi Arabia	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92283678
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261627
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261632



Pays	Référence	Désignation Produit	Lot
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268213
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
South Africa	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292934
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261627
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261631
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261632
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92262011
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268329
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268332
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271859
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271861
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281428



Pays	Référence	Désignation Produit	Lot
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286825
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92290475
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92291433
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
Spain	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268215
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268341
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268345
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281430
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92283354
Spain	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286829
Spain	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290770
Spain	701053444	BO-VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000, SOFTLINE	92290771
Spain	701070440	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92301536
Spain	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403
Spain	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305404
Spain	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
Switzerland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268329

Pays	Référence	Désignation Produit	Lot
Switzerland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271860
Switzerland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271862
Switzerland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292934
Switzerland	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92294673
Switzerland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261640
Switzerland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271866
Switzerland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
Switzerland	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292945
Switzerland	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403
Taiwan	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
Taiwan	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292934
Taiwan	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305357
Taiwan	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669
Thailand	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261789
Turkey	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305669
United Arab Emirates	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268342
United Arab Emirates	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
United Arab Emirates	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92306910
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92260982
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261618

Pays	Référence	Désignation Produit	Lot
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261630
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92268326
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92275269
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92286825
United Kingdom	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92291433
United Kingdom	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261635
United Kingdom	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92261788
United Kingdom	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92274833
United Kingdom	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92281430
United Kingdom	701050109	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92292945
United Kingdom	701070440	VKMO 10000 QUADROX-i Neonatal HMO 10000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92301536
United Kingdom	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92305403
United Kingdom	701070444	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92309362
Uruguay	701049279	VKMO 11000 QUADROX-i Neonatal HMO 11000 with Neonatal Venous Hardshell Cardiotomy Reservoir VHK 11000	92271857