

FDA Home³ Medical Devices⁴ Databases⁵

Class 2 Device Recall Temperature Sensor Catheter

|Recalls11|PMA12|HDE13|Classification14|Standards15 Registration & 510(k) | DeNovo⁸ Adverse

Listing⁹ Events¹⁰

CFR Title 21¹⁶|Radiation-Emitting Products¹⁷|X-Ray Assembler¹⁸|Medsun Reports¹⁹|CLIA²⁰|TPLC²¹

New Search

Back to Search Results

Class 2 Device Recall Temperature Sensor Catheter

See Related Information

Date Initiated by Firm

SuperSearch

July 11, 2018

Create Date

September 20, 2018

Recall Status¹

Open³, Classified

Recall Number

Z-3208-2018

Recall Event ID

8074723

510(K)Number

K933400²⁴

Product Classification

Catheter, urological²⁵ - Product Code KOD²⁶

Product

P400 PREM U/M 16FR TMPSNS, Item Code P4P16TS

Code Information

1814906864

Recalling Firm/ Manufacturer

COVIDIEN LLC

15 Hampshire St

Mansfield MA 02048-1113

For Additional Information Contact Customer Advocacy 800-292-9332

Manufacturer Reason

for Recall

The temperature sensor catheters may be defective in that they will show a lower body

temperature measurement than the actual temperature of the patient.

FDA Determined

Cause 2

Nonconforming Material/Component

Action

On 7/11/18, an Urgent Medical Device Recall letter was distributed to customers. The letter instructed customers to do the following: 1. Please quarantine any unused product of the Covidien item codes and lot numbers listed above (see Attachment A to locate). 2. If you have distributed any of the affected Covidien catheters or kits, then please promptly forward a copy of this letter to those recipients. 3. All unused product of affected item codes and lot numbers must be returned. 4. Please return unused affected product as explained below. You must complete the Recalled Product Return Form (Attachment B) regardless of whether

or not you have inventory of affected product.

Quantity in Commerce

27515 total

Distribution

The products were distributed to the following US states: CA, CO, FL, GA, IA, IL, IN, KY, LA,

MA, MD, MO, NC, NJ, NY, OH, OR, PA, RI, SC, TN, TX, VA, and WA.

Total Product Life Cycle

TPLC Device Report²⁷

updated as the status changes.

¹ A record in this database is created when a firm initiates a correction or removal action. The record is updated if the FDA identifies a violation and classifies the action as a recall, and it is updated for a final time when the recall is terminated. Learn more about medical device recalls²⁸

² Per FDA policy, recall cause determinations are subject to modification up to the point of termination of the recall. ³ The manufacturer has initiated the recall and not all products have been corrected or removed. This record will be