



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

Urgent Field Safety Notice Unicryl Surgical Suture

Date: 13/11/2016

Unimed Reference: FSN01

Dear NCMDR Team

This letter to inform you that we have identified that the Suture Breaking Strength in few incidents are non compliant to European Pharmacopeia requirements.

We are initiating a voluntary recall for the following batches:

Serial	Batch#	Item Code	Description
1	10777-06	DW8902U	Unicryl (ud) 19mm curv.cutt. 2(3/0) 75cm
2	10943-41	DW8902RU	Unicryl (ud) 19mm curv.rev.cutt. 2(3/0) 75cm
3	11090-39	DW8937RU	Unicryl (ud) 24mm curv.rev.cutt. 2(3/0) 75cm
4	10705-16	DW8832RSU7	Unicryl (ud) 26mm curv.rev.cutt. 2(3/0) 75cm
5	11026-28	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
6	11130-02	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
7	10705-17	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
8	10831-01	DW8572U	Unicryl (ud) 60mm str.cutt. 2(3/0) 75cm
9	11043-01	DW8572U	Unicryl (ud) 60mm str.cutt. 2(3/0) 75cm
10	11024-01	DW8759R	Unicryl 13mm curv.rev.cutt. 0.7(6/0) 45cm
11	11118-16	SD8214	Unicryl 18mm curv.cutt 1(5/0)75cm
12	10754-01	DW8900RS	Unicryl 19mm curv.rev.cutt. 1(5/0) 45cm
13	11090-37	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
14	10850-49	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
15	10850-33	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
16	10850-49	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
17	10841-23	SE 540	Unicryl 22mm 1/2c.r.b. 1.5(4/0) 75cm
18	10878-01	S 9890	Unicryl 22mm 1/2c.r.b. 3(2/0) 75cm
19	10841-18	DW8937S	Unicryl 24mm curv.cutt. 2(3/0) 75cm
20	10841-14	DW8938	Unicryl 24mm curv.cutt. 3(2/0) 70cm
21	10705-19	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
22	11148-04	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
23	10879-12	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
24	10920-01	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
25	10706-33	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
26	10943-29	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
27	10777-14	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
28	10900-03	SD8300	Unicryl 26mm 5/8c.r.b.3(2/0) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

29	11127-01A	DW8833	Unicryl 26mm curv.cutt. 3(2/0) 75cm
30	10850-47	DW8832RS	Unicryl 26mm curv.rev.cutt. 2(3/0) 45cm
31	10988-18	DW8832RS7	Unicryl 26mm curv.rev.cutt. 2(3/0) 75cm
32	10846-40	DW8952	Unicryl 30mm 1/2c.r.b. 2(3/0) 75cm
33	11090-42	DW8952	Unicryl 30mm 1/2c.r.b. 2(3/0) 75cm
34	10846-39	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm
35	10879-26	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm
36	10994-02	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm
37	10841-02	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm
38	10879-28	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm
39	11090-45	DW8955	Unicryl 30mm 1/2c.r.b. 4(1) 75cm
40	10841-06	DW8955	Unicryl 30mm 1/2c.r.b. 4(1) 75cm
41	10841-15	DW8813S7	Unicryl 30mm curv.cutt. 3(2/0) 75cm
42	11127-03A	SD8668	Unicryl 35mm curv.cutt. 3.5(0) 75cm
43	11222-03	SD8668	Unicryl 35mm curv.cutt. 3.5(0) 75cm
44	10801-08	SD8163	Unicryl 35mm j.r.b. 4(1) 75cm
45	10841-07	S 9893S9	Unicryl 37mm 1/2c.r.b. 4(1) 90cm
46	10841-22	SD8665	Unicryl 37mm curv.cutt. 3(2/0) 75cm
47	10850-19	DW8525	Unicryl 40mm 1/2c.cutt. 4(1) 75cm
48	10890-01	DW8526S	Unicryl 40mm 1/2c.cutt. 5(2) 70cm
49	11154-01	DW8513	Unicryl 40mm 1/2c.r.b. 3(2/0) 75cm
50	10879-09	DW8513	Unicryl 40mm 1/2c.r.b. 3(2/0) 75cm
51	10370-62K	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
52	10370-62J	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
53	10370-62N	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
54	10370-62L	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
55	11101-61E	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
56	11101-61F	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
57	10706-20	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
58	10988-16A	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
59	10988-16	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
60	10777-12	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm
61	11205-05	DW8514S	Unicryl 40mm 1/2c.r.b. 3.5(0) 90cm
62	10706-21	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm
63	10846-43	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm
64	10705-13	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm
65	10887-22A	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm
66	10879-01	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm
67	11044-15	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm
68	10943-31	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm
69	10988-17	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm
70	10850-20	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

71	11116-33	DW8514H	Unicryl 40mm 1/2c.r.b.hvy. 3.5(0) 75cm
72	11175-38	DW8524R	Unicryl 40mm 1/2c.rev.cutt. 3.5(0) 75cm
73	10943-32	DW8525R	Unicryl 40mm 1/2c.rev.cutt. 4(1) 75cm
74	10826-05	DT8415H	Unicryl 40mm 1/2c.taper.cut.hvy. 3(2/0) 75cm
75	11127-02B	SD8189	Unicryl 40mm curv.cutt. 4(1) 75cm
76	10904-01	DW8854	Unicryl 45mm 1/2c.r.b.hvy. 3.5(0) 75cm
77	10826-02	DW8855	Unicryl 45mm 1/2c.r.b.hvy. 4(1) 75cm
78	10826-03	DW8856	Unicryl 45mm 1/2c.r.b.hvy.5(2) 75cm
79	10846-47	DW8856	Unicryl 45mm 1/2c.r.b.hvy.5(2) 75cm
80	10777-09A	SD8070	Unicryl 48mm 1/2c.cutt. 4(1) 75cm
81	11118-22	SD8070	Unicryl 48mm 1/2c.cutt. 4(1) 75cm
82	10850-29	SD8118	Unicryl 48mm 1/2c.r.b. 3.5(0) 90cm
83	10879-07	SD8118	Unicryl 48mm 1/2c.r.b. 3.5(0) 90cm
84	10879-05	S 9849	Unicryl 48mm 1/2c.r.b. 4(1) 75cm
85	10850-32	S 9849	Unicryl 48mm 1/2c.r.b. 4(1) 75cm
86	10841-09	SD8061	Unicryl 48mm 1/2c.r.b. 4(1) 90cm
87	10879-06	DW8554	Unicryl 50mm 1/2c.r.b. 3.5(0) 75cm
88	10841-20	DW8571S	Unicryl 60mm str.cutt. 1.5(4/0) 67cm
89	10841-21	DW8573	Unicryl 60mm str.cutt. 3(2/0) 75cm
90	10777-13	DW8573	Unicryl 60mm str.cutt. 3(2/0) 75cm
91	10841-19	DW8627D	Unicryl 6mm curv.spat.d/a 0.4(8/0) 30cm
92	10839-13	DW8628DS	Unicryl 6mm curv.spat.d/a 0.5(7/0) 45cm
93	11090-28	DW8504	Unicryl cut lengths 3.5(0) 1.5m
94	11174-10	DW8505	Unicryl cut lengths 4(1) 1.5m
95	11044-12	DW8503	Unicryl cut lengths 3(2/0) 1.5 metres
96	10708-04	DW8503	Unicryl cut lengths 3(2/0) 1.5 metres
97	10841-16	DW8503	Unicryl cut lengths 3(2/0) 1.5 metres
98	10423-37	DW8860SU	Unicryl (ud) 16mm curv.cutt. 1(5/0) 45cm
99	10500-36	DW8861SU	Unicryl (ud) 16mm curv.cutt. 1.5(4/0) 45cm
100	10501-01	DW8861SU	Unicryl (ud) 16mm curv.cutt. 1.5(4/0) 45cm
101	10308-11	DW8861RU	Unicryl (ud) 16mm curv.rev.cutt. 1.5(4/0) 75cm
102	10308-10	DW8860RU	Unicryl (ud) 16mm curv.rev.cutt.1(5/0) 75cm
103	10523-01	DW8900U	Unicryl (ud) 19mm curv.cutt. 1(5/0) 75cm
104	10267-73	DW8901SU	Unicryl (ud) 19mm curv.cutt. 1.5(4/0) 45cm
105	10523-02	DW8901U	Unicryl (ud) 19mm curv.cutt.1.5(4/0) 75cm
106	10633-16	DW8901U	Unicryl (ud) 19mm curv.cutt.1.5(4/0) 75cm
107	10140-80	DW8937RU	Unicryl (ud) 24mm curv.rev.cutt. 2(3/0) 75cm
108	10267-54	DW8832RSU7	Unicryl (ud) 26mm curv.rev.cutt. 2(3/0) 75cm
109	10011-07	DW8832RSU7	Unicryl (ud) 26mm curv.rev.cutt. 2(3/0) 75cm
110	10062-15	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
111	10174-68	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

112	10267-74	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
113	10389-14	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
114	10389-14	DW8833RS7U	Unicryl (ud) 26mm curv.rev.cutt. 3(2/0) 75cm
115	10186-01	DW8881R	Unicryl 12mm curv.rev.cutt. 1.5(4/0) 75cm
116	10389-12	DW8800	Unicryl 16mm 1/2c.r.b. 1 (5/0) 75cm
117	10331-01	DW8861	Unicryl 16mm curv.cutt. 1.5(4/0) 75cm
118	10140-50	SD8114S	Unicryl 19mm 1/2c.cutt. 1.5(4/0) 75cm
119	10170-04	DW8902R	Unicryl 19mm curv.rev.cutt. 2(3/0) 70cm
120	10203-43	DW8902R	Unicryl 19mm curv.rev.cutt. 2(3/0) 70cm
121	10319-06	SE 663	Unicryl 20mm 1/2c.cutt. 1.5(4/0) 75cm
122	10112-11	SE 365	Unicryl 20mm 1/2c.cutt. 2(3/0) 75cm
123	10247-03	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
124	10062-17	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
125	10140-53	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
126	10174-74	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
127	10500-26	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
128	10518-02	DW8921	Unicryl 20mm 1/2c.r.b. 1.5(4/0) 75cm
129	10171-06	DW8922	Unicryl 20mm 1/2c.r.b. 2(3/0) 75cm
130	10140-52	DW8922	Unicryl 20mm 1/2c.r.b. 2(3/0) 75cm
131	10500-25	DW8922	Unicryl 20mm 1/2c.r.b. 2(3/0) 75cm
132	10769-03	DW8922	Unicryl 20mm 1/2c.r.b. 2(3/0) 75cm
133	10469-01	DW8922	Unicryl 20mm 1/2c.r.b. 2(3/0) 75cm
134	10062-19	DW8937R	Unicryl 24mm curv.rev.cutt. 2(3/0) 70cm
135	10174-76	DW8937R	Unicryl 24mm curv.rev.cutt. 2(3/0) 70cm
136	10137-03	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
137	10140-55	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
138	10139-04	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
139	10174-78	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
140	10267-53	DW8942	Unicryl 26mm 1/2c.r.b. 2(3/0) 75cm
141	10140-54	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
142	10174-79	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
143	10267-57	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
144	10190-02	DW8943	Unicryl 26mm 1/2c.r.b. 3(2/0) 75cm
145	10062-24	DW8943S9	Unicryl 26mm 1/2c.r.b. 3(2/0) 90cm
146	10319-04	DW8832	Unicryl 26mm curv.cutt. 2(3/0) 75cm
147	10040-07	DW8833	Unicryl 26mm curv.cutt. 3(2/0) 75cm
148	10363-02	DW8833	Unicryl 26mm curv.cutt. 3(2/0) 75cm
149	10319-05	DW8833R	Unicryl 26mm curv.rev.cutt. 3(2/0) 75cm
150	10171-05	SD8295	Unicryl 27mm 1/2c.r.b. 3(2/0) 75cm
151	10320-01	SE 206	Unicryl 30mm 1/2c.r.b. 3(2/0) 90cm
152	10140-59	DW8952	Unicryl 30mm 1/2c.r.b. 2(3/0) 75cm
153	10137-05	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

154	10666-05	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
155	10062-27	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
156	10061-19	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
157	10139-07	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
158	10140-58	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
159	10262-02	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
160	10500-30	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
161	10488-02	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
162	10091-01	DW8953	Unicryl 30mm 1/2c.r.b. 3(2/0) 75cm	
163	10112-01	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm	
164	10140-56	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm	
165	10140-56	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm	
166	10174-82	DW8954	Unicryl 30mm 1/2c.r.b. 3.5(0) 75cm	
167	10174-83	DW8955	Unicryl 30mm 1/2c.r.b. 4(1) 75cm	
168	10689-04	DW8993	Unicryl 35mm 1/2c.r.b. 3(2/0) 75cm	
169	10062-29	DW8993	Unicryl 35mm 1/2c.r.b. 3(2/0) 75cm	
170	10140-60	DW8993	Unicryl 35mm 1/2c.r.b. 3(2/0) 75cm	
171	10171-07	DW8995	Unicryl 35mm 1/2c.r.b. 4(1) 75cm	
172	10062-30	DW8563	Unicryl 35mm curv.cutt. 3(2/0) 75cm	
173	10205-03	SD8668	Unicryl 35mm curv.cutt. 3.5(0) 75cm	
174	10303-01	SE 130	Unicryl 37mm 1/2c.r.b. 3(2/0) 75cm	
175	10307-01	SE 130	Unicryl 37mm 1/2c.r.b. 3(2/0) 75cm	
176	10509-01	SE 130	Unicryl 37mm 1/2c.r.b. 3(2/0) 75cm	
177	10633-29	SE 130	Unicryl 37mm 1/2c.r.b. 3(2/0) 75cm	
178	10303-02	S 9892	Unicryl 37mm 1/2c.r.b. 3.5(0) 75cm	
179	10194-03	S 9892S	Unicryl 37mm 1/2c.r.b. 3.5(0) 90cm	
180	10247-02	DW8513	Unicryl 40mm 1/2c.r.b. 3(2/0) 75cm	
181	10140-67	DW8513	Unicryl 40mm 1/2c.r.b. 3(2/0) 75cm	
182	10172-17	DW8513	Unicryl 40mm 1/2c.r.b. 3(2/0) 75cm	
183	10564-08	DW8513	Unicryl 40mm 1/2c.r.b. 3(2/0) 75cm	
184	9755-58G	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
185	9755-58J	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
186	9755-58M	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
187	10370-62	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
188	10370-62C	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
189	10370-62B	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
190	10139-09	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
191	10140-65	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
192	10488-01	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
193	10190-01A	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
194	10203-15	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	
195	10633-02	DW8514	Unicryl 40mm 1/2c.r.b. 3.5(0) 75cm	



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

196	10194-02	DW8514S	Unicryl 40mm 1/2c.r.b. 3.5(0) 90cm	
197	10666-03	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm	
198	10062-33	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm	
199	10140-66	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm	
200	10172-18	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm	
201	10500-33	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm	
202	10762-01	DW8515	Unicryl 40mm 1/2c.r.b. 4(1) 75cm	
203	10338-03	DW8515S	Unicryl 40mm 1/2c.r.b. 4(1) 90cm	
204	10194-06	DW8515S	Unicryl 40mm 1/2c.r.b. 4(1) 90cm	
205	10172-20	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm	
206	10174-61	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm	
207	10267-61	DW8516	Unicryl 40mm 1/2c.r.b. 5(2) 75cm	
208	10140-68	DW8514H	Unicryl 40mm 1/2c.r.b.hvy. 3.5(0) 75cm	
209	10174-58	DW8514H	Unicryl 40mm 1/2c.r.b.hvy. 3.5(0) 75cm	
210	10267-62	DW8514H	Unicryl 40mm 1/2c.r.b.hvy. 3.5(0) 75cm	
211	10139-12	DW8515H	Unicryl 40mm 1/2c.r.b.hvy. 4(1) 75cm	
212	10172-19	DW8515H	Unicryl 40mm 1/2c.r.b.hvy. 4(1) 75cm	
213	10140-63	DW8524R	Unicryl 40mm 1/2c.rev.cutt. 3.5(0) 75cm	
214	10174-62	DW8524R	Unicryl 40mm 1/2c.rev.cutt. 3.5(0) 75cm	
215	10140-62	DW8525R	Unicryl 40mm 1/2c.rev.cutt. 4(1) 75cm	
216	10174-63	DW8525R	Unicryl 40mm 1/2c.rev.cutt. 4(1) 75cm	
217	10101-01	DW8525R	Unicryl 40mm 1/2c.rev.cutt. 4(1) 75cm	
218	10140-64	DW8526R	Unicryl 40mm 1/2c.rev.cutt. 5(2) 75cm	
219	10689-03	DT8415H	Unicryl 40mm 1/2c.taper.cut.hvy. 3(2/0) 75cm	
220	10205-01A	SD8189	Unicryl 40mm curv.cutt. 4(1) 75cm	
221	10363-01	SD8189	Unicryl 40mm curv.cutt. 4(1) 75cm	
222	10633-39	DW8855S7	Unicryl 45mm 1/2c.r.b. 4(1) 75cm	
223	10452-04	DW8856S7	Unicryl 45mm 1/2c.r.b. 5(2) 75cm	
224	10399-01	DW8854	Unicryl 45mm 1/2c.r.b.hvy. 3.5(0) 75cm	
225	10452-01	DW8854	Unicryl 45mm 1/2c.r.b.hvy. 3.5(0) 75cm	
226	10376-02	DW3855	Unicryl 45mm 1/2c.r.b.hvy. 4(1) 75cm	
227	10689-01	DW8855	Unicryl 45mm 1/2c.r.b.hvy. 4(1) 75cm	
228	10502-02	DW8856	Unicryl 45mm 1/2c.r.b.hvy.5(2) 75cm	
229	10399-03	DW8856	Unicryl 45mm 1/2c.r.b.hvy.5(2) 75cm	
230	10267-64	DW8856	Unicryl 45mm 1/2c.r.b.hvy.5(2) 75cm	
231	10102-01	DW8856	Unicryl 45mm 1/2c.r.b.hvy.5(2) 75cm	
232	10098-01	DT8536	Unicryl 45mm 1/2c taper.cut 5(2) 75cm	
233	10062-37	DW8556	Unicryl 50mm 1/2c.r.b. 5(2) 75cm	
234	10423-27	DW8055R	Unicryl 50mm 1/2c.rev.cutt. 4(1) 75cm	
235	10633-38	DW8572	Unicryl 60mm str.cutt. 2(3/0) 75cm	
236	10138-01	DW8628DS	Unicryl 6mm curv.spat.d/a.0.5(7/0) 45cm	
237	10423-29	DW8379D	Unicryl 8mm 1/4c.spat.d/a 0.7(6/0) 45cm	



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

238	10712-05	DW8502	Unicryl cut lengths 2(3/0) 1.5m
239	10308-06	DW8504	Unicryl cut lengths 3.5(0) 1.5m
240	10308-07	DW8505	Unicryl cut lengths 4(1) 1.5m
241	10203-47	DW8505	Unicryl cut lengths 4(1) 1.5m
242	10338-01	DW8503	Unicryl cut lengths 3(2/0) 1.5 metres
243	10423-33	DW8503	Unicryl cut lengths 3(2/0) 1.5 metres

Customer Action

Please review your inventory for products affected by this issue and perform the following actions:

1. Immediately identify and quarantine all unused, affected products in your inventory.
2. Return all unused, affected products in your inventory to Unimed. Our Representatives will assist you in the return of affected product as necessary.

Please share this notification with others in your organization as appropriate. If product within scope of this recall has been forwarded to another facility, please notify the facility of the issue and assist with the inconvenience that it may cause

Please be assured that patient safety and products quality remain our primary concern.

If you have any questions, please contact your Unimed Representatives.

Sincerely,

Nader Hamdan Assiri





United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

Urgent Field Safety Notice
Unicaprone Surgical Suture Recall

Date: 13/11/2016

Unimed Reference: FSN02

Dear NCMDR Team

This letter to inform you that we have identified that the Suture Breaking Strength in few incidents are non compliant to European Pharmacopeia requirements

We are initiating a voluntary recall for the following batches:

Serial	Batch#	Item Code	Description
1	10370-112C	DW13029V	Unicaprone (violet) 22mm 1/2c.r.b. 2(3/0) 70cm
2	10370-111A	DW13058V	Unicaprone (violet) 30mm 1/2c.r.b. 3(2/0) 70cm
3	10370-32A	DW13503V	Unicaprone (violet) standard lengths 3(2/0) 1.5 metres
4	10370-32C	DW13503V	Unicaprone (violet) standard lengths 3(2/0) 1.5 metres
5	10850-44	DW13013	Unicaprone 17mm curv.r.b. 1.5(4/0) 70cm
6	10860-02	DP13024	Unicaprone 19mm curv.P-type.rev.cutt. 2(3/0) 45cm
7	10846-13	DW13901RU	Unicaprone 19mm curv.rev.cutt. 1.5(4/0) 45cm
8	10370-107C	DW13902RU	Unicaprone 19mm curv.rev.cutt. 2(3/0) 45cm
9	10370-107D	DW13902RU	Unicaprone 19mm curv.rev.cutt. 2(3/0) 45cm
10	10846-14	DW13902RU	Unicaprone 19mm curv.rev.cutt. 2(3/0) 45cm
11	10801-07	DW13902RS	Unicaprone 19mm curv.rev.cutt. 2(3/0) 75cm
12	10850-45	DW13029	Unicaprone 22mm 1/2c.r.b. 2(3/0) 70cm
13	10370-106B	DW13833RS	Unicaprone 26mm curv.rev.cutt. 3(2/0) 70cm
14	10370-106C	DW13833RS	Unicaprone 26mm curv.rev.cutt. 3(2/0) 70cm
15	10370-106D	DW13833RS	Unicaprone 26mm curv.rev.cutt. 3(2/0) 70cm
16	10370-106F	DW13833RS	Unicaprone 26mm curv.rev.cutt. 3(2/0) 70cm
17	10712-03	DW13572U	Unicaprone 60mm str.cutt. 2(3/0) 70cm
18	10370-33	DW13573	Unicaprone 60mm str.cutt. 3(2/0) 70cm
19	9755-113C	DW13801VS	Unicaprone (violet) 16mm 1/2c.r.b. 1.5(4/0) 75cm
20	9755-112B	DW13029V	Unicaprone (violet) 22mm 1/2c.r.b. 2(3/0) 70cm
21	9755-110B	DW13058V	Unicaprone (violet) 30mm 1/2c.r.b. 3(2/0) 70cm
22	9755-109C	DW13084V	Unicaprone (violet) 36mm 1/2c.r.b. 3.5(0) 70cm
23	9755-29C	DW13503V	Unicaprone (violet) standard lengths 3(2/0) 1.5 metres
24	9755-108F	DW13760RU	Unicaprone 13mm curv.rev.cutt.1(5/0) 45cm
25	10370-109B	DW13760RU	Unicaprone 13mm curv.rev.cutt.1(5/0) 45cm
26	10339-03	DP13008	Unicaprone 16mm curv.P-type.rev.cutt. 1.5(4/0) 45cm
27	10423-05	DW13861RSU	Unicaprone 16mm curv.rev.cutt. 1.5(4/0) 45cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

28	10267-28	DW13861RSU	Unicaprone 16mm curv.rev.cutt. 1.5(4/0) 45cm
29	10240-15	DW13861RS7	Unicaprone 16mm curv.rev.cutt. 1.5(4/0) 75cm
30	10283-01	DW13019	Unicaprone 17mm 1/2c.r.b. 2(3/0) 70cm
31	10290-01	DP13023	Unicaprone 19mm curv.P-type.rev.cutt. 1.5(4/0) 45cm
32	10468-01	DP13023	Unicaprone 19mm curv.P-type.rev.cutt. 1.5(4/0) 45cm
33	9755-107B	DW13901RU	Unicaprone 19mm curv.rev.cutt. 1.5(4/0) 45cm
34	10062-02	DW13901RU	Unicaprone 19mm curv.rev.cutt. 1.5(4/0) 45cm
35	10174-13	DW13901RU	Unicaprone 19mm curv.rev.cutt. 1.5(4/0) 45cm
36	9755-106B	DW13902RU	Unicaprone 19mm curv.rev.cutt. 2(3/0) 45cm
37	9755-106D	DW13902RU	Unicaprone 19mm curv.rev.cutt. 2(3/0) 45cm
38	10249-03	DW13902RS	Unicaprone 19mm curv.rev.cutt. 2(3/0) 75cm
39	10028-32	DW13831S	Unicaprone 26mm curv.cutt. 1.5(4/0) 75cm
40	10249-01	DW13832RS	Unicaprone 26mm curv.rev.cutt. 2(3/0) 75cm
41	10370-106A	DW13833RS	Unicaprone 26mm curv.rev.cutt. 3(2/0) 70cm
42	10608-03	DW13833R	Unicaprone 26mm curv.rev.cutt. 3(2/0) 75cm
43	10377-01	DW13058	Unicaprone 30mm 1/2c.r.b. 3(2/0) 70cm
44	10589-01	DT13079	Unicaprone 35mm 1/2c.taper.cut. 3.5(0) 70cm
45	10062-04	DW13514S	Unicaprone 40mm 1/2c.r.b. 3.5(0) 75cm
46	10334-01	DW13115	Unicaprone 40mm 1/2c.r.b.hvy.4(1) 70cm
47	10370-33A	DW13573	Unicaprone 60mm str.cutt. 3(2/0) 70cm

Customer Action

Please review your inventory for products affected by this issue and perform the following actions:

1. Immediately identify and quarantine all unused, affected products in your inventory.
2. Return all unused, affected products in your inventory to Unimed. Our Representatives will assist you in the return of affected product as necessary.

Please share this notification with others in your organization as appropriate. If product within scope of this recall has been forwarded to another facility, please notify the facility of the issue and assist with the inconvenience that it may cause

Please be assured that patient safety and products quality remain our primary concern.

If you have any questions, please contact your Unimed Representatives.

Sincerely,

Nader Hamdan Assiri





United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

Urgent Field Safety Notice
Unicryl Quick Surgical Suture Recall

Date: 13/11/2016

Unimed Reference: FSN03

Dear NCMDR Team

This letter to inform you that we have identified that the Suture Breaking Strength in few incidents are non compliant to European Pharmacopeia requirements

We are initiating a voluntary recall for the following batches:

Product Name : UniCryl Quick

Serial	Batch#	Item Code	Description
1	10936-01	DP12219R	Unicryl quick 12mm curv.P-type.rev.cutt. 0.7(6/0) 45cm
2	10370-102B	DW12220R	Unicryl quick 12mm curv.rev.cutt. 1(5/0) 45cm
3	10370-101A	DW12801	Unicryl quick 16mm 1/2c.r.b. 1.5(4/0) 75cm
4	10850-42	DW12860	Unicryl quick 16mm curv.cutt. 1(5/0) 75cm
5	10706-02	DW12861	Unicryl quick 16mm curv.cutt. 1.5(4/0) 75cm
6	11023-07	DW12902	Unicryl quick 19mm curv.cutt. 2(3/0) 75cm
7	10943-20	DW12902R	Unicryl quick 19mm curv.rev.cutt. 2(3/0) 75cm
8	10370-100C	DW12922	Unicryl quick 20mm 1/2c.r.b. 2(3/0) 75cm
9	10370-100B	DW12922	Unicryl quick 20mm 1/2c.r.b. 2(3/0) 75cm
10	10370-100D	DW12922	Unicryl quick 20mm 1/2c.r.b. 2(3/0) 75cm
11	10854-01	SD12000	Unicryl quick 22mm 1/2c.r.b. 2(3/0) 75cm
12	10846-23	DW12833	Unicryl quick 26mm curv.cutt. 3(2/0) 75cm
13	11025-02	DW12833	Unicryl quick 26mm curv.cutt. 3(2/0) 75cm
14	10370-98C	DW12833RU	Unicryl quick 26mm curv.rev.cutt. 3(2/0) 75cm
15	10370-98D	DW12833RU	Unicryl quick 26mm curv.rev.cutt. 3(2/0) 75cm
16	10370-99C	DW12953	Unicryl quick 30mm 1/2c.r.b. 3(2/0) 75cm
17	10825-03	DW12953	Unicryl quick 30mm 1/2c.r.b. 3(2/0) 75cm
18	10929-01	DT12873US9	Unicryl quick 35mm 1/2c.taper.cut. 3(2/0) 90cm
19	10881-01	DT12874S9	Unicryl quick 35mm 1/2c.taper.cut. 3.5(0) 90cm
20	10370-95C	SD12030HS	Unicryl quick 36mm 1/2c.r.b.hvy. 3(2/0) 90cm
21	11101-106	SD12030HS	Unicryl quick 36mm 1/2c.r.b.hvy. 3(2/0) 90cm
22	10774-01	DW12513	Unicryl quick 40mm 1/2c.r.b. 3(2/0) 75cm
23	10370-94D	SD12008H	Unicryl quick 50mm 1/2.c.r.b.hvy. 3(2/0) 75cm
24	10370-94C	SD12008H	Unicryl quick 50mm 1/2.c.r.b.hvy. 3(2/0) 75cm
25	11101-105	SD12008H	Unicryl quick 50mm 1/2.c.r.b.hvy. 3(2/0) 75cm
26	9755-103A	DW12801	Unicryl quick 16mm 1/2c.r.b. 1.5(4/0) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

27	10267-39	DW12861	Unicryl quick 16mm curv.cutt. 1.5(4/0) 75cm
28	10200-01	DW12860R	Unicryl quick 16mm curv.rev.cutt. 1(5/0) 75cm
29	10183-03	DW12901R	Unicryl quick 19mm curv.rev.cutt. 1.5(4/0) 75cm
30	9978-03	DW12902R	Unicryl quick 19mm curv.rev.cutt. 2(3/0) 75cm
31	10140-40	DW12902R	Unicryl quick 19mm curv.rev.cutt. 2(3/0) 75cm
32	10139-17	DW12902R	Unicryl quick 19mm curv.rev.cutt. 2(3/0) 75cm
33	10172-04	DW12902R	Unicryl quick 19mm curv.rev.cutt. 2(3/0) 75cm
34	9755-102B	DW12922	Unicryl quick 20mm 1/2c.r.b. 2(3/0) 75cm
35	10191-04	DW12922	Unicryl quick 20mm 1/2c.r.b. 2(3/0) 75cm
36	10062-45	DW12832S	Unicryl quick 26mm curv.cutt. 2(3/0) 45cm
37	9979-33	DW12833	Unicryl quick 26mm curv.cutt. 3(2/0) 75cm
38	10062-46	DW12833	Unicryl quick 26mm curv.cutt. 3(2/0) 75cm
39	10139-03	DW12833	Unicryl quick 26mm curv.cutt. 3(2/0) 75cm
40	10062-47	DP12832RS7	Unicryl quick 26mm curv.P-type.rev.cutt. 2(3/0) 75cm
41	10174-03	DP12832RS7	Unicryl quick 26mm curv.P-type.rev.cutt. 2(3/0) 75cm
42	10612-01	DP12832RS7	Unicryl quick 26mm curv.P-type.rev.cutt. 2(3/0) 75cm
43	10502-03	DW12832RS	Unicryl quick 26mm curv.rev.cutt. 2(3/0) 45cm
44	9755-100B	DW12833RU	Unicryl quick 26mm curv.rev.cutt. 3(2/0) 75cm
45	10370-98	DW12833RU	Unicryl quick 26mm curv.rev.cutt. 3(2/0) 75cm
46	9755-101A	DW12953	Unicryl quick 30mm 1/2c.r.b. 3(2/0) 75cm
47	10721-01	DT12873US9	Unicryl quick 35mm 1/2c.taper.cut. 3(2/0) 90cm
48	10254-01	DT12874S9	Unicryl quick 35mm 1/2c.taper.cut. 3.5(0) 90cm
49	10370-96	SD12062H	Unicryl quick 36mm 1/2c.r.b.hvy 2(3/0) 90cm
50	10370-96A	SD12062H	Unicryl quick 36mm 1/2c.r.b.hvy 2(3/0) 90cm
51	9755-97A	SD12030HS	Unicryl quick 36mm 1/2c.r.b.hvy. 3(2/0) 90cm
52	10370-95	SD12030HS	Unicryl quick 36mm 1/2c.r.b.hvy. 3(2/0) 90cm
53	10370-95A	SD12030HS	Unicryl quick 36mm 1/2c.r.b.hvy. 3(2/0) 90cm
54	10370-95B	SD12030HS	Unicryl quick 36mm 1/2c.r.b.hvy. 3(2/0) 90cm
55	10700-27	ST12059	Unicryl quick 36mm 1/2c.taper.cut. 3(2/0) 90cm
56	10370-94	SD12008H	Unicryl quick 50mm 1/2.c.r.b.hvy. 3(2/0) 75cm
57	10608-01	DW12573	Unicryl quick 60mm str.cutt. 3(2/0) 75cm



United Medical Industries Co. Ltd.
Customer Action

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

Please review your inventory for products affected by this issue and perform the following actions:

1. Immediately identify and quarantine all unused, affected products in your inventory.
2. Return all unused, affected products in your inventory to Unimed. Our Representatives will assist you in the return of affected product as necessary.

Please share this notification with others in your organization as appropriate. If product within scope of this recall has been forwarded to another facility, please notify the facility of the issue and assist with the inconvenience that it may cause

Please be assured that patient safety and products quality remain our primary concern.
If you have any questions, please contact your Unimed Representatives.

Sincerely,
Nader Hamdan Assiri





United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

Urgent Field Safety Notice
Uniglactin Surgical Suture Recall

Date:13/11/2016

Unimed Reference: FSN04

Dear NCMDR Team

This letter to inform you that we have identified that the Suture Breaking Strength in few incidents are non compliant to European Pharmacopeia requirements

We are initiating a voluntary recall for the following batches:

Product Name : Uniglactin

serial	Batch#	Item Code	Description
1	10370-60C	DW11220RSU	Uniglactin (ud) 12mm curv.rev.cutt. 1(5/0) 75cm
2	10760-01A	DW11860U	Uniglactin (ud) 16mm curv.cutt. 1(5/0) 75cm
3	11001-01A	DW11832RS7U	Uniglactin (ud) 26mm curv.rev.cutt. 2(3/0) 75cm
4	10370-89L	DW11219S4	Uniglactin 12mm curv.cutt. 0.7(6/0) 45cm
5	10370-89J	DW11219S4	Uniglactin 12mm curv.cutt. 0.7(6/0) 45cm
6	10370-89N	DW11219S4	Uniglactin 12mm curv.cutt. 0.7(6/0) 45cm
7	10370-89R	DW11219S4	Uniglactin 12mm curv.cutt. 0.7(6/0) 45cm
8	10370-83D	DW11388S4	Uniglactin 13mm curv.r.b. 1.5(4/0) 45cm
9	10370-83E	DW11388S4	Uniglactin 13mm curv.r.b. 1.5(4/0) 45cm
10	10370-83I	DW11388S4	Uniglactin 13mm curv.r.b. 1.5(4/0) 45cm
11	10370-81H	DW11389	Uniglactin 13mm curv.r.b. 2(3/0) 45cm
12	10370-81F	DW11389	Uniglactin 13mm curv.r.b. 2(3/0) 45cm
13	10775-01	DW11398S	Uniglactin 13mm str.cutt. 1.5(4/0) 75cm
14	10370-82F	DW11862	Uniglactin 16mm curv.cutt. 2(3/0) 75cm
15	10370-82K	DW11862	Uniglactin 16mm curv.cutt. 2(3/0) 75cm
16	10370-82J	DW11862	Uniglactin 16mm curv.cutt. 2(3/0) 75cm
17	10370-86B	DW11860	Uniglactin 16mm curv.cutt. 1(5/0) 75cm
18	10370-86D	DW11860	Uniglactin 16mm curv.cutt. 1(5/0) 75cm
19	10370-86E	DW11860	Uniglactin 16mm curv.cutt. 1(5/0) 75cm
20	11101-96	DW11860	Uniglactin 16mm curv.cutt. 1(5/0) 75cm
21	10370-84B	DW11841S	Uniglactin 16mm curv.r.b. 1.5(4/0) 75cm
22	10370-80B	DW11922	Uniglactin 20mm 1/2c.r.b. 2(3/0) 75cm
23	10370-80A	DW11922	Uniglactin 20mm 1/2c.r.b. 2(3/0) 75cm
24	10825-04	DW11922	Uniglactin 20mm 1/2c.r.b. 2(3/0) 75cm
25	11101-89	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm
26	11101-89G	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm
27	11101-89F	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

28	11101-83	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
29	11101-83I	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
30	11101-83D	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
31	11101-79	DT11823	Uniglactin 25mm 1/2c.taper.cut. 3(2/0) 75cm
32	10370-75B	DT11823H	Uniglactin 25mm 1/2c.taper.cut.hvy. 3(2/0) 75cm
33	10370-75C	DT11823H	Uniglactin 25mm 1/2c.taper.cut.hvy. 3(2/0) 75cm
34	10370-65C	SD11170H	Uniglactin 25mm 5/8c.r.b.hvy. 3.5(0) 75cm
35	10370-65E	SD11170H	Uniglactin 25mm 5/8c.r.b.hvy. 3.5(0) 75cm
36	10700-05	DW11941	Uniglactin 26mm 1/2c.r.b. 1.5(4/0) 75cm
37	10949-01	DW11943	Uniglactin 26mm 1/2c.r.b. 3(2/0) 75cm
38	11107-01	DW11943	Uniglactin 26mm 1/2c.r.b. 3(2/0) 75cm
39	10700-02A	DW11945	Uniglactin 26mm 1/2c.r.b. 4(1) 75cm
40	11107-04	SD11130	Uniglactin 26mm 1/2c.r.b.1(5/0) 75cm
41	10370-74A	DW11943SA	Uniglactin 26mm(Black) 1/2c.r.b. 3(2/0) 70cm
42	10961-02	DW11943SA	Uniglactin 26mm(Black) 1/2c.r.b. 3(2/0) 70cm
43	11101-84B	DW11953	Uniglactin 30mm 1/2c.r.b. 3(2/0) 75cm
44	11101-84	DW11953	Uniglactin 30mm 1/2c.r.b. 3(2/0) 75cm
45	10700-07	DW11953	Uniglactin 30mm 1/2c.r.b. 3(2/0) 75cm
46	10370-66I	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
47	10370-66H	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
48	10370-66J	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
49	10370-66N	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
50	10370-66M	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
51	10370-66O	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
52	11101-66	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
53	11101-66E	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
54	10370-69H	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
55	10370-69G	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
56	10370-69J	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
57	10370-69L	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
58	10370-69M	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
59	10370-69N	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
60	10370-69O	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
61	10370-69P	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
62	11101-70A	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm
63	10370-77A	DW11953D	Uniglactin 30mm 1/2c.r.b.d/a 3(2/0) 75cm
64	11101-81A	DW11953D	Uniglactin 30mm 1/2c.r.b.d/a 3(2/0) 75cm
65	10370-76E	SD11085A	Uniglactin 30mm curv.r.b. 3(2/0) 75cm
66	10370-76F	SD11085A	Uniglactin 30mm curv.r.b. 3(2/0) 75cm
67	11101-65A	SD11171	Uniglactin 30mm curv.r.b.3.5(0) 75cm
68	11101-65	SD11171	Uniglactin 30mm curv.r.b.3.5(0) 75cm
69	11094-01	DW11993	Uniglactin 35mm 1/2c.r.b. 3(2/0) 75cm



®

United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

70	11101-71	DT11565	Uniglactin 35mm 1/2c.taper.cut.hvy. 4(1) 75cm
71	11101-71D	DT11565	Uniglactin 35mm 1/2c.taper.cut.hvy. 4(1) 75cm
72	10370-61B	SD11156S	Uniglactin 36mm j.r.b. 3.5(0) 70cm
73	10370-61C	SD11156S	Uniglactin 36mm j.r.b. 3.5(0) 70cm
74	10370-67V	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
75	10370-67Y	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
76	10370-67AH	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
77	10370-67AI	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
78	10370-67AF	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
79	10370-67AE	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
80	10370-67AJ	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
81	10370-67AD	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
82	10370-67AN	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
83	10370-67AK	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
84	11101-67A	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
85	11101-67J	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm
86	11101-72M	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm
87	11101-72A	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm
88	10834-01	DW11515S	Uniglactin 40mm 1/2c.r.b. 4(1) 90cm
89	11177-01A	DW11515S	Uniglactin 40mm 1/2c.r.b. 4(1) 90cm
90	10700-15	DW11516	Uniglactin 40mm 1/2c.r.b. 5(2) 75cm
91	10370-78D	SD11175	Uniglactin 40mm curv.r.b. 3(2/0) 75cm
92	10370-78H	SD11175	Uniglactin 40mm curv.r.b. 3(2/0) 75cm
93	10370-78F	SD11175	Uniglactin 40mm curv.r.b. 3(2/0) 75cm
94	11101-68	DW11855S	Uniglactin 45mm 1/2c.r.b. 4(1) 75cm
95	11101-68B	DW11855S	Uniglactin 45mm 1/2c.r.b. 4(1) 75cm
96	10913-01	DW11855BS	Uniglactin 45mm 1/2c.r.b.blunt. 4(1) 90cm
97	10370-71N	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
98	10370-71P	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
99	10370-71Q	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
100	10370-71V	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
101	10370-71U	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
102	10370-71AA	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
103	10370-71Y	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
104	10370-71AB	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
105	10370-71X	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
106	10370-71AD	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
107	10370-71AC	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
108	10370-71Z	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
109	10370-71AI	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
110	10370-71AG	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
111	10370-71AH	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

112	10370-71AE	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
113	11101-74E	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
114	11101-74A	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
115	11101-74C	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
116	11101-74B	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm
117	10986-01	DT11536H	Uniglactin 45mm 1/2c.taper.cut.hvy. 5(2) 75cm
118	10700-18	SD11155	Uniglactin 48mm 1/2c.r.b. 5(2) 75cm
119	10370-63I	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
120	10370-63H	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
121	10370-63J	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
122	10370-63K	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
123	11101-62A	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
124	10370-59C	DW11555	Uniglactin 50mm 1/2c.r.b. 4(1) 75cm
125	10370-59D	DW11555	Uniglactin 50mm 1/2c.r.b. 4(1) 75cm
126	10370-59F	DW11555	Uniglactin 50mm 1/2c.r.b. 4(1) 75cm
127	10370-59H	DW11555	Uniglactin 50mm 1/2c.r.b. 4(1) 75cm
128	11101-58B	DW11555	Uniglactin 50mm 1/2c.r.b. 4(1) 75cm
129	10370-68B	SD11172BH	Uniglactin 50mm 1/2c.r.b.blnt.hvy.4(1) 75cm
130	10370-68D	SD11172BH	Uniglactin 50mm 1/2c.r.b.blnt.hvy.4(1) 75cm
131	10370-79I	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
132	10370-79H	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
133	10370-79L	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
134	10370-79O	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
135	10370-79Q	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
136	10370-79P	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
137	10370-79R	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
138	11101-85	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
139	10370-72B	SD11174	Uniglactin 65mm curv.r.b. 3(2/0) 75cm
140	10370-64C	SD11169	Uniglactin 65mm curv.r.b. 3.5(0) 75cm
141	10370-64D	SD11169	Uniglactin 65mm curv.r.b. 3.5(0) 75cm
142	11101-63	SD11169	Uniglactin 65mm curv.r.b. 3.5(0) 75cm
143	10370-93F	SD11053S	Uniglactin 8mm 1/2c.rev.cutt. d/a 0.4(8/0) 45cm
144	10370-91F	SD11121RD	Uniglactin 8mm 1/2c.rev.cutt.d/a 0.5(7/0) 45cm
145	10370-92F	DW11667S	Uniglactin 8mm curv.r.b. 0.4(8/0) 30cm
146	10370-105A	DW11462	Uniglactin cut lengths 2(3/0) 10x50cm
147	10370-105B	DW11462	Uniglactin cut lengths 2(3/0) 10x50cm
148	11101-116	DW11462	Uniglactin cut lengths 2(3/0) 10x50cm
149	10370-104B	DW11463	Uniglactin cut lengths 3(2/0) 10x50cm
150	11101-115A	DW11463	Uniglactin cut lengths 3(2/0) 10x50cm
151	10370-73C	DW11503	Uniglactin cut lengths 3(2/0) 1.5 metres
152	10370-73E	DW11503	Uniglactin cut lengths 3(2/0) 1.5 metres
153	11101-76	DW11503	Uniglactin cut lengths 3(2/0) 1.5 metres



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

154	10918-01	DW11503	Uniglactin cut lengths 3(2/0) 1.5 metres
155	10370-103B	DW11464	Uniglactin cut lengths 3.5(0) 10x50cm
156	10760-01	DW11860U	Uniglactin (ud) 16mm curv.cutt. 1(5/0) 75cm
157	10358-01	DP11902SU	Uniglactin (ud) 19mm curv.P-type.cutt. 2(3/0) 45cm
158	10578-01	SD11091S4UE	Uniglactin (ud) 20mm curv.extra.cutt. 2(3/0) 45cm
159	9755-93N	SD11001S4	Uniglactin 13mm 1/2c.r.b. 0.7(6/0) 45cm
160	9755-93Q	SD11001S4	Uniglactin 13mm 1/2c.r.b. 0.7(6/0) 45cm
161	9755-89F	DW11387S4	Uniglactin 13mm curv.r.b. 1(5/0) 45cm
162	9755-89E	DW11387S4	Uniglactin 13mm curv.r.b. 1(5/0) 45cm
163	9755-89G	DW11387S4	Uniglactin 13mm curv.r.b. 1(5/0) 45cm
164	9755-89J	DW11387S4	Uniglactin 13mm curv.r.b. 1(5/0) 45cm
165	9755-85E	DW11388S4	Uniglactin 13mm curv.r.b. 1.5(4/0) 45cm
166	9755-85H	DW11388S4	Uniglactin 13mm curv.r.b. 1.5(4/0) 45cm
167	9755-82H	DW11389	Uniglactin 13mm curv.r.b. 2(3/0) 45cm
168	9755-82G	DW11389	Uniglactin 13mm curv.r.b. 2(3/0) 45cm
169	10555-01	DW11800	Uniglactin 16mm 1/2c.r.b. 1 (5/0) 75cm
170	9755-83F	DW11862	Uniglactin 16mm curv.cutt. 2(3/0) 75cm
171	9755-83L	DW11862	Uniglactin 16mm curv.cutt. 2(3/0) 75cm
172	9755-90G	DW11860	Uniglactin 16mm curv.cutt. 1(5/0) 75cm
173	9755-87J	DW11861	Uniglactin 16mm curv.cutt. 1.5(4/0) 75cm
174	9755-87i	DW11861	Uniglactin 16mm curv.cutt. 1.5(4/0) 75cm
175	9755-87L	DW11861	Uniglactin 16mm curv.cutt. 1.5(4/0) 75cm
176	9755-84K	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm
177	9755-84G	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm
178	9755-84N	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm
179	9755-84O	SD11008	Uniglactin 25mm 1/2c.r.b. 2(3/0) 75cm
180	9755-78Z	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
181	9755-78O	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
182	9755-78X	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
183	9755-78Y	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
184	9755-78Q	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
185	9755-78W	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
186	10431-01	SD11021	Uniglactin 25mm 1/2c.r.b. 3(2/0) 75cm
187	9755-74E	DT11823H	Uniglactin 25mm 1/2c.taper.cut.hvy. 3(2/0) 75cm
188	9755-74D	DT11823H	Uniglactin 25mm 1/2c.taper.cut.hvy. 3(2/0) 75cm
189	10370-75	DT11823H	Uniglactin 25mm 1/2c.taper.cut.hvy. 3(2/0) 75cm
190	9755-79F	DW11953	Uniglactin 30mm 1/2c.r.b. 3(2/0) 75cm
191	10390-04	DW11953	Uniglactin 30mm 1/2c.r.b. 3(2/0) 75cm
192	9755-63K	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
193	9755-63J	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
194	10370-66A	DW11954	Uniglactin 30mm 1/2c.r.b. 3.5(0) 75cm
195	9755-67L	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

196	9755-67K	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm	
197	9755-67M	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm	
198	9755-67O	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm	
199	10370-69C	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm	
200	10370-69D	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm	
201	10370-69A	DW11955	Uniglactin 30mm 1/2c.r.b. 4(1) 75cm	
202	9755-76D	DW11953D	Uniglactin 30mm 1/2c.r.b.d/a 3(2/0) 75cm	
203	9755-75B	SD11085A	Uniglactin 30mm curv.r.b. 3(2/0) 75cm	
204	9755-75E	SD11085A	Uniglactin 30mm curv.r.b. 3(2/0) 75cm	
205	10370-76A	SD11085A	Uniglactin 30mm curv.r.b. 3(2/0) 75cm	
206	9755-62B	SD11171	Uniglactin 30mm curv.r.b.3.5(0) 75cm	
207	9755-62C	SD11171	Uniglactin 30mm curv.r.b.3.5(0) 75cm	
208	9755-68J	DT11565	Uniglactin 35mm 1/2c.taper.cut.hvy. 4(1) 75cm	
209	9755-68K	DT11565	Uniglactin 35mm 1/2c.taper.cut.hvy. 4(1) 75cm	
210	9755-68M	DT11565	Uniglactin 35mm 1/2c.taper.cut.hvy. 4(1) 75cm	
211	10370-61A	SD11156S	Uniglactin 36mm j.r.b. 3.5(0) 70cm	
212	9755-64M	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
213	9755-64T	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
214	9755-64R	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
215	9755-64AG	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
216	9755-64AA	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
217	9755-64Z	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
218	10370-67C	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
219	10370-67G	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
220	10370-67F	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
221	10370-67H	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
222	10370-67D	DW11514	Uniglactin 40mm 1/2c.r.b. 3.5(0) 75cm	
223	9755-69S	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
224	9755-69W	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
225	9755-69V	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
226	9755-69T	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
227	9755-69X	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
228	9755-69Q	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
229	9755-69U	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
230	9755-69AG	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
231	10305-02	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
232	10390-01B	DW11515	Uniglactin 40mm 1/2c.r.b. 4(1) 75cm	
233	10357-06	DW11515S	Uniglactin 40mm 1/2c.r.b. 4(1) 90cm	
234	9755-77F	SD11175	Uniglactin 40mm curv.r.b. 3(2/0) 75cm	
235	9755-65H	DW11855S	Uniglactin 45mm 1/2c.r.b. 4(1) 75cm	
236	9755-65L	DW11855S	Uniglactin 45mm 1/2c.r.b. 4(1) 75cm	
237	9755-71i	DT11536	Uniglactin 45mm 1/2c.taper.cut. 5(2) 75cm	



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

238	9755-71H	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
239	9755-71M	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
240	9755-71N	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
241	9755-71U	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
242	9755-71S	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
243	9755-71K	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
244	9755-71W	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
245	9755-71V	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
246	10370-71A	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
247	10370-71i	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
248	10370-71H	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
249	10370-71G	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
250	10370-71F	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
251	10731-02	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
252	10305-01	DT11536	Uniglactin 45mm 1/2c.taper.cutt. 5(2) 75cm
253	10700-17	SD11075	Uniglactin 48mm 1/2c.r.b. 4(1) 75cm
254	9755-59G	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
255	9755-59i	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
256	10370-63A	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
257	10370-63	DW11554	Uniglactin 50mm 1/2c.r.b. 3.5(0) 75cm
258	10370-59	DW11555	Uniglactin 50mm 1/2c.r.b. 4(1) 75cm
259	10633-34	DW11556S9	Uniglactin 50mm 1/2c.r.b. 5(2) 90cm
260	9755-80i	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
261	9755-80J	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
262	9755-80K	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
263	9755-80U	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
264	9755-80S	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
265	10370-79	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
266	10370-79E	DW11573	Uniglactin 60mm str.cutt. 3(2/0) 75cm
267	10370-72A	SD11174	Uniglactin 65mm curv.r.b. 3(2/0) 75cm
268	10506-02	DW11627D	Uniglactin 6mm curv.spat.d/a 0.4(8/0) 30cm
269	9755-96P	SD11053S	Uniglactin 8mm 1/2c.rev.cutt. d/a 0.4(8/0) 45cm
270	9755-94H	SD11121RD	Uniglactin 8mm 1/2c.rev.cutt.d/a 0.5(7/0) 45cm
271	9755-88G	DW11380D	Uniglactin 8mm 1/4c.spat.d/a 1(5/0) 45cm
272	9755-88N	DW11380D	Uniglactin 8mm 1/4c.spat.d/a 1(5/0) 45cm
273	10375-02	DW11379DS	Uniglactin 8mm 1/4c.spat.d/a. 0.7(6/0) 20cm
274	9755-95E	DW11667S	Uniglactin 8mm curv.r.b. 0.4(8/0) 30cm
275	9755-95F	DW11667S	Uniglactin 8mm curv.r.b. 0.4(8/0) 30cm
276	9755-95L	DW11667S	Uniglactin 8mm curv.r.b. 0.4(8/0) 30cm
277	9755-91E	DW11369	Uniglactin 8mm curv.r.b. 0.7(6/0) 30cm
278	10370-73	DW11503	Uniglactin cut lengths 3(2/0) 1.5 metres
279	10370-103	DW11464	Uniglactin cut lengths 3.5(0) 10x50cm



United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

Customer Action

Please review your inventory for products affected by this issue and perform the following actions:

1. Immediately identify and quarantine all unused, affected products in your inventory.
2. Return all unused, affected products in your inventory to Unimed. Our Representatives will assist you in the return of affected product as necessary.

Please share this notification with others in your organization as appropriate. If product within scope of this recall has been forwarded to another facility, please notify the facility of the issue and assist with the inconvenience that it may cause

Please be assured that patient safety and products quality remain our primary concern.

If you have any questions, please contact your Unimed Representatives.

Sincerely,

Nader Hamdan Assir:





United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 11 265 0665
Fax: +966 11 265 0569
www.unimed.com.sa
general@unimed.com.sa

**Urgent Field Safety Notice
Unicryl M Surgical Suture Recall**

Date: 13/11/2016

Unimed Reference: FSN05

Dear NCMDR Team

This letter to inform you that we have identified that the Suture Breaking Strength in few incidents are non compliant to European Pharmacopeia requirements

We are initiating a voluntary recall for the following batches:

Serial	Batch#	Item Code	Description
1	10370-38B	DW10800S7	Unicryl M 16mm 1/2c.r.b. 1(5/0) 70cm
2	10971-01	DW10800S7	Unicryl M 16mm 1/2c.r.b. 1(5/0) 70cm
3	10370-37D	DW10801	Unicryl M 16mm 1/2c.r.b. 1.5(4/0) 70cm
4	10925-01	DW10861	Unicryl M 16mm curv.cutt. 1.5(4/0) 70cm
5	10370-49B	SD10040	Unicryl M 25mm 1/2c.r.b. 1(5/0) 75cm
6	10370-49C	SD10040	Unicryl M 25mm 1/2c.r.b. 1(5/0) 75cm
7	10370-48B	SD10034S7	Unicryl M 25mm 1/2c.r.b. 1.5(4/0) 75cm
8	10370-46	SD10036	Unicryl M 25mm 1/2c.r.b. 2(3/0) 75cm
9	10370-46A	SD10036	Unicryl M 25mm 1/2c.r.b. 2(3/0) 75cm
10	10370-46H	SD10036	Unicryl M 25mm 1/2c.r.b. 2(3/0) 75cm
11	10370-47A	DW10952S	Unicryl M 30mm 1/2c.r.b. 2(3/0) 75cm
12	11101-44A	DW10952S	Unicryl M 30mm 1/2c.r.b. 2(3/0) 75cm
13	10829-01	DW10953	Unicryl M 30mm 1/2c.r.b. 3(2/0) 70cm
14	10370-44C	DW10953S	Unicryl M 30mm 1/2c.r.b. 3(2/0) 75cm
15	10370-44D	DW10953S	Unicryl M 30mm 1/2c.r.b. 3(2/0) 75cm
16	11214-01	DW10954	Unicryl M 30mm 1/2c.r.b. 3.5(0) 70cm
17	11169-01	DW10954	Unicryl M 30mm 1/2c.r.b. 3.5(0) 70cm
18	10370-41D	SD10149	Unicryl M 36mm 1/2c.r.b. 3.5(0) 90cm
19	10370-41E	SD10149	Unicryl M 36mm 1/2c.r.b. 3.5(0) 90cm
20	10846-07	DW10514	Unicryl M 40mm 1/2c.r.b. 3.5(0) 70cm
21	11044-03	DW10515LS1	Unicryl M 40mm 1/2c.r.b. 4(1) 1.5m loop
22	10839-02	DW10515	Unicryl M 40mm 1/2c.r.b. 4(1) 90cm
23	11076-01	DT10417S9	Unicryl M 40mm 1/2c.taper.cut.4(1) 90cm
24	10370-07D	DW10064LS	Unicryl M 50mm 1/2c.r.b. 3.5(0) 1.5m loop
25	10370-07E	DW10064LS	Unicryl M 50mm 1/2c.r.b. 3.5(0) 1.5m loop
26	9755-35A	DW10800S7	Unicryl M 16mm 1/2c.r.b. 1(5/0) 70cm
27	9755-46D	SD10040	Unicryl M 25mm 1/2c.r.b. 1(5/0) 75cm
28	10370-49	SD10040	Unicryl M 25mm 1/2c.r.b. 1(5/0) 75cm
29	10370-49	SD10040	Unicryl M 25mm 1/2c.r.b. 1(5/0) 75cm

30	9755-45A	SD10034S7	Unicryl M 25mm 1/2c.r.b. 1.5(4/0) 75cm
31	9755-45B	SD10034S7	Unicryl M 25mm 1/2c.r.b. 1.5(4/0) 75cm
32	9755-43B	SD10036	Unicryl M 25mm 1/2c.r.b. 2(3/0) 75cm
33	9755-43C	SD10036	Unicryl M 25mm 1/2c.r.b. 2(3/0) 75cm
34	10467-01	SD10010	Unicryl M 25mm 1/2c.r.b. 3(2/0) 70cm
35	10657-03	DW10943S	Unicryl M 26mm 1/2c.r.b. 3(2/0) 75cm
36	10370-45	DW10943S9A	Unicryl M 26mm(Black) 1/2c.r.b. 3(2/0) 90cm
37	9755-44C	DW10952S	Unicryl M 30mm 1/2c.r.b. 2(3/0) 75cm
38	10306-04	DW10953	Unicryl M 30mm 1/2c.r.b. 3(2/0) 70cm
39	9755-41C	DW10953S	Unicryl M 30mm 1/2c.r.b. 3(2/0) 75cm
40	10370-44	DW10953S	Unicryl M 30mm 1/2c.r.b. 3(2/0) 75cm
41	10370-44A	DW10953S	Unicryl M 30mm 1/2c.r.b. 3(2/0) 75cm
42	9755-38C	SD10149	Unicryl M 36mm 1/2c.r.b. 3.5(0) 90cm
43	10370-41A	SD10149	Unicryl M 36mm 1/2c.r.b. 3.5(0) 90cm
44	9755-40B	SD10151	Unicryl M 36mm 1/2c.r.b. 4(1) 90cm
45	9755-40C	SD10151	Unicryl M 36mm 1/2c.r.b. 4(1) 90cm
46	10370-43	SD10151	Unicryl M 36mm 1/2c.r.b. 4(1) 90cm
47	10370-43A	SD10151	Unicryl M 36mm 1/2c.r.b. 4(1) 90cm
48	10284-05	DW10515	Unicryl M 40mm 1/2c.r.b. 4(1) 90cm
49	10685-04	DW10515	Unicryl M 40mm 1/2c.r.b. 4(1) 90cm
50	10486-01	DW10514LHS	Unicryl M 40mm 1/2c.r.b.hvy. 3.5(0) 1.5m loop
51	10562-01	DT10417S9	Unicryl M 40mm 1/2c.taper.cut.4(1) 90cm

Customer Action

Please review your inventory for products affected by this issue and perform the following actions:

1. Immediately identify and quarantine all unused, affected products in your inventory.
2. Return all unused, affected products in your inventory to Unimed. Our Representatives will assist you in the return of affected product as necessary.

Please share this notification with others in your organization as appropriate. If product within scope of this recall has been forwarded to another facility, please notify the facility of the issue and assist with the inconvenience that it may cause

Please be assured that patient safety and products quality remain our primary concern.

If you have any questions, please contact your Unimed Representatives.

Sincerely,

Nader Hamdan Assiri

