

Urgent Field Safety Notice **Covidien DAR™ airway products** **Customer Notification**

April 2021

Medtronic Reference FA975

Dear Customer,

The purpose of this letter is to advise you of an important Field Safety Notice for certain production lots of Medtronic's Covidien DAR™ airway products.

Issue Description:

Medtronic was notified by our supplier of sterilization services, Steril Milano, of potential deviations from validated parameters for ethylene oxide sterilization. The deviations affect certain production lots of Medtronic's Covidien DAR™ airway products and occurred between March 2018 and February 2021. Some of these impacted lots were distributed to your facility.

These deviations were the subject of Circulars by the Italian Ministry of Health, dated 11 March and 30 March 2021. In accordance with the 30 March 2021, Circular, Medtronic is notifying customers who have received potentially impacted lots and is requesting that lots be quarantined and held at their facility pending further notice by Medtronic.

Impacted lots that are within Medtronic control have been quarantined and sterilization activities with Steril Milano have ceased.

The full impact of the deviations is still being investigated, and Medtronic is analyzing data to assess whether sterilization has been achieved despite the deviations. We are working with our notified body and competent authorities to review our approach, and whether lots under quarantine might be released if sufficient data is available. Medtronic is aware that adequate supply of these devices is critical in the pandemic and we are working to identify solutions.

This Field Safety Notice affects only the item codes and lot numbers listed in Attachment A.

Actions you should take:

1. Immediately quarantine and discontinue use of item codes with associated lot numbers listed in Attachment A.
2. If you have distributed the devices listed in Attachment A, please promptly forward the information from this letter to those recipients.

Actions being taken by Medtronic:

1. Medtronic is investigating the sterilization records for affected production lots to determine whether sterilization has been achieved despite the deviations.
2. Medtronic will communicate with affected customers regarding the finding of our investigation once complete and will provide direction on what customers should do with the devices quarantined because of this Field Safety Notice.

Medtronic has notified the Competent Authority of your country of this action.

We regret any inconvenience this may cause. We are committed to patient safety and appreciate your prompt attention to this matter. If you have any questions regarding this communication, please contact your Medtronic representative.

Sincerely,

Asi, Ala

Identifying Affected Product

Item code

Lot number

COVIDIEN positive results for life®
DAR™
 Adult-Pediatric Electrostatic Filter HME
 Small
 Vt 150-1200 mL

REF 352/5877Z
LOT XXXXXXXXXX
 Use-by YYYY-MM-DD

COVIDIEN positive results for life®
DAR™
 Adult-Pediatric Electrostatic Filter HME
 Small
 Vt 150-1200 mL

REF 352/5877Z
LOT XXXXXXXXXX
 Use-by YYYY-MM-DD

Electrostatisk filter-fugtvarmeveksler til voksne og børn Lille
 Adult-pediatrik elektrostatisk filter HME Liten
 Ηλεκτροστατικό φίλτρο HME ενήλικων/ παιδιατρικό Μικρό
 Felnőt-gyermek elektrostatikus szűrő HME Kicsi
 Filtr elektrostatyczny dla dzieci i dorosłych, z wymiennikiem ciepła/wilgoci Mały
 Elektrostatiský filtr HME pro dospěle a děti Malá
 FBO электростатического фильтра, для взрослых и детей Малый
 Yetişkin-Pediyatrik Elektrostatik Filtre HME (Isı ve Nem Değiştirici) Küçük

STERILE EO **Rx ONLY**
 Not made with natural rubber latex
 Single use
 Not made with DEHP
 Caution, consult accompanying documents
 CE 0123

© 2011 Covidien. Made in Italy.
 Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA
 Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.
 10071182 Rev C

Item Code	Description	Affected Lot Numbers					
031122 / 031122-25	STERIVENT MINI X MTP FILTER	18L0473FAX	18L1287FAX	19E0682FAX	20A0688FAX	20E1017FAX	20F1129FAX
		18L0474FAX	19C0145FAX	19J0567FAX	20A0689FAX	20E1018FAX	20F1131FAX
		18L0475FAX	19D0638FAX	19K0524FAX	20B0623FAX	20F0942FAX	
		18L1286FAX	19E0681FAX	19K1052FAX	20C0679FAX	20F0943FAX	
285/25035	Adult Breathing Circuit with Exhalation Valve	20L0611FAX	20M1116FAX				
285/25038	Adult Breathing Circuit with Exhalation Valve	20L0612FAX	20M1118FAX				
285/25390	Adult Breathing Circuit	21A0011FAX					
285/25408	3-in-1 Coaxial Circuit	19E0725FAX	20G0975FAX	20G1024FAX			
285/25417	3-in-1 Coaxial Circuit	19A1305FAX	20G1079FAX	20L1104FAX			
285/25446	3-in-1 Coaxial Circuit	20D0415FAX	20G1077FAX	20K1133FAX	20E0204FAX	20J0557FAX	
285/25517	Adult Breathing Circuit with Exhalation Valve	20G0426FAX					
285/25540	Adult Breathing Circuit	20G0425FAX	20L1155FAX	20L1156FAX	20L1157FAX	21A0015FAX	21A0014FAX
		21A0016FAX	21A1135FAX				
285/25541	Adult Filtered Breathing Circuit	20M0037FAX					
285/25557	Adult Coaxial Breathing Circuit	19D0464FAX					
285/5058	Adult Breathing Circuit	20L0020FAX					
285/5059	Adult Breathing Circuit	20H0850FAX	20J0692FAX	20M0437FAX			
285/5060	Adult Breathing Circuit	20J0120FAX					
285/5061	Adult Breathing Circuit	20G0080FAX					
285/5062	Adult Breathing Circuit	20M0028FAX					
285/5063	Adult Breathing Circuit	18B0369FAX	20F0923FAX	20J0029FAX	20K1255FAX	20M1047FAX	20M1048FAX
		20M1049FAX	20M1091FAX	20M1092FAX	21A0844FAX		
285/5063Z	Adult Breathing Circuit	20G1120FAX					
285/5064	Adult Breathing Circuit	18B0683FAX	20K1253FAX	20L1124FAX	20L1126FAX	20M1050FAX	20M1052FAX
		18B1074FAX	20K1254FAX	20L1125FAX	20L1154FAX	20M1051FAX	20M1053FAX
		20F0122FAX					
285/7030	Adult Breathing Circuit	20J0691FAX					
285/7032	Adult Breathing Circuit	19A1212FAX	19E0477FAX	20A0313FAX			
285/7056	Adult Breathing Circuit with Exhalation Valve	21A0013FAX	20L0616FAX				
285/7147	Adult Breathing Circuit with Exhalation Valve	20L0615FAX					
285/7185	Adult Breathing Circuit	20G0421FAX					
285/7227	Adult Breathing Circuit	20L0617FAX					
285/7283	Adult Breathing Circuit	20C0459FAX					
285/7413	Adult Breathing Circuit	20C0466FAX					
285/7515	Adult Breathing Circuit	18E1108FAX	20J0123FAX	20L1049FAX			
285/7627	Adult Mapleson E Ayres Circuit	20G0935FAX					
285/7770	Adult Breathing Circuit	20L1128FAX	20L1129FAX				
285/7785	Adult Breathing Circuit	20M0479FAX					
285/7786	Adult Breathing Circuit	20M0480FAX					
285/7945	Adult Bain Coaxial Circuit	20M0035FAX					
285S25640	Adult Breathing Circuit	19E0480FAX	20L1052FAX	20L1255FAX			
285S25726	Adult Breathing Circuit with Aerosol O2 Mask	20M1100FAX					
285S25735	Adult Breathing Circuit	20G0936FAX					
285S25762	Adult Breathing Circuit	20L1287FAX	20M0607FAX				
285S25790	Adult Breathing Circuit	21A0052FAX					
285S25791	Adult Breathing Circuit	20L0625FAX					
285S5059Z	Adult Breathing Circuit	20K0254FAX					
285S5064Z	Adult Breathing Circuit	19A1415FAX					
285S7031Z	Adult Breathing Circuit	20J0039FAX					
285SX0005	Adult Breathing Circuit with Exhalation Valve	20L0970FAX	20M1120FAX				
286/25190	Pediatric Breathing Circuit	20K0944FAX					
286/5065	Pediatric Breathing Circuit	20F0402FAX					
286/5066	Pediatric Breathing Circuit	20F0401FAX	20J1175FAX				
286/5068	Pediatric Breathing Circuit	20E0435FAX	20M0131FAX				
286/5069	Pediatric Breathing Circuit	20K1187FAX	20K1188FAX				
286/5070	Pediatric Breathing Circuit	20K0828FAX					
286/5071	Pediatric Breathing Circuit	20H0392FAX					
286/7359	Pediatric Bain Coaxial Circuit	21A0023FAX					
286/7669	Infant Mapleson E Ayres Circuit	20L0627FAX					
287/25010	Guedel Airway 3 CM M.000	19A0646FAX	19G0580FAX	19L1322FAX			
287/25011	Guedel Airway 4 CM M.00	19G0581FAX	19L1320FAX				
287/25012	Guedel Airway 5 CM M.0	19L1319FAX					
287/25013	Guedel Airway 6 CM M.1	19L1321FAX					
287/25014	Guedel Airway 7 CM M.2	19A0647FAX	19A0648FAX	19A1373FAX	19A1374FAX	19G0582FAX	19L1323FAX
287/25015	Guedel Airway 8 CM M.3	18L0308FAX	19D0773FAX	19E1075FAX	19L1324FAX		
287/25016	Guedel Airway 9 CM M.4	19G0583FAX	19G0630FAX	19L1325FAX			
287/25017	Guedel Airway 10 CM M.5	19L1326FAX					
287/25018	Guedel Airway 11 CM M.6	19G0584FAX	19L1327FAX				
287/7631	Guedel Airway MIS. 10 CM	18K1255FAX	19B0385FAX	19J0371FAX	19L0879FAX	20G1268FAX	20K0457FAX
		18K1256FAX	19D0834FAX	19J0372FAX	20B0455FAX	20H0535FAX	20K0458FAX
		18M0557FAX	19D1340FAX	19K0437FAX	20B0456FAX	20H0536FAX	20K1065FAX
		18M0558FAX	19D1341FAX	19K0438FAX	20B0662FAX	20J0481FAX	20K1066FAX
		18M0922FAX	19D1342FAX	19K0439FAX	20C0535FAX	20J0482FAX	20L0402FAX

Item Code	Description	Affected Lot Numbers					
		19A0171FAX	19E1481FAX	19K1259FAX	20C0536FAX	20J0483FAX	20M0614FAX
		19A0172FAX	19E1482FAX	19K1260FAX	20C0929FAX	20J1288FAX	21A0881FAX
		19A0173FAX	19G1400FAX	19K1261FAX	20G1265FAX	20K0377FAX	21A0882FAX
		19B0383FAX	19J0369FAX	19K1262FAX	20G1266FAX	20K0378FAX	21A1001FAX
		19B0384FAX	19J0370FAX	19L0878FAX	20G1267FAX	20K0456FAX	21A1002FAX
		18K1156FAX	19B1017FAX	19F0313FAX	19L0871FAX	20D0456FAX	20K0983FAX
		18L0368FAX	19B1018FAX	19F0666FAX	19M0311FAX	20D0457FAX	20K1572FAX
		18L0369FAX	19B1019FAX	19F0667FAX	19M0312FAX	20D0811FAX	20L0123FAX
		18L0370FAX	19B1152FAX	19F1241FAX	20A0356FAX	20G0673FAX	20L0722FAX
		18L0371FAX	19B1153FAX	19F1249FAX	20A0358FAX	20G0674FAX	20L0723FAX
		18L0697FAX	19C0117FAX	19H0870FAX	20A0571FAX	20G0675FAX	20L0724FAX
		19A0167FAX	19D0127FAX	19H0871FAX	20A0572FAX	20G0676FAX	20L0725FAX
		19A0168FAX	19D0128FAX	19H0872FAX	20A0573FAX	20G1177FAX	20L0726FAX
		19A0169FAX	19E0780FAX	19H0873FAX	20B0110FAX	20H0534FAX	20L0727FAX
		19A0292FAX	19E0781FAX	19J1314FAX	20B0111FAX	20H0597FAX	20L0728FAX
		19A1361FAX	19E0782FAX	19J1315FAX	20B1041FAX	20J0680FAX	20M0366FAX
		19A1362FAX	19E0783FAX	19J1316FAX	20C0883FAX	20J1339FAX	20M1492FAX
		19A1363FAX	19F0088FAX	19L0249FAX	20C0884FAX	20J1340FAX	20M1493FAX
		19A1364FAX	19F0089FAX	19L0250FAX	20C0885FAX	20J1341FAX	20M1494FAX
		19A1365FAX	19F0090FAX	19L0424FAX	20C0886FAX	20K0981FAX	20M1497FAX
		19B1016FAX	19F0233FAX	19L0425FAX	20D0455FAX	20K0982FAX	21A0977FAX
		18L0755FAX	19A1289FAX	19E0444FAX	19J0634FAX	20D0060FAX	20K0376FAX
		18L0756FAX	19A1290FAX	19E0445FAX	19J1084FAX	20D0618FAX	20L0720FAX
		18L0757FAX	19B0002FAX	19E0446FAX	19J1085FAX	20D0619FAX	20M0817FAX
		18L0758FAX	19B0168FAX	19E0447FAX	20A0964FAX	20D0620FAX	20M0818FAX
		18L0759FAX	19C0430FAX	19E1128FAX	20A0965FAX	20G0671FAX	20M0819FAX
		18L0760FAX	19C0431FAX	19E1133FAX	20A0966FAX	20H0439FAX	21A0464FAX
		18M0923FAX	19C0432FAX	19E1134FAX	20A0967FAX	20H0440FAX	21A0465FAX
		18M0924FAX	19C0526FAX	19E1135FAX	20A0968FAX	20H0459FAX	21A0466FAX
		19A0641FAX	19D0536FAX	19G0785FAX	20B0029FAX	20J0974FAX	21A0467FAX
		19A0642FAX	19D0537FAX	19H0751FAX	20C0495FAX	20J0975FAX	21A0867FAX
		19A0916FAX	19D1160FAX	19H0752FAX	20D0059FAX	20J0976FAX	21A1004FAX
		19A1288FAX	19E0443FAX				
		18L0825FAX	19B0724FAX	19J0241FAX	20A0962FAX	20G0994FAX	20M1235FAX
		18M0596FAX	19C1390FAX	19L0500FAX	20B0173FAX	20H1044FAX	20M1310FAX
		19B0167FAX	19D0180FAX	19L0501FAX	20C0097FAX	20J1110FAX	
		19B0607FAX	19G1116FAX	19L1175FAX	20C0098FAX	20M0889FAX	
		19B0166FAX	19D1164FAX	19L0148FAX	20B0142FAX	20G0995FAX	20M1371FAX
		19B1099FAX	19E1397FAX	20A0398FAX	20B0587FAX	20H0801FAX	21A0107FAX
		19B1242FAX	19K0580FAX	20A0574FAX			
		18M0595FAX	19F0883FAX	20A0359FAX	20A1431FAX	20G0996FAX	20M0616FAX
		19A0175FAX	19G0577FAX	20A0397FAX	20A1432FAX	20L0401FAX	20M1311FAX
		19E1396FAX	19L0187FAX	20A1157FAX	20B1345FAX		
		18M0594FAX	19F0761FAX	19K0001FAX	19L1224FAX	20A0799FAX	20H0710FAX
		19A0644FAX	19G0579FAX	19L0517FAX	20A0651FAX	20A0969FAX	20J1109FAX
		19B0916FAX					
		19A0174FAX	19B1238FAX	19J0143FAX	19L0943FAX	20B0453FAX	20G0997FAX
		19B0771FAX	19G0578FAX	19L0144FAX	20B0333FAX	20C1089FAX	20H0802FAX
		19F0315FAX	19L0248FAX	20A1028FAX	20C0954FAX	20H0448FAX	20M0890FAX
		19J0448FAX	20A0699FAX	20C0080FAX	20D0616FAX		
291/25488	CO2 Sampling Line	20L0747FAX					
291/7617	Suction Valve	20J1258FAX	20M0387FAX				
291/9250	Bifurcated Connector	20D1496FAX					
300/13532	Adult Breathing Circuit	21A0292FAX					
300/13620	Adult Breathing Circuit	18E0943FAX	19C1483FAX	20A0318FAX	20B0401FAX	20C0337FAX	
300/13774	Adult Breathing Circuit	20L1313FAX					
300/13783	Adult Breathing Circuit	20H0098FAX	20H0100FAX	20H0824FAX	20J0018FAX		
300/13815	Adult Breathing Circuit	19C1486FAX					
300/13957	Adult Breathing Circuit	20G0954FAX					
300/14147	Adult Breathing Circuit	19B0186FAX	20H0089FAX	20H0090FAX	20H1228FAX	20H1229FAX	20K0048FAX
300/6037	Adult Breathing Circuit	20M0002FAX					
300/6042	Adult Breathing Circuit	19M0810FAX	19M0811FAX	20K0052FAX	20K0053FAX	20K0055FAX	20K1513FAX
		20M1108FAX					
		18K0565FAX	19C0995FAX	19F1163FAX	20A0332FAX	20E1185FAX	20L0203FAX
		19A0976FAX	19C1012FAX	19F1192FAX	20D0638FAX	20J0672FAX	20L1148FAX
		19A1164FAX	19C1013FAX	19F1193FAX	20D0639FAX	20J0673FAX	20L1149FAX
		19A1165FAX	19C1027FAX	20A0329FAX	20D0640FAX	20J0675FAX	20L1152FAX
		19A1166FAX	19D0099FAX	20A0330FAX	20D0641FAX	20J0676FAX	20L1318FAX
		19B0906FAX	19D0100FAX	20A0331FAX	20D0644FAX	20J0772FAX	20L1319FAX
		19B0907FAX	19D0167FAX				
		18B0296FAX	19B0910FAX	19F1232FAX	20C0282FAX	20F0644FAX	20M0005FAX
		18E0612FAX	19C0867FAX	19F1233FAX	20C0322FAX	20F0750FAX	20M0006FAX
		18F0182FAX	19C0868FAX	19G0007FAX	20C0324FAX	20J1143FAX	20M0007FAX
		18K0566FAX	19F1209FAX	19G0008FAX	20C0325FAX	20J1144FAX	20M0008FAX
		19A1414FAX	19F1210FAX	19M0787FAX	20D0513FAX	20K0488FAX	20M0536FAX

Item Code	Description	Affected Lot Numbers					
		19A1421FAX	19F1211FAX	19M0788FAX	20D0514FAX	20L0198FAX	20M0537FAX
		19A1422FAX	19F1212FAX	19M0789FAX	20D0612FAX	20L0200FAX	20M1417FAX
		19B0909FAX	19F1231FAX	20B0423FAX	20D0637FAX	20L0201FAX	21A0001FAX
300/6049	Adult Breathing Circuit	20A0208FAX	20A0209FAX	20M0018FAX			
300/6081	Adult Breathing Circuit	19K1436FAX					
300/6099	Adult Breathing Circuit	20F0351FAX	20K0333FAX	21A0287FAX			
300/6211	Adult Breathing Circuit	20M0038FAX					
300/6235	Adult Breathing Circuit	19C1078FAX	20L1158FAX	20L1159FAX	20L1160FAX		
300/6259	Adult Breathing Circuit	20L1245FAX					
300/6365	Adult Breathing Circuit	17C0486FAX	20M1123FAX				
300/6449	Adult Breathing Circuit	19B0174FAX	19E0620FAX	20A1229FAX	20K0335FAX	20M0016FAX	21A0288FAX
		19D0312FAX	20A1228FAX	20K0334FAX			
300/6484	Adult Breathing Circuit	19D0313FAX	20H1273FAX				
300/6567	Adult Breathing Circuit	20H1275FAX	20M0011FAX	21A0290FAX			
300/6644	Adult Breathing Circuit with Filter HME	20D1108FAX					
300S14224	Adult Breathing Circuit	20E0589FAX	20L1256FAX				
300S14269	Adult Breathing Circuit	18B0180FAX	18E1555FAX	20B0131FAX	20E0450FAX	20G1598FAX	20M0844FAX
		18B0181FAX	19K0892FAX	20C0291FAX	20E0453FAX	20H0082FAX	20M0845FAX
		18B0763FAX	20A0118FAX	20C0292FAX	20G1596FAX	20H0083FAX	20M0846FAX
		18B0764FAX	20B0068FAX	20C0293FAX	20G1597FAX	20L1072FAX	20M0951FAX
		18B0765FAX	20B0069FAX	20C0302FAX			
300S14353	Adult Breathing Circuit	19B0346FAX					
300S14385	Adult Breathing Circuit	19K1077FAX	20A0701FAX	20L1257FAX			
300S14425	Adult Breathing Circuit	18B0781FAX					
300S6043Z	Adult Breathing Circuit	19C1302FAX	20G1121FAX				
300S6044Z	Adult Breathing Circuit	19A1419FAX	20G1122FAX				
301/13404	Adult Breathing Circuit	20K0315FAX	20M0012FAX	21A0291FAX			
301/14103	Adult Breathing Circuit	20C0288FAX	20M1307FAX				
301/14105	Adult Breathing Circuit	18J1171FAX	20H1236FAX	20K1257FAX			
301/6326	Adult Breathing Circuit	20J0020FAX	20K1511FAX	20L1315FAX			
301/6378	Adult Breathing Circuit	20A1236FAX					
301S14401Z	Adult Breathing Circuit with Filter HME	20J0670FAX					
303/5140	Adult Breathing Circuit	19B0489FAX	20B0408FAX				
303/6057	Adult Breathing Circuit	20D0476FAX	20L0773FAX	20M0605FAX			
303/6062	Adult Breathing Circuit	20H0207FAX	20M1308FAX				
303/6090	Adult Breathing Circuit	19D0141FAX					
303/6112	Adult Breathing Circuit	20L1046FAX	20L1047FAX				
303/6117	Adult Breathing Circuit	20L1345FAX					
303/6298	Adult Breathing Circuit	20A0756FAX	20K1119FAX	21A0098FAX			
303/6376	Adult Breathing Circuit	20L1073FAX					
303S14127Z	Adult Breathing Circuit	20J1382FAX					
303S14424	Adult Breathing Circuit with Water Traps	20A0705FAX	20E0590FAX	20M1304FAX			
304/13031	Adult Single-Limb Heated Circuit	19B0625FAX	19D0650FAX	20B0410FAX	20L1174FAX	20L1342FAX	
304/13070	Adult Single-Limb Heated Circuit	20L0866FAX					
304/13358	Adult Single-Limb Heated Circuit	19C0643FAX					
304/13482	Adult Breathing Circuit	19B0515FAX					
304/13567	Adult Breathing Circuit	21A0111FAX					
304/13686	Adult Breathing Circuit	18B0766FAX	18B0767FAX	18E0658FAX	18E0659FAX	18E0660FAX	19F1151FAX
		19F1152FAX	19F1183FAX	19F1184FAX	20A0722FAX		
304/13693	Adult Breathing Circuit	18B0182FAX	20A0720FAX	20K1123FAX	20M0032FAX		
304/14053	Adult Single-Limb Heated Circuit	20K1499FAX	20L0865FAX	21A1025FAX			
304/14054	Adult Single-Limb Heated Circuit	19B0306FAX	20M0979FAX				
304/14145	Adult Single-Limb Heated Circuit	20K0360FAX					
304/6085	Adult Breathing Circuit	18F0247FAX	19D0398FAX	19D0617FAX	19L0070FAX	20K1260FAX	20L1103FAX
		18F0248FAX	19D0399FAX	19E0347FAX	19L0071FAX	20L1101FAX	20L1167FAX
		19B0486FAX	19D0589FAX	19E0348FAX	20K1015FAX	20L1102FAX	21A0054FAX
304/6400	Adult Breathing Circuit	20A0721FAX	20K1266FAX	20M1423FAX			
304/6404	Adult Breathing Circuit	20L1166FAX					
304/6486	Adult Breathing Circuit	21A0099FAX					
304/6518	Adult Breathing Circuit	19C1080FAX					
304S14262	Adult Single-Limb Heated Circuit	19C0644FAX	20A0561FAX	20M1176FAX			
304S14263	Adult Dual-Limb Heated Circuit	20A1330FAX	20D0537FAX	21A0352FAX			
304S14291	Adult Single-Limb Heated Circuit	20M1179FAX					
304S14292	Adult Single-Limb Heated Circuit	20M1330FAX					
304S14292Z	Adult Single-Limb Heated Circuit	20L0872FAX					
304S14300	Adult Dual-Limb Heated Circuit	20M0621FAX	20M0622FAX	20M0623FAX			
304S14313	Adult Single-Limb Heated Circuit	20K0382FAX					
304S14365	Adult Breathing Circuit	19K1166FAX					
304S14402Z	Adult Single-Limb Heated Circuit	19B0630FAX	20L0979FAX				
304S14416	Adult Single-Limb Heated Circuit	20L0861FAX					
305/6712	Pediatric Breathing Circuit	20J0127FAX	20J0128FAX				
305/6713	Pediatric Breathing Circuit	20A0304FAX	20K1504FAX				
305/6714	Pediatric Breathing Circuit	20A0305FAX	20J0151FAX	20M0707FAX			

Item Code	Description	Affected Lot Numbers					
305/6722	Pediatric Breathing Circuit	19E0133FAX	20M0706FAX				
305/6852	Infant Breathing Circuit	20G0924FAX					
305/6907	Pediatric Breathing Circuit	20L1248FAX					
305/6980	Infant Breathing Circuit	20L0247FAX					
305/8173	Infant Breathing Circuit	20M0759FAX	20M1127FAX	20M1128FAX			
305/8195	Infant Breathing Circuit	19K1439FAX	20A0593FAX	20J0202FAX	20M0608FAX		
305/8483	Infant Breathing Circuit	19K1440FAX					
305/8773	Pediatric Breathing Circuit	18J1127FAX	18J1128FAX	18J1129FAX	20B0541FAX	20B0542FAX	
305/8804	Pediatric Breathing Circuit	20K1510FAX	20M0704FAX				
305G8926	Pediatric Breathing Circuit with Breathing Bag	20B0562FAX					
305S8934	Infant Breathing Circuit	20M0556FAX					
305S8935	Pediatric Breathing Circuit	20D0344FAX	20D0345FAX				
305S9928	Infant Breathing Circuit	20M0753FAX					
306/6741	Pediatric Breathing Circuit	20M0048FAX					
306/6843	Pediatric Single-Limb Heated Circuit	20B0412FAX	20J0892FAX				
306/6948	Pediatric Breathing Circuit	18E1552FAX	19C0730FAX	19E0502FAX	20M0840FAX		
306G8925	Pediatric Breathing Circuit	20M0426FAX					
306S9904	Pediatric Single-Limb Heated Circuit	20B0414FAX					
307/6747	Infant Breathing Circuit	20G1575FAX	20H0765FAX				
307/8447	Infant Single-Limb Heated Circuit	20J0781FAX					
307/8448	Infant Breathing Circuit	20M0748FAX					
307/8682	Infant Single-Limb Heated Circuit	20M1062FAX					
309/13432	Adult Mapleson C Manual Breathing Unit	20M0855FAX					
309/6310	Adult Mapleson C Manual Breathing Unit	18B0012FAX	18E1462FAX	19D0114FAX	19E0709FAX	20J0874FAX	20M0110FAX
		18B1156FAX	19D0113FAX	19E0708FAX	20D0471FAX	20L0392FAX	20M0520FAX
		18E1042FAX					
309/6351	Infant Mapleson C Manual Breathing Unit	20C0312FAX	20C0313FAX				
309/6353	Pediatric Mapleson C Manual Breathing Unit	20M0950FAX					
309/6529	Adult Mapleson C Manual Breathing Unit	20J1395FAX					
309/6557	Adult Mapleson C Manual Breathing Unit	18E0849FAX					
309/6563	Pediatric Mapleson C Manual Breathing Unit	20G0769FAX					
309/6612	Adult Mapleson C Manual Breathing Unit	18E1257FAX	20C0387FAX	20K0565FAX	20L0426FAX	20M0521FAX	20M0538FAX
		19C1033FAX	20C0565FAX	20L0424FAX	20L0427FAX	20M0522FAX	20M0712FAX
		19E0483FAX	20G1006FAX	20L0425FAX	20L0428FAX	20M0523FAX	20M0713FAX
309S14165	Adult Mapleson C Manual Breathing Unit	19E0482FAX	20G1014FAX	20K0809FAX			
309S14166	Adult Mapleson C Manual Breathing Unit	20J1191FAX					
309S14217	Adult Mapleson C Manual Breathing Unit	19E0363FAX	20B0224FAX	20G1017FAX			
309S14219	Pediatric Mapleson C Manual Breathing Unit	19E0361FAX	20G1016FAX	20L1105FAX			
330/5110	Catheter Mount 15CM 22M/15F-22F	20H0411FAX					
330/5337	Catheter Mount 22M/15F	20G1638FAX	20G1639FAX				
331/5113	Catheter Mount CM.15 22M/15F-22F	19B1334FAX					
331/5305	Catheter Mount CM.10 22M/15F-15M	19C0245FAX					
331/5353	Catheter Mount CM.5 22M/15F-15M	20H0410FAX					
331/5373	Catheter Mount CM.20 15M.RAC.GIR.	20H0427FAX					
331/5385	Catheter Mount CM.15 22M/15F-22F	20H0421FAX	20H0422FAX	20H0423FAX	20J1311FAX	20J1313FAX	20L0668FAX
		20L0692FAX	20L1010FAX	20M0449FAX	20M0450FAX		
331/5389	Catheter Mount RAC.GIR. CM.10	18E1567FAX	20B0603FAX	20B0604FAX	20H0680FAX	20H0681FAX	20H0682FAX
		20J1014FAX	20J1015FAX	20J1017FAX	20J1208FAX	20J1209FAX	
		17E1483FAX	19L1401FAX	20H0417FAX	20J1216FAX	20K1585FAX	20L1306FAX
331/5390	Catheter Mount GIREV.T. GOM.15 CM	17E1484FAX	19L1402FAX	20H0418FAX	20K0639FAX	20L0703FAX	20L1307FAX
		17E1485FAX	20H0415FAX	20H0547FAX	20K0973FAX	20L0704FAX	20M0454FAX
		17E1486FAX	20H0416FAX	20J1215FAX	20K1584FAX	20L0705FAX	20M0456FAX
		19L1400FAX					
331/5391	Catheter Mount GIREV.T.GOM. 10CM.	20J1277FAX	20J1279FAX	20K0652FAX	20L0707FAX	20L0708FAX	20L0709FAX
		20J1278FAX	20J1280FAX	20L0706FAX			
331/5395	Catheter Mount CM.20 22M/15F-22F	19L0318FAX	20C0962FAX	20L0308FAX	20L0310FAX	20M0663FAX	20M0664FAX
		20A0825FAX	20G1235FAX	20L0309FAX			
331/5661	Catheter Mount Extendable Swivel	20C0774FAX	20H0419FAX	20H0420FAX	20H0696FAX	20H0697FAX	20H0698FAX
		20H0699FAX	20H1185FAX	20J1318FAX	20M0860FAX		
331/5662	Catheter Mount Extendable	20B1105FAX	20B1106FAX	20H0404FAX	20J0416FAX		
331/5669	Catheter Mount Extendable Swivel	19B1191FAX	19M0117FAX	20D0594FAX	20H0499FAX	20J1001FAX	20J1003FAX
		20J1004FAX	20K1575FAX	20K1576FAX	20K1577FAX		
332/5115	Catheter Mount 15CM 15F 15M	20H0721FAX					
332/5310	Catheter Mount CM.15 15F-22F	20D0723FAX	20H0413FAX	20L0698FAX			
332/5332	Catheter Mount 20CM 15F 15M	20A0803FAX					
332/5663	Catheter Mount Extendable	19L1385FAX	20D0599FAX	20J1040FAX	20J1041FAX	20J1273FAX	20J1274FAX
350/19003	Electrostatic Filter	18L1143FAX	19D0326FAX	20C0197FAX	20D0364FAX	20D1619FAX	20G0319FAX
		19C0790FAX	19D0327FAX	20C0198FAX	20D0365FAX	20G0176FAX	20G0320FAX
		19C0791FAX	19G0468FAX	20D0246FAX	20D0647FAX	20G0177FAX	20G0321FAX
		19C0792FAX	19G0469FAX	20D0247FAX	20D0864FAX	20G0220FAX	20G0322FAX
		19C0793FAX	19G0470FAX	20D0248FAX	20D1052FAX	20G0221FAX	20G0323FAX
		19C0794FAX	20A0575FAX	20D0249FAX	20D1053FAX	20G0222FAX	20G0324FAX
		19C0795FAX	20A1013FAX	20D0250FAX	20D1530FAX	20G0317FAX	20G0325FAX
		19D0325FAX	20A1014FAX	20D0363FAX	20D1618FAX	20G0318FAX	

Item Code	Description	Affected Lot Numbers					
350/5420	Electrostatic Filter with Connector	19A1295FAX	19E0602FAX	19F0768FAX	19K0056FAX	20A0477FAX	20F0374FAX
		20G0178FAX					
350/5422	Electrostatic Filter	18M0400FAX	19B0903FAX	19G0690FAX	20D0866FAX	20G1520FAX	20J0047FAX
		18M0401FAX	19B0904FAX	19G0920FAX	20D0867FAX	20G1521FAX	20J0335FAX
		18M0402FAX	19C0491FAX	19G0921FAX	20D0868FAX	20G1522FAX	20J0336FAX
		19A0058FAX	19C0492FAX	19G0922FAX	20E1399FAX	20G1523FAX	20J0337FAX
		19A0160FAX	19C0493FAX	19K0433FAX	20E1400FAX	20G1602FAX	20J0338FAX
		19A0161FAX	19D0084FAX	19K0434FAX	20E1423FAX	20G1603FAX	20J0339FAX
		19A0162FAX	19E1263FAX	19K0435FAX	20E1424FAX	20G1604FAX	20J0340FAX
		19A0434FAX	19E1264FAX	19K0436FAX	20F0180FAX	20G1612FAX	20J0502FAX
		19B0726FAX	19E1266FAX	19K1040FAX	20F0181FAX	20G1613FAX	20J0503FAX
		19B0727FAX	19E1299FAX	20A0568FAX	20F0182FAX	20H0511FAX	20J0687FAX
		19B0728FAX	19E1300FAX	20A0569FAX	20F0224FAX	20H0512FAX	21A0212FAX
		19B0729FAX	19E1301FAX	20C0443FAX	20F0225FAX	20H0513FAX	21A0213FAX
		19B0730FAX	19G0687FAX	20C0444FAX	20F0226FAX	20H0514FAX	
		19B0731FAX	19G0688FAX	20C0445FAX	20F0603FAX	20H1279FAX	
		19B0902FAX	19G0689FAX	20C1292FAX	20G1507FAX	20H1280FAX	
		350/5422TC	Electrostatic Filter	20A0567FAX	20G1506FAX	20J0258FAX	20M0399FAX
350/5845	Electrostatic Filter	18L1257FAX	19G0042FAX	19L1101FAX	20D1540FAX	20G0638FAX	20J0094FAX
		18L1258FAX	19G0043FAX	20A0261FAX	20D1541FAX	20G0639FAX	20J0095FAX
		19B1374FAX	19G0044FAX	20A0262FAX	20D1542FAX	20G0698FAX	20J0096FAX
		19B1375FAX	19G0701FAX	20A1183FAX	20E0155FAX	20G1115FAX	20L0408FAX
		19B1376FAX	19G0702FAX	20A1184FAX	20E0156FAX	20H0198FAX	20M0278FAX
		19D0160FAX	19G0703FAX	20A1185FAX	20E0157FAX	20H0199FAX	20M0279FAX
		19D0161FAX	19G0704FAX	20D0379FAX	20F0533FAX	20H0200FAX	20M0983FAX
		19D1127FAX	19G0705FAX	20D0380FAX	20F0858FAX	20H0201FAX	21A0474FAX
		19D1128FAX	19H0700FAX	20D0381FAX	20F0859FAX	20H0202FAX	21A0475FAX
		19E0118FAX	19H0701FAX	20D0382FAX	20F0860FAX	20H0203FAX	21A0476FAX
		19E0119FAX	19J1325FAX	20D0383FAX	20F0861FAX	20H0204FAX	21A0477FAX
		19E0120FAX	19L0113FAX	20D0384FAX	20F0862FAX	20H0205FAX	21A1130FAX
		19E0121FAX	19L0114FAX	20D1537FAX	20G0635FAX	20J0091FAX	
		19E0267FAX	19L0115FAX	20D1538FAX	20G0636FAX	20J0092FAX	
		19E0268FAX	19L0116FAX	20D1539FAX	20G0637FAX	20J0093FAX	
		350/5845TC	Electrostatic Filter	20F0529FAX			
350/5865	Electrostatic Filter 22F	19A1296FAX	19E1655FAX	19L1352FAX	20C0652FAX	20F0253FAX	20G0892FAX
		19A1297FAX	19E1656FAX	19M0270FAX	20C0653FAX	20F1271FAX	20G1672FAX
		19A1298FAX	19G0499FAX	19M0272FAX	20C0707FAX	20F1272FAX	20G1673FAX
		19A1299FAX	19G0500FAX	19M0416FAX	20C0742FAX	20G0194FAX	20G1674FAX
		19B0045FAX	19G0501FAX	19M0925FAX	20C0743FAX	20G0195FAX	20G1675FAX
		19B0046FAX	19G0685FAX	20A0110FAX	20C0744FAX	20G0196FAX	20G1676FAX
		19B0047FAX	19H0728FAX	20A0671FAX	20C0777FAX	20G0197FAX	20H0328FAX
		19B0048FAX	19H0729FAX	20A1021FAX	20C1116FAX	20G0198FAX	20H0329FAX
		19B0425FAX	19H0730FAX	20A1022FAX	20C1117FAX	20G0326FAX	20H0330FAX
		19B0426FAX	19H0731FAX	20A1023FAX	20C1118FAX	20G0327FAX	20H0869FAX
		19B0427FAX	19H0732FAX	20A1312FAX	20C1147FAX	20G0328FAX	20H0870FAX
		19B0428FAX	19H0733FAX	20A1313FAX	20C1148FAX	20G0329FAX	20H0871FAX
		19B0429FAX	19H0734FAX	20A1314FAX	20C1288FAX	20G0330FAX	20J0697FAX
		19C0529FAX	19K0025FAX	20A1365FAX	20C1289FAX	20G0331FAX	20J0698FAX
		19C0910FAX	19K0026FAX	20A1366FAX	20D1132FAX	20G0342FAX	20J0723FAX
		19D0031FAX	19K0656FAX	20A1367FAX	20D1133FAX	20G0497FAX	20J0724FAX
		19D0032FAX	19K1043FAX	20A1465FAX	20D1413FAX	20G0498FAX	20J0725FAX
		19D0033FAX	19K1044FAX	20B0154FAX	20D1414FAX	20G0499FAX	20J0983FAX
		19D0199FAX	19K1045FAX	20B0155FAX	20D1418FAX	20G0500FAX	20J0984FAX
		19D0200FAX	19K1336FAX	20B0156FAX	20E0251FAX	20G0887FAX	20J1033FAX
19D0201FAX	19K1337FAX	20B0266FAX	20F0248FAX	20G0888FAX	20J1034FAX		
19E0017FAX	19K1338FAX	20B0267FAX	20F0249FAX	20G0889FAX	20J1035FAX		
19E0019FAX	19L1198FAX	20B0268FAX	20F0251FAX	20G0890FAX	20J1036FAX		
19E0020FAX	19L1199FAX	20C0645FAX	20F0252FAX	20G0891FAX	20J1166FAX		
350/5865Z	Electrostatic Filter 22F X JAPAN	18M0095FAX	19E0596FAX	19F0490FAX	19L0926FAX	20E0715FAX	20E0716FAX
		20E0717FAX	20K0281FAX	20L0786FAX	20L0787FAX		
350/5865Z	Electrostatic Filter 22F X JAPAN	18D0305FAX	19C0380FAX	19H0781FAX	20B1023FAX	20F0217FAX	20K0455FAX
		18D0306FAX	19C0387FAX	19H0782FAX	20B1024FAX	20F0218FAX	20K0915FAX
		18J0040FAX	19C0388FAX	19H0964FAX	20B1234FAX	20F0456FAX	20K0916FAX
		18J0041FAX	19C0389FAX	19H0965FAX	20B1293FAX	20F0457FAX	20K0917FAX
		18J0045FAX	19C0549FAX	19H0966FAX	20B1294FAX	20F0458FAX	20K1030FAX
		18J0141FAX	19C0550FAX	19H0967FAX	20B1295FAX	20F0459FAX	20K1031FAX
		18J0157FAX	19C0551FAX	19J0081FAX	20B1296FAX	20F0733FAX	20K1032FAX
		18J0158FAX	19C0552FAX	19J0082FAX	20B1297FAX	20F0734FAX	20K1033FAX
		18J0355FAX	19C0554FAX	19J0083FAX	20C0047FAX	20F0735FAX	20K1034FAX
		18J0356FAX	19C0632FAX	19J0084FAX	20C0048FAX	20F0736FAX	20K1232FAX
		18J1166FAX	19C0977FAX	19J0250FAX	20C0049FAX	20F0737FAX	20K1233FAX
		18J1167FAX	19C0978FAX	19J0251FAX	20C0050FAX	20F0738FAX	20K1234FAX
		18K0553FAX	19C1067FAX	19J0253FAX	20C0063FAX	20F0739FAX	20K1235FAX
		18K0554FAX	19C1068FAX	19J0254FAX	20C0164FAX	20F0821FAX	20K1480FAX
		18K0600FAX	19C1069FAX	19J0722FAX	20C0165FAX	20F1009FAX	20K1481FAX
		18K0742FAX	19C1070FAX	19J1022FAX	20C0221FAX	20F1300FAX	20K1482FAX
		18K0743FAX	19C1434FAX	19J1023FAX	20C0222FAX	20G0054FAX	20K1483FAX

Item Code	Description	Affected Lot Numbers					
		18K0770FAX	19C1435FAX	19J1024FAX	20C0223FAX	20G0055FAX	20K1484FAX
		18K0771FAX	19C1436FAX	19J1176FAX	20C0224FAX	20G0056FAX	20K1485FAX
		18K1021FAX	19C1437FAX	19J1334FAX	20C0225FAX	20G0200FAX	20K1486FAX
		18K1022FAX	19D0021FAX	19J1335FAX	20C0226FAX	20G0201FAX	20L0282FAX
		18K1023FAX	19D0022FAX	19J1336FAX	20C0448FAX	20G0202FAX	20L0283FAX
		18K1024FAX	19D0190FAX	19K0201FAX	20C0449FAX	20G0203FAX	20L0284FAX
		18K1297FAX	19D0191FAX	19K0202FAX	20C0638FAX	20G0351FAX	20L0385FAX
		18K1298FAX	19D0370FAX	19K0388FAX	20C0639FAX	20G0352FAX	20L0386FAX
		18K1299FAX	19D0371FAX	19K0414FAX	20C0640FAX	20G0353FAX	20L0387FAX
		18K1341FAX	19D0372FAX	19K0415FAX	20C0641FAX	20G0354FAX	20L0495FAX
		18L0007FAX	19D0373FAX	19K0482FAX	20C0657FAX	20G0355FAX	20L0496FAX
		18L0008FAX	19D0469FAX	19K0483FAX	20C0658FAX	20G0775FAX	20L0497FAX
		18L0009FAX	19D0470FAX	19K0562FAX	20C0659FAX	20G0777FAX	20L0498FAX
		18L0215FAX	19D0644FAX	19K0563FAX	20C0660FAX	20G0778FAX	20L0541FAX
		18L0216FAX	19D0645FAX	19K0696FAX	20C0661FAX	20G0779FAX	20L0542FAX
		18L0217FAX	19D1068FAX	19K1046FAX	20C0662FAX	20G0865FAX	20L0543FAX
		18L0218FAX	19E0224FAX	19K1047FAX	20C0663FAX	20G0866FAX	20L0544FAX
		18L0376FAX	19E0225FAX	19K1048FAX	20C0745FAX	20G0867FAX	20L0832FAX
		18L0377FAX	19E0226FAX	19K1049FAX	20C0746FAX	20G0868FAX	20L0833FAX
		18L0378FAX	19E0227FAX	19K1107FAX	20C0747FAX	20G0869FAX	20L0834FAX
		18L0707FAX	19E0357FAX	19K1108FAX	20C0748FAX	20G0870FAX	20L0835FAX
		18L0708FAX	19E0358FAX	19K1109FAX	20C0749FAX	20G0871FAX	20L0836FAX
		18L0709FAX	19E0407FAX	19K1110FAX	20C0890FAX	20G1322FAX	20L0837FAX
		18L0710FAX	19E0408FAX	19L0219FAX	20C0891FAX	20G1323FAX	20L0838FAX
		18L0881FAX	19E0409FAX	19L0223FAX	20C0892FAX	20G1324FAX	20L0839FAX
		18L0882FAX	19E0410FAX	19L0224FAX	20C0893FAX	20G1325FAX	20L1094FAX
		18L0883FAX	19E0924FAX	19L0225FAX	20C0949FAX	20G1326FAX	20L1095FAX
		18L1033FAX	19E0925FAX	19L0372FAX	20C0950FAX	20G1475FAX	20L1096FAX
		18L1034FAX	19E0996FAX	19L0373FAX	20C0951FAX	20G1476FAX	20L1097FAX
		18L1035FAX	19E0997FAX	19L0374FAX	20C0952FAX	20G1477FAX	20L1098FAX
		18L1262FAX	19E0999FAX	19L0625FAX	20C0953FAX	20G1478FAX	20L1110FAX
		18L1263FAX	19E1109FAX	19L0626FAX	20C1127FAX	20G1479FAX	20L1111FAX
		18L1264FAX	19E1110FAX	19L0627FAX	20C1128FAX	20G1480FAX	20L1112FAX
		18L1265FAX	19E1111FAX	19L0628FAX	20C1129FAX	20G1549FAX	20L1328FAX
		18L1441FAX	19E1112FAX	19L0629FAX	20C1130FAX	20G1550FAX	20L1329FAX
		18L1442FAX	19E1136FAX	19L0872FAX	20C1131FAX	20G1551FAX	20L1330FAX
		18L1443FAX	19E1137FAX	19L0873FAX	20D0282FAX	20G1552FAX	20L1331FAX
		18L1444FAX	19E1439FAX	19L0874FAX	20D0283FAX	20G1553FAX	20L1332FAX
		18M0056FAX	19E1441FAX	19L0875FAX	20D0284FAX	20G1554FAX	20L1333FAX
		18M0057FAX	19E1443FAX	19L1131FAX	20D0285FAX	20G1555FAX	20L1386FAX
		18M0058FAX	19E1578FAX	19L1132FAX	20D0317FAX	20G1556FAX	20L1387FAX
		18M0103FAX	19E1579FAX	19L1133FAX	20D0318FAX	20H0053FAX	20L1388FAX
		18M0104FAX	19E1580FAX	19L1134FAX	20D0319FAX	20H0054FAX	20L1389FAX
		18M0327FAX	19E1686FAX	19L1176FAX	20D0320FAX	20H0055FAX	20M0069FAX
		18M0328FAX	19E1687FAX	19L1177FAX	20D0321FAX	20H0056FAX	20M0070FAX
		18M0329FAX	19F0151FAX	19L1178FAX	20D0322FAX	20H0057FAX	20M0072FAX
		18M0426FAX	19F0152FAX	19L1179FAX	20D0323FAX	20H0058FAX	20M0073FAX
		18M0427FAX	19F0153FAX	19L1377FAX	20D0693FAX	20H0350FAX	20M0074FAX
		18M0517FAX	19F0356FAX	19L1378FAX	20D0694FAX	20H0351FAX	20M0075FAX
		18M0518FAX	19F0358FAX	19L1379FAX	20D0695FAX	20H0352FAX	20M0464FAX
		18M0519FAX	19F0359FAX	19L1380FAX	20D0696FAX	20H0353FAX	20M0465FAX
		18M0575FAX	19F0492FAX	19L1381FAX	20D0769FAX	20H0354FAX	20M0466FAX
		18M0576FAX	19F0493FAX	19M0256FAX	20D0770FAX	20H0355FAX	20M0467FAX
		18M0577FAX	19F0494FAX	19M0257FAX	20D0771FAX	20H0504FAX	20M0468FAX
		18M0958FAX	19F0495FAX	19M0258FAX	20D0772FAX	20H0505FAX	20M0469FAX
		18M0959FAX	19F0496FAX	19M0259FAX	20D0773FAX	20H0506FAX	20M0470FAX
		18M0960FAX	19F0668FAX	19M0279FAX	20D0815FAX	20H0538FAX	20M0567FAX
		18M0965FAX	19F0669FAX	19M0280FAX	20D0816FAX	20H0541FAX	20M0568FAX
		18M0966FAX	19F0670FAX	19M0408FAX	20D0817FAX	20H0542FAX	20M0569FAX
		18M0967FAX	19F0671FAX	19M0409FAX	20D0826FAX	20H0768FAX	20M0570FAX
		19A0053FAX	19F0672FAX	19M0464FAX	20D0827FAX	20H0769FAX	20M0571FAX
		19A0054FAX	19F0673FAX	19M0465FAX	20D0828FAX	20H0770FAX	20M0892FAX
		19A0116FAX	19F0748FAX	19M0467FAX	20D0895FAX	20H0771FAX	20M0893FAX
		19A0117FAX	19F0749FAX	19M0620FAX	20D0896FAX	20H0853FAX	20M0941FAX
		19A0401FAX	19F0750FAX	19M0621FAX	20D0897FAX	20H0854FAX	20M0942FAX
		19A0402FAX	19F0872FAX	19M0920FAX	20D0898FAX	20H0855FAX	20M0985FAX
		19A0403FAX	19F0873FAX	19M0921FAX	20D0899FAX	20H0856FAX	20M0986FAX
		19A0404FAX	19F0874FAX	19M0922FAX	20D0900FAX	20H1080FAX	20M0987FAX
		19A0476FAX	19F0875FAX	19M0934FAX	20D1072FAX	20H1081FAX	20M0988FAX
		19A0477FAX	19F0907FAX	20A0159FAX	20D1073FAX	20H1082FAX	20M1234FAX
		19A0478FAX	19F0908FAX	20A0160FAX	20D1074FAX	20H1083FAX	20M1331FAX
		19A0609FAX	19F0909FAX	20A0161FAX	20D1075FAX	20H1141FAX	20M1448FAX
		19A0610FAX	19F0965FAX	20A0162FAX	20D1232FAX	20H1142FAX	20M1449FAX
		19A0611FAX	19F0966FAX	20A0379FAX	20D1233FAX	20H1143FAX	20M1450FAX
		19A0612FAX	19F0967FAX	20A0380FAX	20D1234FAX	20H1144FAX	20M1451FAX
		19A0678FAX	19F0968FAX	20A0621FAX	20D1235FAX	20H1145FAX	20M1452FAX
		19A0679FAX	19F1167FAX	20A0622FAX	20D1236FAX	20H1146FAX	20M1453FAX

350/5879 Electrostatic Filter

Item Code	Description	Affected Lot Numbers					
		19A0680FAX	19F1203FAX	20A0623FAX	20D1237FAX	20H1147FAX	20M1509FAX
		19A0986FAX	19F1204FAX	20A0624FAX	20D1238FAX	20J0205FAX	20M1510FAX
		19A0987FAX	19F1205FAX	20A0714FAX	20D1239FAX	20J0206FAX	20M1511FAX
		19A1146FAX	19F1206FAX	20A0715FAX	20D1597FAX	20J0207FAX	20M1512FAX
		19A1147FAX	19G0163FAX	20A0716FAX	20D1598FAX	20J0208FAX	20M1513FAX
		19A1148FAX	19G0164FAX	20A0970FAX	20D1599FAX	20J0466FAX	21A0200FAX
		19A1149FAX	19G0165FAX	20A1075FAX	20D1600FAX	20J0467FAX	21A0201FAX
		19A1150FAX	19G0166FAX	20A1076FAX	20E0005FAX	20J0468FAX	21A0202FAX
		19A1278FAX	19G0167FAX	20A1077FAX	20E0006FAX	20J0469FAX	21A0387FAX
		19A1279FAX	19G0168FAX	20A1078FAX	20E0007FAX	20J0638FAX	21A0388FAX
		19A1449FAX	19G0169FAX	20A1255FAX	20E0008FAX	20J0639FAX	21A0389FAX
		19A1450FAX	19G0170FAX	20A1256FAX	20E0009FAX	20J0640FAX	21A0390FAX
		19B0280FAX	19G0644FAX	20A1257FAX	20E0010FAX	20J0641FAX	21A0391FAX
		19B0281FAX	19G0645FAX	20A1407FAX	20E0011FAX	20J0936FAX	21A0412FAX
		19B0282FAX	19G0646FAX	20A1409FAX	20E0012FAX	20J0937FAX	21A0413FAX
		19B0460FAX	19G0647FAX	20B0149FAX	20E0173FAX	20J0939FAX	21A0414FAX
		19B0462FAX	19G0648FAX	20B0150FAX	20E0174FAX	20J0940FAX	21A0416FAX
		19B0463FAX	19G0680FAX	20B0151FAX	20E0175FAX	20J1066FAX	21A0417FAX
		19B0464FAX	19G0681FAX	20B0152FAX	20E0176FAX	20J1067FAX	21A0418FAX
		19B0546FAX	19G0682FAX	20B0153FAX	20E0179FAX	20J1068FAX	21A0419FAX
		19B0547FAX	19G0683FAX	20B0236FAX	20E0180FAX	20J1091FAX	21A0420FAX
		19B0703FAX	19G0770FAX	20B0237FAX	20E0658FAX	20J1300FAX	21A0421FAX
		19B0704FAX	19G0771FAX	20B0238FAX	20E0659FAX	20J1301FAX	21A0422FAX
		19B0705FAX	19G0772FAX	20B0445FAX	20E0660FAX	20J1342FAX	21A0423FAX
		19B0706FAX	19G0826FAX	20B0446FAX	20E0661FAX	20J1343FAX	21A0424FAX
		19B0707FAX	19G0827FAX	20B0447FAX	20E0662FAX	20J1344FAX	21A0425FAX
		19B0708FAX	19G0828FAX	20B0448FAX	20E0718FAX	20J1345FAX	21A0426FAX
		19B0975FAX	19G0837FAX	20B0449FAX	20E0810FAX	20J1346FAX	21A0427FAX
		19B0976FAX	19G0838FAX	20B0582FAX	20E0920FAX	20K0029FAX	21A0571FAX
		19B0977FAX	19G0839FAX	20B0583FAX	20E0921FAX	20K0030FAX	21A0572FAX
		19B0978FAX	19G0923FAX	20B0584FAX	20E0922FAX	20K0031FAX	21A0573FAX
		19B1129FAX	19G1109FAX	20B0585FAX	20E1004FAX	20K0032FAX	21A0574FAX
		19B1130FAX	19G1110FAX	20B0630FAX	20E1005FAX	20K0033FAX	21A0575FAX
		19B1131FAX	19G1203FAX	20B0631FAX	20E1006FAX	20K0034FAX	21A0576FAX
		19B1381FAX	19G1204FAX	20B0632FAX	20E1058FAX	20K0246FAX	21A0791FAX
		19B1382FAX	19G1205FAX	20B0690FAX	20E1059FAX	20K0247FAX	21A0792FAX
		19B1383FAX	19G1206FAX	20B0691FAX	20E1060FAX	20K0248FAX	21A0795FAX
		19B1384FAX	19H0685FAX	20B0692FAX	20E1061FAX	20K0249FAX	21A0796FAX
		19B1516FAX	19H0686FAX	20B0693FAX	20E1062FAX	20K0448FAX	21A1111FAX
		19B1517FAX	19H0687FAX	20B0864FAX	20F0030FAX	20K0450FAX	21A1112FAX
		19B1518FAX	19H0688FAX	20B0865FAX	20F0032FAX	20K0451FAX	21B0250FAX
		19B1519FAX	19H0689FAX	20B0992FAX	20F0033FAX	20K0452FAX	
		19B1520FAX	19H0779FAX	20B0993FAX	20F0186FAX	20K0453FAX	
		19C0379FAX	19H0780FAX	20B0994FAX	20F0188FAX	20K0454FAX	
350/5879TC	Electrostatic Filter	19A1263FAX	20D1103FAX	20H0617FAX	20J0811FAX	20K0077FAX	20M0386FAX
		19G1306FAX	20G0117FAX	20J0098FAX	20J1126FAX	20K0459FAX	21A0680FAX
		19J1114FAX	20G1359FAX	20J0099FAX	20K0015FAX	20K0460FAX	21A0681FAX
		19M0091FAX	20G1582FAX	20J0810FAX	20K0076FAX	20M0385FAX	
350/5879Z	Electrostatic Filter	18M0591FAX	19G0897FAX	20D0862FAX	20D0863FAX	20G0640FAX	20G0641FAX
		20H0716FAX					
350/5882	Electrostatic Filter with Catheter Mount	20J0204FAX					
350S19006	Electrostatic Filter Angled	18M0322FAX	19C0530FAX	19G0134FAX	20C0005FAX	20F0837FAX	20L0411FAX
		18M0323FAX	19D0474FAX	19G0992FAX	20C0006FAX	20G1302FAX	20L0412FAX
		18M0598FAX	19D0475FAX	19K0185FAX	20D0108FAX	20G1303FAX	20M0167FAX
		19A0918FAX	19D0476FAX	19L0414FAX	20D0109FAX	20G1304FAX	20M0168FAX
		19A0919FAX	19E0274FAX	19L0415FAX	20D0968FAX	20G1305FAX	20M0169FAX
		19A0920FAX	19E0275FAX	19L0416FAX	20E1096FAX	20L0026FAX	20M0170FAX
		19A0923FAX	19G0130FAX	20A0324FAX	20F0835FAX	20L0027FAX	20M0531FAX
		19A0924FAX	19G0131FAX	20C0004FAX	20F0836FAX	20L0410FAX	20M1337FAX
		19A0925FAX	19G0133FAX				
350S5807	Electrostatic Filter	18L0194FAX	18L1252FAX	19D0294FAX	19E1691FAX	19G0097FAX	19H0884FAX
		19J0636FAX	19K1219FAX	20A0478FAX	20D1102FAX	20H0163FAX	20J0050FAX
351/19002	Mechanical Filter with Catheter Mount	20H0564FAX	20L0440FAX				
351/5410	Mechanical Filter	18K0757FAX	19C0759FAX	19J0485FAX	20C0058FAX	20E0903FAX	20L1023FAX
		18K0758FAX	19C0760FAX	19J0631FAX	20C0059FAX	20F0890FAX	20M0617FAX
		18K0759FAX	19D0171FAX	19K0500FAX	20C0060FAX	20F0891FAX	20M0618FAX
		18K0760FAX	19D0172FAX	19K0501FAX	20C0692FAX	20F0892FAX	20M0619FAX
		18K1333FAX	19D0173FAX	19K0502FAX	20C0693FAX	20F1232FAX	20M0620FAX
		18K1334FAX	19D0174FAX	19K0605FAX	20C0694FAX	20F1233FAX	20M0681FAX
		18K1335FAX	19D0175FAX	19K0720FAX	20C0708FAX	20F1234FAX	21A0215FAX
		18L0118FAX	19D0176FAX	19K0731FAX	20C0709FAX	20G0199FAX	21A0216FAX
		18L0119FAX	19D0177FAX	19L0534FAX	20C0757FAX	20G1184FAX	21A0217FAX
		18L0120FAX	19D0542FAX	19L0538FAX	20C0758FAX	20G1185FAX	21A0218FAX
		18L0121FAX	19E0973FAX	19L0539FAX	20C0759FAX	20H0602FAX	21A0219FAX
		18L1136FAX	19F0807FAX	19L0683FAX	20C0760FAX	20H0603FAX	21A0255FAX
		18L1137FAX	19F0808FAX	19L0904FAX	20C1158FAX	20H0799FAX	21A0256FAX
		18L1138FAX	19F0809FAX	19L0920FAX	20C1159FAX	20H0800FAX	21A0257FAX

Item Code	Description	Affected Lot Numbers							
351/5410	Mechanical Filter	18L1139FAX	19F0810FAX	19L0921FAX	20C1160FAX	20J1220FAX	21A0258FAX		
		18L1283FAX	19G1041FAX	19L0922FAX	20C1161FAX	20J1221FAX	21A0469FAX		
		18L1284FAX	19G1376FAX	19M0903FAX	20C1162FAX	20J1222FAX	21A0470FAX		
		18L1285FAX	19G1377FAX	19M0904FAX	20C1163FAX	20K0500FAX	21A0471FAX		
		19B0066FAX	19G1378FAX	19M0905FAX	20C1286FAX	20K0501FAX	21A1185FAX		
		19B0067FAX	19J0283FAX	20B0486FAX	20D1400FAX	20K0502FAX	21A1186FAX		
		19B0068FAX	19J0284FAX	20B0487FAX	20D1401FAX	20K0503FAX	21A1189FAX		
		19B0069FAX	19J0392FAX	20B0488FAX	20D1402FAX	20K0504FAX	21B0017FAX		
		19B0070FAX	19J0438FAX	20B0489FAX	20D1403FAX	20K0505FAX	21B0018FAX		
		19B0071FAX	19J0439FAX	20B0638FAX	20E0090FAX	20K0506FAX	21B0019FAX		
		19B0072FAX	19J0440FAX	20B0639FAX	20E0091FAX	20K0507FAX	21B0020FAX		
		19B0740FAX	19J0482FAX	20B0640FAX	20E0092FAX	20K0853FAX	21B0021FAX		
		19C0757FAX	19J0483FAX	20B0641FAX	20E0093FAX	20L1021FAX			
		19C0758FAX	19J0484FAX	20B0702FAX	20E0094FAX	20L1022FAX			
		351/5410TC	Mechanical Filter	18L1120FAX	19B0073FAX	20D1404FAX	20H0604FAX	20K0013FAX	21A0569FAX
		351/5410Z	Mechanical Filter X JAPAN	18L1140FAX	19F0811FAX	19L0684FAX	20D1399FAX	20K0499FAX	20L1025FAX
				19C0761FAX	19K0604FAX	20B0490FAX			
351/5834	Mechanical Filter 'S'	18L0826FAX	18M0275FAX	19B1503FAX	19G0625FAX	19G1486FAX	19H0783FAX		
		19J1074FAX	19K0416FAX	19L0045FAX	20A1188FAX	20B1306FAX	20E0376FAX		
		20E0480FAX	20E0481FAX	20E0482FAX	20K1051FAX	20L0219FAX	20M0554FAX		
351/5848	Mechanical Filter with Catheter Mount	21A0779FAX							
351/5856	Mechanical Filter 22M-22F	18K1195FAX	19E1077FAX	20B0784FAX	20D1086FAX	20H0452FAX	20K0107FAX		
		18K1196FAX	19E1593FAX	20B0785FAX	20F0560FAX	20J1302FAX	20K0108FAX		
		18L1084FAX	19E1594FAX	20B0786FAX	20F0561FAX	20J1303FAX	20K0109FAX		
		18L1085FAX	19E1595FAX	20B1400FAX	20F0562FAX	20J1304FAX	20K0110FAX		
		18M0017FAX	19K0885FAX	20C0094FAX	20F0563FAX	20J1305FAX	20K0111FAX		
		19B0075FAX	19K1000FAX	20C0509FAX	20G0858FAX	20J1306FAX	20K0222FAX		
		19C1063FAX	19L0532FAX	20C0510FAX	20G0859FAX	20J1307FAX	20L1273FAX		
		19C1064FAX	19L0924FAX	20C0511FAX	20G1216FAX	20J1350FAX	20M0682FAX		
		19C1253FAX	19L1023FAX	20C0572FAX	20H0428FAX	20J1351FAX	21A0214FAX		
		19C1254FAX	19M0524FAX	20D1084FAX	20H0429FAX	20J1352FAX			
		19E1076FAX	19M0753FAX	20D1085FAX	20H0451FAX	20J1353FAX			
351/5856Z	Mechanical Filter 22M-22F	18L0117FAX	19B0074FAX	19C0536FAX	20L1275FAX	20L1276FAX	20L1277FAX		
		20L1278FAX							
351/5878	Mechanical Filter	18L0586FAX	19A0074FAX	19K0341FAX	20A1192FAX	20F1230FAX	20K1049FAX		
		18L0587FAX	19A1035FAX	19K0342FAX	20A1193FAX	20F1231FAX	20K1050FAX		
		18L0588FAX	19A1036FAX	19K0343FAX	20A1194FAX	20G1664FAX	20L0223FAX		
		18L0589FAX	19A1037FAX	19K1071FAX	20A1388FAX	20G1665FAX	20L0224FAX		
		18M0277FAX	19B1504FAX	19L0048FAX	20A1389FAX	20G1666FAX	20L0225FAX		
		18M0278FAX	19C0328FAX	19L0049FAX	20A1390FAX	20G1667FAX	20L0226FAX		
		18M0279FAX	19C0533FAX	19L0050FAX	20A1482FAX	20H0515FAX	20L0951FAX		
		18M0280FAX	19D0685FAX	19L0051FAX	20B1311FAX	20J0506FAX	20M0812FAX		
		18M0407FAX	19E0053FAX	19L0052FAX	20B1312FAX	20J0507FAX	20M0813FAX		
		18M0408FAX	19E0055FAX	19L0420FAX	20D0078FAX	20J0508FAX	20M0853FAX		
		18M0409FAX	19F0204FAX	19L0421FAX	20D0079FAX	20J0509FAX	20M1514FAX		
		18M0873FAX	19F0674FAX	19M0166FAX	20E1308FAX	20J1018FAX	20M1515FAX		
		18M0874FAX	19G0411FAX	19M0167FAX	20E1309FAX	20K1047FAX	20M1516FAX		
		18M0875FAX	19J1264FAX	20A0340FAX	20E1310FAX	20K1048FAX			
		351/5878TC	Mechanical Filter	20F0002FAX	20L0398FAX				
		351/5878Z	Mechanical Filter	18L0585FAX	19A1034FAX	19K0345FAX	20A1191FAX	20E1311FAX	20K1046FAX
18L1151FAX	19C0338FAX			19L0046FAX	20D0071FAX	20G1668FAX	20L0222FAX		
18M0276FAX	19F0205FAX			20A1190FAX	20E0377FAX	20J0480FAX	21A0780FAX		
19A0015FAX									
351/5969	Mechanical Filter with Connector	20K0835FAX							
351/5970	Mechanical Filter	18K1120FAX	19C1221FAX	19J0244FAX	20C0780FAX	20E0996FAX	20H0866FAX		
		18K1121FAX	19C1222FAX	19J0245FAX	20C0781FAX	20E1338FAX	20H0867FAX		
		18K1122FAX	19C1542FAX	19J0246FAX	20C0799FAX	20E1339FAX	20H0868FAX		
		18K1358FAX	19C1543FAX	19J0247FAX	20C0800FAX	20E1340FAX	20I0699FAX		
		18K1359FAX	19C1544FAX	19J0248FAX	20C0801FAX	20F0042FAX	20J0700FAX		
		18K1360FAX	19D0034FAX	19J0249FAX	20C0802FAX	20F0163FAX	20J0701FAX		
		18K1361FAX	19D0035FAX	19J1152FAX	20C1199FAX	20F0164FAX	20J0702FAX		
		18L0115FAX	19D0036FAX	19J1153FAX	20C1200FAX	20F0219FAX	20J0800FAX		
		18L0116FAX	19D0037FAX	19J1338FAX	20C1201FAX	20F0220FAX	20J0801FAX		
		18L0711FAX	19D0554FAX	19J1339FAX	20C1202FAX	20F0223FAX	20J0802FAX		
		18L0712FAX	19D0555FAX	19K0120FAX	20C1203FAX	20F0545FAX	20J0803FAX		
		18L0713FAX	19D0556FAX	19K0122FAX	20D0146FAX	20F0546FAX	20J0943FAX		
		18L0714FAX	19D1028FAX	19K0123FAX	20D0147FAX	20F1027FAX	20J0944FAX		
		18L0715FAX	19D1071FAX	19K1329FAX	20D0148FAX	20G0021FAX	20J1132FAX		
		18L1092FAX	19D1072FAX	19K1330FAX	20D0149FAX	20G0022FAX	20J1133FAX		
		18L1093FAX	19D1073FAX	19K1331FAX	20D0150FAX	20G0023FAX	20J1134FAX		
		18L1094FAX	19D1074FAX	19L0091FAX	20D0324FAX	20G0138FAX	20J1135FAX		
		18L1095FAX	19E0475FAX	19L0092FAX	20D0325FAX	20G0139FAX	20J1136FAX		
		18L1096FAX	19E0549FAX	19L0093FAX	20D0819FAX	20G0204FAX	20J1137FAX		
		18M0270FAX	19E0550FAX	19L0094FAX	20D0820FAX	20G0205FAX	20J1401FAX		
		18M0271FAX	19E0551FAX	19L0095FAX	20D0999FAX	20G0206FAX	20J1402FAX		
		18M0272FAX	19E0712FAX	19L1215FAX	20D1001FAX	20G0371FAX	20J1403FAX		

Item Code	Description	Affected Lot Numbers					
351/5979	Mechanical Filter MINI	18M0380FAX	19E1138FAX	19L1216FAX	20D1002FAX	20G0372FAX	20J1404FAX
		18M0387FAX	19E1410FAX	19L1217FAX	20D1253FAX	20G0373FAX	20J1405FAX
		18M0388FAX	19E1411FAX	19L1218FAX	20D1254FAX	20G0491FAX	20K0404FAX
		18M0389FAX	19E1413FAX	19M0208FAX	20D1255FAX	20G0492FAX	20K0405FAX
		18M0544FAX	19E1452FAX	19M0209FAX	20D1256FAX	20G0493FAX	20K0406FAX
		18M0545FAX	19E1453FAX	19M0263FAX	20D1421FAX	20G0494FAX	20K0407FAX
		18M0546FAX	19E1454FAX	19M0264FAX	20D1422FAX	20G0495FAX	20K0408FAX
		18M0547FAX	19F0034FAX	19M0265FAX	20D1529FAX	20G0496FAX	20K0409FAX
		18M0548FAX	19F0038FAX	19M0716FAX	20D1673FAX	20G0862FAX	20K0669FAX
		18M0807FAX	19G0189FAX	19M0896FAX	20D1674FAX	20G0863FAX	20K0670FAX
		18M0808FAX	19G0190FAX	19M0897FAX	20D1675FAX	20G0864FAX	20K0671FAX
		18M1009FAX	19G0191FAX	19M0898FAX	20E0401FAX	20H0060FAX	20K0672FAX
		19A0692FAX	19G0192FAX	19M0899FAX	20E0402FAX	20H0061FAX	20K0673FAX
		19B0027FAX	19G0193FAX	20A0115FAX	20E0403FAX	20H0071FAX	20K0674FAX
		19B0029FAX	19G0194FAX	20A0265FAX	20E0404FAX	20H0072FAX	20K1458FAX
		19B0794FAX	19G0676FAX	20A1495FAX	20E0405FAX	20H0324FAX	20K1459FAX
		19B1022FAX	19G0677FAX	20A1496FAX	20E0561FAX	20H0325FAX	20K1460FAX
		19B1044FAX	19G0678FAX	20B0160FAX	20E0562FAX	20H0326FAX	20K1461FAX
		19B1435FAX	19G1001FAX	20B0161FAX	20E0617FAX	20H0327FAX	20K1462FAX
		19B1436FAX	19G1002FAX	20B0162FAX	20E0618FAX	20H0590FAX	20M1398FAX
		19B1437FAX	19H0581FAX	20B0163FAX	20E0619FAX	20H0593FAX	20M1399FAX
		19B1506FAX	19H0753FAX	20B0791FAX	20E0720FAX	20H0780FAX	20M1400FAX
		19B1531FAX	19H0754FAX	20B0792FAX	20E0721FAX	20H0781FAX	20M1401FAX
		19B1532FAX	19H0755FAX	20B0793FAX	20E0815FAX	20H0782FAX	20M1498FAX
		19B1533FAX	19H0756FAX	20B0947FAX	20E0816FAX	20H0783FAX	20M1499FAX
		19B1534FAX	19H0757FAX	20B0948FAX	20E0817FAX	20H0861FAX	20M1500FAX
		19C0951FAX	19H0951FAX	20C0501FAX	20E0992FAX	20H0862FAX	
		19C0952FAX	19H0952FAX	20C0507FAX	20E0993FAX	20H0863FAX	
		19C1219FAX	19H0953FAX	20C0508FAX	20E0994FAX	20H0864FAX	
		19C1220FAX	19H0954FAX	20C0564FAX	20E0995FAX	20H0865FAX	
		351/5979TC	Mechanical Filter with Connector Small, Tethered Cap	19E0379FAX	19E1360FAX	19F1051FAX	19G1004FAX
		19M0713FAX	20D1262FAX	20G1357FAX	20G1358FAX	20K0014FAX	
351/5980	Mechanical Filter with Catheter Mount	19L1213FAX					
351/5984	Mechanical Filter with Connector	18M0984FAX	19E0719FAX	19J0052FAX	20D1596FAX	20H0067FAX	20H0070FAX
		19B1021FAX	19G0199FAX	20A1460FAX	20G0860FAX	20H0068FAX	20J1223FAX
		19B1434FAX	19G0994FAX	20D1248FAX	20G0861FAX	20H0069FAX	20J1225FAX
		19E0718FAX	19G0995FAX	20D1250FAX	20H0066FAX		
351S19032	Mechanical Filter with Catheter Mount	20A0566FAX	20J0530FAX	20J0531FAX			
352/5802	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	19B0661FAX					
352/5805	Adult-Pediatric Electrostatic Filter HME	18L0628FAX	19G0365FAX	20A0703FAX	20D0523FAX	20E1012FAX	20H1089FAX
		18L0629FAX	19G0366FAX	20A0704FAX	20D0524FAX	20E1080FAX	20H1095FAX
		18L0630FAX	19G0367FAX	20A0774FAX	20D0525FAX	20E1081FAX	20H1096FAX
		18L0631FAX	19G0368FAX	20A0775FAX	20D0526FAX	20E1082FAX	20H1097FAX
		18M0774FAX	19G0369FAX	20A0839FAX	20D0665FAX	20E1257FAX	20H1098FAX
		18M0775FAX	19G0370FAX	20A0840FAX	20D0666FAX	20E1258FAX	20H1099FAX
		18M0776FAX	19G0371FAX	20A1020FAX	20D0667FAX	20E1259FAX	20J1044FAX
		18M0777FAX	19G0372FAX	20B0460FAX	20D0668FAX	20E1397FAX	20J1045FAX
		19A0973FAX	19G0373FAX	20B0461FAX	20D0680FAX	20E1398FAX	20J1046FAX
		19B1386FAX	19J0263FAX	20B0462FAX	20E0169FAX	20E1417FAX	20J1227FAX
		19D0182FAX	19J0264FAX	20B0463FAX	20E0170FAX	20E1418FAX	20J1228FAX
		19D0183FAX	19J0265FAX	20B0464FAX	20E0171FAX	20F0822FAX	20J1229FAX
		19D0184FAX	19J0354FAX	20B0465FAX	20E0172FAX	20F0823FAX	20K0089FAX
		19D0186FAX	19J0355FAX	20B0571FAX	20E0259FAX	20F0824FAX	20K0090FAX
		19D0187FAX	19J0356FAX	20B0572FAX	20E0260FAX	20F0825FAX	20K1490FAX
		19D0543FAX	19J0408FAX	20B0626FAX	20E0261FAX	20F0826FAX	20K1491FAX
		19D0544FAX	19J0409FAX	20B0627FAX	20E0262FAX	20F0827FAX	20K1493FAX
		19D0545FAX	19J0410FAX	20B0700FAX	20E0263FAX	20F1048FAX	20K1494FAX
		19D0546FAX	19J0575FAX	20B0715FAX	20E0264FAX	20F1049FAX	20K1495FAX
		19D0547FAX	19K1008FAX	20C0664FAX	20E0265FAX	20F1052FAX	20K1496FAX
		19D0677FAX	19K1341FAX	20C0665FAX	20E0475FAX	20F1053FAX	20L0358FAX
		19D1333FAX	19K1342FAX	20C0666FAX	20E0476FAX	20F1114FAX	20L0359FAX
		19D1335FAX	19K1343FAX	20C0667FAX	20E0477FAX	20F1115FAX	20L0360FAX
		19D1336FAX	19L0067FAX	20C0668FAX	20E0478FAX	20F1116FAX	20L0361FAX
		19D1337FAX	19L0109FAX	20C0995FAX	20E0479FAX	20F1117FAX	20L0362FAX
		19E0359FAX	19L0110FAX	20D0306FAX	20E0669FAX	20F1118FAX	20L1308FAX
		19E0360FAX	19L0111FAX	20D0307FAX	20E0670FAX	20F1119FAX	20L1309FAX
		19E0362FAX	19L0370FAX	20D0326FAX	20E0671FAX	20G1380FAX	20L1310FAX
		19E0476FAX	19L0371FAX	20D0327FAX	20E0999FAX	20G1381FAX	20M0081FAX
		19E0683FAX	19L0527FAX	20D0408FAX	20E1000FAX	20H1085FAX	20M0082FAX
		19E0684FAX	19L0528FAX	20D0409FAX	20E1001FAX	20H1086FAX	20M1376FAX
19E1040FAX	19L0529FAX	20D0410FAX	20E1002FAX	20H1087FAX	20M1377FAX		
19E1041FAX	19L1330FAX	20D0522FAX	20E1003FAX	20H1088FAX			
352/5805TC	Adult-Pediatric Electrostatic Filter HME	20E0907FAX	20E0908FAX	20G0886FAX			
352/5944	Adult-Pediatric Electrostatic Filter HME	18L0419FAX	19D0285FAX	19E0791FAX	20B1084FAX	20D1533FAX	20K0586FAX
		18L0420FAX	19D1123FAX	19H0702FAX	20D0311FAX	20D1534FAX	20K1139FAX
		18L1259FAX	19D1124FAX	19J1326FAX	20D0312FAX	20D1535FAX	20L0407FAX

Item Code	Description	Affected Lot Numbers					
352/5844	Adult-Pediatric Electrostatic Filter HME	19A0289FAX	19E0265FAX	19J1327FAX	20D1531FAX	20G1111FAX	20M1385FAX
		19A1268FAX	19E0266FAX	19M0770FAX	20D1532FAX	20G1112FAX	20M1386FAX
		19D0284FAX					
352/5844Z	Adult-Pediatric Electrostatic Filter HME	18L1260FAX	19E0264FAX	19J1328FAX	19M0768FAX	20E0802FAX	20F0717FAX
		19A0290FAX	19E1495FAX	19J1329FAX	20B1026FAX	20E0803FAX	20F0718FAX
		19A1269FAX	19E1496FAX	19L0105FAX	20B1027FAX	20E0937FAX	20K0587FAX
		19B0382FAX	19E1497FAX	19L0106FAX	20B1028FAX	20E0938FAX	20K0588FAX
		19D0287FAX	19F0162FAX	19L0107FAX	20C0138FAX	20E0939FAX	20K0589FAX
		19D0288FAX	19F0163FAX	19L1102FAX	20E0603FAX	20E1011FAX	20K0590FAX
		19D0459FAX	19F0164FAX	19L1241FAX	20E0604FAX	20F0713FAX	20K0591FAX
		19D0460FAX	19F0165FAX	19L1242FAX	20E0607FAX	20F0714FAX	20K0592FAX
		19D0461FAX	19H0703FAX	19L1243FAX	20E0609FAX	20F0715FAX	20L0519FAX
		19D0462FAX	19H0704FAX	19M0123FAX	20E0800FAX	20F0716FAX	20L0520FAX
		19D0463FAX	19H0705FAX	19M0124FAX	20E0801FAX		
352/5855	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	17E1489FAX	20A0996FAX	20D0121FAX	20D0521FAX	20J0850FAX	20K0522FAX
		18E1564FAX	20A1026FAX	20D0122FAX	20D0660FAX	20J1098FAX	20K0523FAX
		18E1565FAX	20A1027FAX	20D0151FAX	20D0661FAX	20J1102FAX	20M0171FAX
		18E1566FAX	20A1351FAX	20D0152FAX	20D0662FAX	20J1103FAX	20M0172FAX
		19L0178FAX	20A1352FAX	20D0153FAX	20D1289FAX	20J1202FAX	20M0173FAX
		19L0179FAX	20A1353FAX	20D0428FAX	20D1290FAX	20J1203FAX	20M1039FAX
		19L0180FAX	20A1354FAX	20D0429FAX	20D1291FAX	20J1204FAX	20M1040FAX
		20A0473FAX	20C0682FAX	20D0430FAX	20D1292FAX	20J1205FAX	20M1387FAX
		20A0707FAX	20C0683FAX	20D0431FAX	20G1630FAX	20K0041FAX	20M1388FAX
		20A0971FAX	20D0117FAX	20D0517FAX	20G1631FAX	20K0042FAX	20M1389FAX
		20A0972FAX	20D0118FAX	20D0518FAX	20H0186FAX	20K0043FAX	20M1390FAX
		20A0994FAX	20D0119FAX	20D0519FAX	20J0848FAX	20K0044FAX	20M1391FAX
		20A0995FAX	20D0120FAX	20D0520FAX	20J0849FAX	20K0045FAX	
		352/5863	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	19L1355FAX	20D0528FAX	20J1271FAX	20K0083FAX
352/5867	Adult-Pediatric Electrostatic Filter HME with Connector	18L1245FAX	19D0195FAX	19H0750FAX	20A1204FAX	20E0101FAX	20K0542FAX
		18L1246FAX	19D0196FAX	19H0900FAX	20A1205FAX	20E0102FAX	20K0543FAX
		18L1247FAX	19D0197FAX	19H0901FAX	20A1206FAX	20E0384FAX	20K0997FAX
		18L1248FAX	19D0198FAX	19H0902FAX	20A1218FAX	20E0385FAX	20K0998FAX
		18L1466FAX	19D0253FAX	19H0903FAX	20B0907FAX	20E0386FAX	20K0999FAX
		18L1467FAX	19D0254FAX	19H0904FAX	20B0908FAX	20F0366FAX	20K1002FAX
		18M0937FAX	19E0251FAX	19J0604FAX	20B0909FAX	20F0589FAX	20K1004FAX
		18M0938FAX	19E0252FAX	19J1172FAX	20B1098FAX	20F0992FAX	20K1094FAX
		19A0566FAX	19E0255FAX	19J1173FAX	20B1099FAX	20G0484FAX	20K1514FAX
		19A0567FAX	19E0256FAX	19K0136FAX	20B1100FAX	20G0486FAX	20K1515FAX
		19A1243FAX	19E1491FAX	19K0373FAX	20B1276FAX	20G0487FAX	20K1516FAX
		19B1359FAX	19E1493FAX	19K0374FAX	20B1277FAX	20G0488FAX	20K1517FAX
		19B1360FAX	19E1494FAX	19K0389FAX	20B1402FAX	20G0489FAX	20K1518FAX
		19B1361FAX	19F1224FAX	19K0411FAX	20B1403FAX	20J0279FAX	20K1519FAX
		19B1362FAX	19F1225FAX	19K1117FAX	20B1404FAX	20J0280FAX	20L0330FAX
		19C0122FAX	19F1262FAX	19K1118FAX	20C0259FAX	20J0281FAX	20L0331FAX
		19C0123FAX	19F1263FAX	19K1282FAX	20C1157FAX	20J0282FAX	20L0332FAX
		19C0124FAX	19F1264FAX	19K1283FAX	20C1190FAX	20J0283FAX	20L0334FAX
		19C0125FAX	19G0173FAX	19K1284FAX	20D0252FAX	20J0637FAX	20L0355FAX
		19C0390FAX	19G0174FAX	19L1143FAX	20D1305FAX	20J1253FAX	20L0356FAX
		19C0391FAX	19G0175FAX	19L1144FAX	20D1325FAX	20K0270FAX	20L0357FAX
		19C0392FAX	19G0176FAX	19L1145FAX	20D1326FAX	20K0271FAX	20M1383FAX
		19C0393FAX	19G0177FAX	19M0365FAX	20D1327FAX	20K0272FAX	21A1173FAX
		19C1208FAX	19G0178FAX	19M0457FAX	20D1328FAX	20K0273FAX	21A1174FAX
		19C1438FAX	19H0678FAX	19M0458FAX	20E0095FAX	20K0274FAX	21A1175FAX
		19C1439FAX	19H0679FAX	19M0459FAX	20E0096FAX	20K0275FAX	21A1176FAX
		19C1440FAX	19H0680FAX	19M0460FAX	20E0097FAX	20K0538FAX	21A1177FAX
		19C1441FAX	19H0747FAX	20A0845FAX	20E0098FAX	20K0539FAX	21A1178FAX
		19D0038FAX	19H0748FAX	20A0992FAX	20E0099FAX	20K0540FAX	
		19D0039FAX	19H0749FAX	20A0993FAX	20E0100FAX	20K0541FAX	
		352/5867TC	Adult-Pediatric Electrostatic Filter HME with Connector	18L1124FAX	19A0312FAX	19A1244FAX	19F0168FAX
19M0782FAX	20C0303FAX			20D0245FAX	20G1222FAX	20G1578FAX	20H0367FAX
20J0284FAX	20J1165FAX			20M0656FAX	20M1384FAX		
352/5867Z	Adult-Pediatric Electrostatic Filter HME with Connector	18L1261FAX	19E0554FAX	19J0318FAX	19M0364FAX	20F0885FAX	20K0871FAX
		18M0592FAX	19E0584FAX	19K0210FAX	20B0123FAX	20F0886FAX	20K0872FAX
		18M0593FAX	19F0157FAX	19K0211FAX	20B0124FAX	20F0887FAX	20L0413FAX
		19A0310FAX	19F0484FAX	19K0212FAX	20B0355FAX	20F0888FAX	20L0414FAX
		19A0311FAX	19F0485FAX	19K0390FAX	20B0356FAX	20F0889FAX	20L0415FAX
		19B1368FAX	19F0615FAX	19K0391FAX	20B0358FAX	20F1149FAX	20L1019FAX
		19B1369FAX	19F0616FAX	19K0417FAX	20B0359FAX	20F1150FAX	20L1020FAX
		19B1370FAX	19F0650FAX	19K0418FAX	20B0524FAX	20G0721FAX	20L1092FAX
		19B1371FAX	19F0651FAX	19K0419FAX	20B0936FAX	20G0722FAX	20L1093FAX
		19B1372FAX	19F0652FAX	19K0420FAX	20C0159FAX	20G0723FAX	20L1258FAX
		19C0137FAX	19F0755FAX	19K0569FAX	20C0160FAX	20G0724FAX	20L1259FAX
		19C0138FAX	19F0756FAX	19K0570FAX	20C0161FAX	20H0034FAX	20L1260FAX
		19C0139FAX	19F0757FAX	19K0627FAX	20C0304FAX	20H0172FAX	20L1261FAX
		19C0140FAX	19F1220FAX	19K0628FAX	20C0305FAX	20H0173FAX	20L1262FAX
		19D0423FAX	19F1221FAX	19K1406FAX	20C0306FAX	20J0310FAX	20L1263FAX
		19D0600FAX	19F1222FAX	19K1407FAX	20C0307FAX	20J0311FAX	20L1348FAX

Item Code	Description	Affected Lot Numbers					
19D0601FAX	19F1223FAX	19K1408FAX	20C0308FAX	20J0312FAX	20L1349FAX		
19D0602FAX	19F1258FAX	19L0043FAX	20C0309FAX	20J0313FAX	20L1381FAX		
19D0674FAX	19F1259FAX	19L0082FAX	20D1192FAX	20J1073FAX	20L1382FAX		
19D0675FAX	19F1260FAX	19L0085FAX	20D1193FAX	20J1074FAX	21A1179FAX		
19D0676FAX	19G1009FAX	19L0086FAX	20D1194FAX	20J1075FAX	21A1180FAX		
19D1255FAX	19G1010FAX	19L0616FAX	20D1195FAX	20K0626FAX	21A1181FAX		
19D1256FAX	19G1012FAX	19L1293FAX	20D1196FAX	20K0627FAX	21A1182FAX		
19D1327FAX	19G1013FAX	19L1294FAX	20D1197FAX	20K0628FAX	21A1183FAX		
19D1328FAX	19J0172FAX	19L1295FAX	20D1543FAX	20K0629FAX	21A1184FAX		
19D1329FAX	19J0173FAX	19L1296FAX	20E0580FAX	20K0849FAX			
19D1330FAX	19J0174FAX	19M0199FAX	20E0582FAX	20K0850FAX			
17E1488FAX	19B1128FAX	19J0234FAX	20B1103FAX	20F0902FAX	20J1131FAX		
18D0302FAX	19B1316FAX	19J0235FAX	20B1104FAX	20F0903FAX	20J1195FAX		
18D0303FAX	19B1317FAX	19J0271FAX	20B1231FAX	20F0904FAX	20J1196FAX		
18D0304FAX	19B1319FAX	19J0272FAX	20B1232FAX	20F0905FAX	20J1197FAX		
18J0042FAX	19B1320FAX	19J0273FAX	20B1233FAX	20F0906FAX	20J1198FAX		
18J0043FAX	19B1321FAX	19J0274FAX	20B1414FAX	20F0907FAX	20J1199FAX		
18J0044FAX	19B1322FAX	19J0275FAX	20B1415FAX	20F1003FAX	20J1200FAX		
18J0046FAX	19B1323FAX	19J0401FAX	20B1416FAX	20F1004FAX	20J1201FAX		
18J0047FAX	19B1510FAX	19J0402FAX	20B1417FAX	20F1005FAX	20J1293FAX		
18J1168FAX	19B1511FAX	19J0403FAX	20C0037FAX	20F1006FAX	20J1294FAX		
18J1169FAX	19B1514FAX	19J0540FAX	20C0038FAX	20F1007FAX	20J1295FAX		
18J1170FAX	19B1515FAX	19J0541FAX	20C0039FAX	20F1043FAX	20J1408FAX		
18J1234FAX	19C0129FAX	19J0542FAX	20C0173FAX	20F1044FAX	20K0017FAX		
18J1235FAX	19C0130FAX	19J0652FAX	20C0174FAX	20F1045FAX	20K0018FAX		
18J1236FAX	19C0131FAX	19J0653FAX	20C0175FAX	20F1046FAX	20K0019FAX		
18J1237FAX	19C0132FAX	19J0654FAX	20C0189FAX	20F1176FAX	20K0024FAX		
18J1238FAX	19C0133FAX	19J0967FAX	20C0190FAX	20F1177FAX	20K0025FAX		
18K0126FAX	19C0278FAX	19J0968FAX	20C0191FAX	20F1178FAX	20K0026FAX		
18K0127FAX	19C0381FAX	19J0969FAX	20C0192FAX	20F1179FAX	20K0027FAX		
18K0128FAX	19C0382FAX	19J1062FAX	20C0193FAX	20F1295FAX	20K0028FAX		
18K0129FAX	19C0383FAX	19J1063FAX	20C0194FAX	20F1296FAX	20K0213FAX		
18K0130FAX	19C0384FAX	19J1064FAX	20C0195FAX	20F1297FAX	20K0214FAX		
18K0237FAX	19C0416FAX	19J1066FAX	20C0196FAX	20F1298FAX	20K0215FAX		
18K0238FAX	19C0417FAX	19J1161FAX	20C0371FAX	20G0058FAX	20K0216FAX		
18K0239FAX	19C0418FAX	19J1162FAX	20C0372FAX	20G0059FAX	20K0217FAX		
18K0240FAX	19C0419FAX	19J1163FAX	20C0434FAX	20G0060FAX	20K0240FAX		
18K0241FAX	19C0420FAX	19J1164FAX	20C0435FAX	20G0061FAX	20K0241FAX		
18K0354FAX	19C0421FAX	19J1330FAX	20C0436FAX	20G0127FAX	20K0242FAX		
18K0355FAX	19C0984FAX	19J1331FAX	20C0437FAX	20G0128FAX	20K0243FAX		
18K0356FAX	19C0985FAX	19J1332FAX	20C0438FAX	20G0129FAX	20K0244FAX		
18K0357FAX	19C0990FAX	19J1333FAX	20C0439FAX	20G0130FAX	20K0410FAX		
18K0358FAX	19C1020FAX	19K0032FAX	20C0500FAX	20G0228FAX	20K0411FAX		
18K0448FAX	19C1021FAX	19K0098FAX	20C0538FAX	20G0244FAX	20K0412FAX		
18K0449FAX	19C1119FAX	19K0099FAX	20C0539FAX	20G0245FAX	20K0413FAX		
18K0450FAX	19C1120FAX	19K0100FAX	20C0540FAX	20G0246FAX	20K0414FAX		
18K0451FAX	19C1121FAX	19K0101FAX	20C0541FAX	20G0247FAX	20K0415FAX		
18K0452FAX	19C1122FAX	19K0102FAX	20C0542FAX	20G0248FAX	20K0416FAX		
18K0547FAX	19C1123FAX	19K0103FAX	20C0556FAX	20G0249FAX	20K0481FAX		
18K0548FAX	19C1209FAX	19K0292FAX	20C0557FAX	20G0251FAX	20K0482FAX		
18K0549FAX	19C1210FAX	19K0293FAX	20C0558FAX	20G0252FAX	20K0486FAX		
18K0550FAX	19C1211FAX	19K0294FAX	20C0559FAX	20G0253FAX	20K0487FAX		
18K0551FAX	19C1233FAX	19K0295FAX	20C0620FAX	20G0356FAX	20K0654FAX		
18K0555FAX	19C1234FAX	19K0296FAX	20C0621FAX	20G0357FAX	20K0655FAX		
18K0556FAX	19C1353FAX	19K0297FAX	20C0622FAX	20G0358FAX	20K0656FAX		
18K0557FAX	19C1354FAX	19K0383FAX	20C0810FAX	20G0359FAX	20K0657FAX		
18K0603FAX	19C1355FAX	19K0384FAX	20C0811FAX	20G0360FAX	20K0658FAX		
18K0604FAX	19C1356FAX	19K0385FAX	20C0812FAX	20G0361FAX	20K0885FAX		
18K0744FAX	19C1357FAX	19K0386FAX	20C0813FAX	20G0362FAX	20K0886FAX		
18K0769FAX	19C1358FAX	19K0387FAX	20C0814FAX	20G0363FAX	20K0887FAX		
18K0772FAX	19C1359FAX	19K0558FAX	20C0815FAX	20G0365FAX	20K0888FAX		
18K0773FAX	19C1360FAX	19K0559FAX	20C0816FAX	20G0366FAX	20K0889FAX		
18K0774FAX	19C1527FAX	19K0560FAX	20C0817FAX	20G0367FAX	20K0890FAX		
18K0832FAX	19C1528FAX	19K0561FAX	20C0818FAX	20G0708FAX	20K0891FAX		
18K0833FAX	19C1529FAX	19K0646FAX	20C0819FAX	20G0709FAX	20K0892FAX		
18K0834FAX	19C1530FAX	19K0647FAX	20C1132FAX	20G0710FAX	20K0893FAX		
18K0835FAX	19D0016FAX	19K0648FAX	20C1133FAX	20G0711FAX	20K0941FAX		
18K0836FAX	19D0017FAX	19K0649FAX	20C1134FAX	20G0712FAX	20K0942FAX		
18K0912FAX	19D0018FAX	19K0663FAX	20C1135FAX	20G0713FAX	20K0943FAX		
18K1009FAX	19D0019FAX	19K0664FAX	20C1136FAX	20G0748FAX	20K1026FAX		
18K1010FAX	19D0020FAX	19K0665FAX	20C1137FAX	20G0749FAX	20K1027FAX		
18K1011FAX	19D0299FAX	19K0975FAX	20C1138FAX	20G0750FAX	20K1028FAX		
18K1012FAX	19D0300FAX	19K1030FAX	20C1139FAX	20G0751FAX	20K1029FAX		
18K1013FAX	19D0303FAX	19K1031FAX	20C1140FAX	20G0753FAX	20K1164FAX		
18K1025FAX	19D0363FAX	19K1039FAX	20C1141FAX	20G0755FAX	20K1165FAX		
18K1157FAX	19D0364FAX	19K1058FAX	20C1193FAX	20G0840FAX	20K1166FAX		
18K1158FAX	19D0365FAX	19K1059FAX	20C1194FAX	20G0841FAX	20K1167FAX		

Item Code	Description	Affected Lot Numbers					
		18K1186FAX	19D0471FAX	19K1060FAX	20C1195FAX	20G0842FAX	20K1168FAX
		18K1187FAX	19D0472FAX	19K1103FAX	20C1211FAX	20G0843FAX	20K1169FAX
		18K1188FAX	19D0473FAX	19K1104FAX	20C1212FAX	20G0844FAX	20K1415FAX
		18K1189FAX	19D0494FAX	19K1105FAX	20C1213FAX	20G0845FAX	20K1416FAX
		18K1190FAX	19D0495FAX	19K1111FAX	20C1214FAX	20G0846FAX	20K1417FAX
		18K1300FAX	19D0623FAX	19K1195FAX	20C1215FAX	20G0847FAX	20K1418FAX
		18K1301FAX	19D0624FAX	19K1196FAX	20C1216FAX	20G0874FAX	20K1419FAX
		18K1302FAX	19D0625FAX	19K1225FAX	20D0036FAX	20G0875FAX	20K1420FAX
		18K1303FAX	19D0705FAX	19K1243FAX	20D0037FAX	20G0876FAX	20K1561FAX
		18L0010FAX	19D0706FAX	19K1244FAX	20D0038FAX	20G0877FAX	20K1562FAX
		18L0011FAX	19D0707FAX	19K1245FAX	20D0039FAX	20G0878FAX	20K1563FAX
		18L0012FAX	19D0736FAX	19K1346FAX	20D0040FAX	20G0879FAX	20K1564FAX
		18L0013FAX	19D0737FAX	19K1347FAX	20D0041FAX	20G0880FAX	20K1565FAX
		18L0014FAX	19D0738FAX	19K1348FAX	20D0138FAX	20G0881FAX	20K1566FAX
		18L0015FAX	19D0739FAX	19K1349FAX	20D0139FAX	20G1069FAX	20K1567FAX
		18L0219FAX	19D0942FAX	19K1350FAX	20D0140FAX	20G1070FAX	20K1569FAX
		18L0220FAX	19D0943FAX	19K1351FAX	20D0141FAX	20G1071FAX	20L0127FAX
		18L0221FAX	19D0944FAX	19K1352FAX	20D0142FAX	20G1075FAX	20L0131FAX
		18L0222FAX	19D1257FAX	19K1353FAX	20D0143FAX	20G1076FAX	20L0132FAX
		18L0223FAX	19D1258FAX	19K1354FAX	20D0144FAX	20G1187FAX	20L0133FAX
		18L0224FAX	19D1259FAX	19K1355FAX	20D0145FAX	20G1188FAX	20L0134FAX
		18L0260FAX	19D1260FAX	19K1356FAX	20D0275FAX	20G1189FAX	20L0135FAX
		18L0261FAX	19D1307FAX	19L0149FAX	20D0276FAX	20G1190FAX	20L0276FAX
		18L0262FAX	19D1308FAX	19L0150FAX	20D0277FAX	20G1191FAX	20L0277FAX
		18L0263FAX	19D1309FAX	19L0151FAX	20D0278FAX	20G1192FAX	20L0278FAX
		18L0350FAX	19D1310FAX	19L0152FAX	20D0279FAX	20G1193FAX	20L0279FAX
		18L0351FAX	19E0468FAX	19L0153FAX	20D0280FAX	20G1194FAX	20L0280FAX
		18L0352FAX	19E0469FAX	19L0154FAX	20D0432FAX	20G1195FAX	20L0281FAX
		18L0353FAX	19E0470FAX	19L0155FAX	20D0433FAX	20G1264FAX	20L0333FAX
		18L0354FAX	19E0471FAX	19L0234FAX	20D0434FAX	20G1292FAX	20L0334FAX
		18L0355FAX	19E0544FAX	19L0236FAX	20D0435FAX	20G1293FAX	20L0335FAX
		18L0555FAX	19E0545FAX	19L0238FAX	20D0436FAX	20G1294FAX	20L0336FAX
		18L0556FAX	19E0546FAX	19L0239FAX	20D0437FAX	20G1295FAX	20L0337FAX
		18L0557FAX	19E0547FAX	19L0269FAX	20D0438FAX	20G1296FAX	20L0370FAX
		18L0558FAX	19E0567FAX	19L0270FAX	20D0439FAX	20G1297FAX	20L0457FAX
		18L0559FAX	19E0568FAX	19L0271FAX	20D0440FAX	20G1298FAX	20L0458FAX
		18L0704FAX	19E0741FAX	19L0354FAX	20D0441FAX	20G1299FAX	20L0459FAX
		18L0705FAX	19E0742FAX	19L0355FAX	20D0651FAX	20G1360FAX	20L0460FAX
		18L0706FAX	19E0743FAX	19L0356FAX	20D0652FAX	20G1361FAX	20L0513FAX
		18L0804FAX	19E0811FAX	19L0496FAX	20D0765FAX	20G1362FAX	20L0514FAX
		18L0805FAX	19E0812FAX	19L0497FAX	20D0766FAX	20G1363FAX	20L0515FAX
		18L0875FAX	19E0813FAX	19L0530FAX	20D0767FAX	20G1364FAX	20L0516FAX
		18L0876FAX	19E0906FAX	19L0531FAX	20D0768FAX	20G1365FAX	20L0517FAX
		18L0877FAX	19E0907FAX	19L0634FAX	20D0822FAX	20G1366FAX	20L0518FAX
		18L0878FAX	19E0908FAX	19L0635FAX	20D0823FAX	20G1367FAX	20L0652FAX
		18L0879FAX	19E1027FAX	19L0636FAX	20D0824FAX	20G1368FAX	20L0653FAX
		18L0880FAX	19E1028FAX	19L0653FAX	20D0825FAX	20G1370FAX	20L0654FAX
		18L1114FAX	19E1029FAX	19L0801FAX	20D0939FAX	20G1371FAX	20L0655FAX
		18L1115FAX	19E1030FAX	19L0805FAX	20D0940FAX	20G1488FAX	20L0656FAX
		18L1116FAX	19E1031FAX	19L0806FAX	20D0941FAX	20G1489FAX	20L0657FAX
		18L1117FAX	19E1120FAX	19L0864FAX	20D0942FAX	20G1490FAX	20L0672FAX
		18L1268FAX	19E1121FAX	19L0865FAX	20D0943FAX	20G1491FAX	20L0673FAX
		18L1269FAX	19E1293FAX	19L0866FAX	20D0944FAX	20G1492FAX	20L0674FAX
		18L1270FAX	19E1294FAX	19L0885FAX	20D1068FAX	20G1493FAX	20L0675FAX
		18L1271FAX	19E1302FAX	19L0886FAX	20D1069FAX	20G1494FAX	20L0676FAX
		18L1272FAX	19E1303FAX	19L0955FAX	20D1070FAX	20G1495FAX	20L0677FAX
		18L1273FAX	19E1304FAX	19L0956FAX	20D1071FAX	20H0044FAX	20L0678FAX
		18L1274FAX	19E1305FAX	19L0957FAX	20D1134FAX	20H0045FAX	20L0679FAX
		18L1437FAX	19E1306FAX	19L0958FAX	20D1135FAX	20H0046FAX	20L0791FAX
		18L1438FAX	19E1435FAX	19L1125FAX	20D1136FAX	20H0047FAX	20L0792FAX
		18L1439FAX	19E1501FAX	19L1126FAX	20D1137FAX	20H0298FAX	20L0793FAX
		18L1440FAX	19E1502FAX	19L1331FAX	20D1138FAX	20H0299FAX	20L1353FAX
		18M0067FAX	19E1503FAX	19L1332FAX	20D1139FAX	20H0300FAX	20L1354FAX
		18M0069FAX	19E1504FAX	19L1333FAX	20D1140FAX	20H0301FAX	20L1355FAX
		18M0072FAX	19E1556FAX	19L1334FAX	20D1424FAX	20H0302FAX	20L1356FAX
		18M0177FAX	19E1561FAX	19M0183FAX	20D1425FAX	20H0303FAX	20L1357FAX
		18M0178FAX	19E1562FAX	19M0187FAX	20D1426FAX	20H0304FAX	20L1358FAX
		18M0179FAX	19E1563FAX	19M0211FAX	20D1427FAX	20H0313FAX	20L1359FAX
		18M0324FAX	19E1564FAX	19M0212FAX	20D1428FAX	20H0314FAX	20L1360FAX
		18M0325FAX	19E1565FAX	19M0213FAX	20D1503FAX	20H0315FAX	20L1361FAX
		18M0326FAX	19E1688FAX	19M0251FAX	20D1504FAX	20H0316FAX	20L1362FAX
		18M0368FAX	19E1689FAX	19M0252FAX	20D1505FAX	20H0317FAX	20L1363FAX
		18M0379FAX	19E1690FAX	19M0253FAX	20D1506FAX	20H0318FAX	20L1364FAX
		18M0422FAX	19F0185FAX	19M0254FAX	20D1508FAX	20H0319FAX	20L1365FAX
		18M0423FAX	19F0186FAX	19M0255FAX	20D1512FAX	20H0320FAX	20L1367FAX
		18M0424FAX	19F0187FAX	19M0300FAX	20D1513FAX	20H0321FAX	20M0062FAX
		18M0425FAX	19F0188FAX	19M0301FAX	20D1514FAX	20H0322FAX	20M0063FAX

352/5877 Adult-Pediatric Electrostatic Filter HME

Item Code	Description	Affected Lot Numbers					
18M0481FAX	19F0346FAX	19M0302FAX	20D1515FAX	20H0323FAX	20M0064FAX		
18M0482FAX	19F0347FAX	19M0303FAX	20D1517FAX	20H0507FAX	20M0157FAX		
18M0484FAX	19F0348FAX	19M0304FAX	20D1601FAX	20H0508FAX	20M0158FAX		
18M0526FAX	19F0349FAX	19M0358FAX	20D1602FAX	20H0509FAX	20M0159FAX		
18M0527FAX	19F0350FAX	19M0359FAX	20D1603FAX	20H0510FAX	20M0160FAX		
18M0528FAX	19F0520FAX	19M0360FAX	20D1604FAX	20H0543FAX	20M0161FAX		
18M0529FAX	19F0524FAX	19M0361FAX	20D1605FAX	20H0544FAX	20M0162FAX		
18M0564FAX	19F0525FAX	19M0547FAX	20E0078FAX	20H0545FAX	20M0163FAX		
18M0565FAX	19F0526FAX	19M0549FAX	20E0079FAX	20H0546FAX	20M0471FAX		
18M0566FAX	19F0527FAX	19M0550FAX	20E0080FAX	20H0608FAX	20M0472FAX		
18M0738FAX	19F0528FAX	19M0551FAX	20E0081FAX	20H0609FAX	20M0473FAX		
18M0758FAX	19F0742FAX	19M0553FAX	20E0082FAX	20H0610FAX	20M0474FAX		
18M0759FAX	19F0743FAX	19M0717FAX	20E0469FAX	20H0611FAX	20M0475FAX		
18M0760FAX	19F0744FAX	19M0718FAX	20E0470FAX	20H0685FAX	20M0476FAX		
18M0761FAX	19F0745FAX	19M0719FAX	20E0471FAX	20H0688FAX	20M0572FAX		
18M0941FAX	19F0746FAX	19M0720FAX	20E0472FAX	20H0689FAX	20M0573FAX		
18M0942FAX	19F0747FAX	19M0721FAX	20E0473FAX	20H0690FAX	20M0574FAX		
18M0943FAX	19F0878FAX	19M0915FAX	20E0474FAX	20H0772FAX	20M0651FAX		
18M0944FAX	19F0879FAX	19M0916FAX	20E0647FAX	20H0773FAX	20M0652FAX		
18M0945FAX	19F0882FAX	19M0917FAX	20E0648FAX	20H0774FAX	20M0653FAX		
18M1027FAX	19F0942FAX	19M0918FAX	20E0649FAX	20H0775FAX	20M0654FAX		
18M1028FAX	19F0943FAX	20A0094FAX	20E0650FAX	20H0776FAX	20M0655FAX		
18M1029FAX	19F0944FAX	20A0095FAX	20E0651FAX	20H0777FAX	20M0685FAX		
18M1030FAX	19F0945FAX	20A0096FAX	20E0652FAX	20H0778FAX	20M0686FAX		
18M1031FAX	19F0946FAX	20A0098FAX	20E0653FAX	20H0779FAX	20M0687FAX		
18M1032FAX	19F0973FAX	20A0156FAX	20E0654FAX	20H0857FAX	20M0688FAX		
19A0180FAX	19F0974FAX	20A0219FAX	20E0655FAX	20H0858FAX	20M0689FAX		
19A0181FAX	19F1164FAX	20A0220FAX	20E0730FAX	20H0859FAX	20M0690FAX		
19A0183FAX	19F1165FAX	20A0221FAX	20E0731FAX	20H0860FAX	20M0691FAX		
19A0184FAX	19F1166FAX	20A0222FAX	20E0732FAX	20H1067FAX	20M0692FAX		
19A0406FAX	19F1171FAX	20A0370FAX	20E0733FAX	20H1068FAX	20M0741FAX		
19A0407FAX	19F1172FAX	20A0371FAX	20E0734FAX	20H1133FAX	20M0742FAX		
19A0451FAX	19F1173FAX	20A0372FAX	20E0735FAX	20H1134FAX	20M0743FAX		
19A0452FAX	19F1174FAX	20A0373FAX	20E0941FAX	20H1135FAX	20M0744FAX		
19A0453FAX	19F1175FAX	20A0404FAX	20E0942FAX	20H1136FAX	20M0745FAX		
19A0487FAX	19F1215FAX	20A0406FAX	20E0943FAX	20H1137FAX	20M0746FAX		
19A0489FAX	19F1216FAX	20A0514FAX	20E0944FAX	20H1138FAX	20M0835FAX		
19A0490FAX	19F1217FAX	20A0515FAX	20E0945FAX	20H1139FAX	20M0836FAX		
19A0649FAX	19F1218FAX	20A0516FAX	20E0946FAX	20H1140FAX	20M0837FAX		
19A0650FAX	19F1219FAX	20A0517FAX	20E0947FAX	20H1223FAX	20M0838FAX		
19A0651FAX	19G0107FAX	20A0625FAX	20E0948FAX	20H1224FAX	20M1378FAX		
19A0652FAX	19G0108FAX	20A0626FAX	20E1063FAX	20H1225FAX	20M1379FAX		
19A0669FAX	19G0109FAX	20A0627FAX	20E1064FAX	20H1226FAX	20M1380FAX		
19A0670FAX	19G0110FAX	20A0628FAX	20E1065FAX	20H1293FAX	20M1381FAX		
19A0688FAX	19G0111FAX	20A0629FAX	20E1066FAX	20H1294FAX	20M1382FAX		
19A0689FAX	19G0112FAX	20A0656FAX	20E1067FAX	20H1295FAX	20M1466FAX		
19A0690FAX	19G0217FAX	20A0657FAX	20E1068FAX	20H1296FAX	20M1467FAX		
19A0691FAX	19G0218FAX	20A0658FAX	20E1069FAX	20J0100FAX	20M1468FAX		
19A0903FAX	19G0219FAX	20A0726FAX	20E1070FAX	20J0101FAX	20M1503FAX		
19A0904FAX	19G0220FAX	20A0727FAX	20E1144FAX	20J0102FAX	20M1504FAX		
19A0905FAX	19G0221FAX	20A0728FAX	20E1145FAX	20J0103FAX	20M1505FAX		
19A0906FAX	19G0222FAX	20A0729FAX	20E1146FAX	20J0104FAX	20M1506FAX		
19A0907FAX	19G0402FAX	20A0730FAX	20E1147FAX	20J0105FAX	20M1507FAX		
19A1131FAX	19G0403FAX	20A0846FAX	20E1148FAX	20J0106FAX	20M1508FAX		
19A1132FAX	19G0404FAX	20A1068FAX	20E1149FAX	20J0107FAX	21A0168FAX		
19A1133FAX	19G0405FAX	20A1069FAX	20E1150FAX	20J0108FAX	21A0169FAX		
19A1134FAX	19G0406FAX	20A1070FAX	20E1151FAX	20J0109FAX	21A0170FAX		
19A1135FAX	19G0407FAX	20A1071FAX	20E1419FAX	20J0291FAX	21A0171FAX		
19A1136FAX	19G0476FAX	20A1072FAX	20E1420FAX	20J0292FAX	21A0172FAX		
19A1280FAX	19G0477FAX	20A1074FAX	20E1421FAX	20J0293FAX	21A0173FAX		
19A1281FAX	19G0657FAX	20A1245FAX	20E1422FAX	20J0294FAX	21A0174FAX		
19A1282FAX	19G0659FAX	20A1246FAX	20F0022FAX	20J0295FAX	21A0175FAX		
19A1352FAX	19G0660FAX	20A1247FAX	20F0023FAX	20J0296FAX	21A0245FAX		
19A1353FAX	19G0661FAX	20A1248FAX	20F0024FAX	20J0297FAX	21A0246FAX		
19A1354FAX	19G0768FAX	20A1303FAX	20F0025FAX	20J0298FAX	21A0247FAX		
19A1355FAX	19G0811FAX	20A1304FAX	20F0026FAX	20J0299FAX	21A0248FAX		
19A1393FAX	19G0812FAX	20A1305FAX	20F0174FAX	20J0300FAX	21A0249FAX		
19A1394FAX	19G0916FAX	20A1306FAX	20F0198FAX	20J0301FAX	21A0250FAX		
19A1395FAX	19G0917FAX	20A1375FAX	20F0199FAX	20J0446FAX	21A0380FAX		
19A1396FAX	19G0918FAX	20A1376FAX	20F0200FAX	20J0447FAX	21A0381FAX		
19B0093FAX	19G1098FAX	20A1378FAX	20F0201FAX	20J0448FAX	21A0382FAX		
19B0149FAX	19G1099FAX	20A1383FAX	20F0202FAX	20J0449FAX	21A0383FAX		
19B0386FAX	19G1239FAX	20A1384FAX	20F0203FAX	20J0450FAX	21A0384FAX		
19B0387FAX	19G1242FAX	20A1406FAX	20F0204FAX	20J0597FAX	21A0385FAX		
19B0388FAX	19G1299FAX	20A1490FAX	20F0205FAX	20J0598FAX	21A0386FAX		
19B0389FAX	19G1300FAX	20A1504FAX	20F0206FAX	20J0599FAX	21A0496FAX		
19B0390FAX	19G1301FAX	20B0174FAX	20F0207FAX	20J0600FAX	21A0497FAX		

Item Code	Description	Affected Lot Numbers					
		19B0391FAX	19G1302FAX	20B0175FAX	20F0394FAX	20J0601FAX	21A0498FAX
		19B0451FAX	19G1303FAX	20B0176FAX	20F0395FAX	20J0602FAX	21A0499FAX
		19B0453FAX	19G1304FAX	20B0177FAX	20F0396FAX	20J0603FAX	21A0500FAX
		19B0455FAX	19G1347FAX	20B0232FAX	20F0397FAX	20J0604FAX	21A0501FAX
		19B0456FAX	19H0681FAX	20B0233FAX	20F0398FAX	20J0606FAX	21A0502FAX
		19B0457FAX	19H0682FAX	20B0234FAX	20F0399FAX	20J0607FAX	21A0578FAX
		19B0523FAX	19H0683FAX	20B0235FAX	20F0400FAX	20J0608FAX	21A0579FAX
		19B0524FAX	19H0684FAX	20B0249FAX	20F0524FAX	20J0609FAX	21A0580FAX
		19B0551FAX	19H0741FAX	20B0346FAX	20F0525FAX	20J0715FAX	21A0581FAX
		19B0552FAX	19H0742FAX	20B0347FAX	20F0564FAX	20J0716FAX	21A0582FAX
		19B0553FAX	19H0743FAX	20B0348FAX	20F0565FAX	20J0717FAX	21A0583FAX
		19B0668FAX	19H0744FAX	20B0438FAX	20F0566FAX	20J0718FAX	21A0584FAX
		19B0669FAX	19H0745FAX	20B0439FAX	20F0567FAX	20J0719FAX	21A0703FAX
		19B0670FAX	19H0746FAX	20B0440FAX	20F0568FAX	20J0720FAX	21A0704FAX
		19B0671FAX	19H0804FAX	20B0441FAX	20F0569FAX	20J0721FAX	21A0784FAX
		19B0672FAX	19H0805FAX	20B0442FAX	20F0570FAX	20J0722FAX	21A0785FAX
		19B0709FAX	19H0806FAX	20B0443FAX	20F0614FAX	20J0890FAX	21A0786FAX
		19B0710FAX	19H0807FAX	20B0444FAX	20F0615FAX	20J0891FAX	21A0787FAX
		19B0711FAX	19H0808FAX	20B0551FAX	20F0616FAX	20J0893FAX	21A0883FAX
		19B0712FAX	19H0960FAX	20B0552FAX	20F0627FAX	20J0894FAX	21A0789FAX
		19B0751FAX	19H0961FAX	20B0553FAX	20F0628FAX	20J0895FAX	21A0790FAX
		19B0911FAX	19H0962FAX	20B0568FAX	20F0629FAX	20J0896FAX	21A0883FAX
		19B0912FAX	19H0963FAX	20B0569FAX	20F0630FAX	20J1026FAX	21A0885FAX
		19B0913FAX	19J0074FAX	20B0636FAX	20F0631FAX	20J1027FAX	21A1126FAX
		19B0914FAX	19J0075FAX	20B0637FAX	20F0632FAX	20J1028FAX	21A1127FAX
		19B0915FAX	19J0076FAX	20B0694FAX	20F0789FAX	20J1029FAX	21A1128FAX
		19B1058FAX	19J0077FAX	20B0695FAX	20F0790FAX	20J1030FAX	21A1129FAX
		19B1059FAX	19J0078FAX	20B0696FAX	20F0791FAX	20J1031FAX	21A1199FAX
		19B1060FAX	19J0228FAX	20B0945FAX	20F0792FAX	20J1032FAX	21A1203FAX
		19B1123FAX	19J0229FAX	20B0998FAX	20F0793FAX	20J1093FAX	21A1204FAX
		19B1124FAX	19J0230FAX	20B1000FAX	20F0794FAX	20J1094FAX	21B0006FAX
		19B1125FAX	19J0231FAX	20B1001FAX	20F0899FAX	20J1127FAX	21B0007FAX
		19B1126FAX	19J0232FAX	20B1101FAX	20F0900FAX	20J1129FAX	
		19B1127FAX	19J0233FAX	20B1102FAX	20F0901FAX	20J1130FAX	
		18L1122FAX	20B0349FAX	20C0626FAX	20D1354FAX	20H0075FAX	20M0287FAX
		18L1123FAX	20B0350FAX	20C1145FAX	20D1355FAX	20H0076FAX	20M0560FAX
		19A1266FAX	20B0362FAX	20C1146FAX	20D1608FAX	20H0717FAX	20M0561FAX
		19A1267FAX	20B0363FAX	20D0360FAX	20D1609FAX	20H1122FAX	20M0944FAX
		19B0936FAX	20B0628FAX	20D0413FAX	20F0495FAX	20H1123FAX	20M0945FAX
		19B1355FAX	20B0629FAX	20D0697FAX	20F0773FAX	20J0218FAX	20M0946FAX
		19C1484FAX	20B1082FAX	20D0891FAX	20G0532FAX	20J0219FAX	20M1236FAX
		19C1485FAX	20B1083FAX	20D0892FAX	20G0533FAX	20J0460FAX	20M1237FAX
		19G0048FAX	20B1175FAX	20D0893FAX	20G0818FAX	20J0750FAX	20M1238FAX
		19H0406FAX	20B1398FAX	20D0960FAX	20G0820FAX	20J0751FAX	20M1454FAX
		19K1068FAX	20B1399FAX	20D0961FAX	20G1135FAX	20L0788FAX	20M1455FAX
		19L0407FAX	20C0144FAX	20D1066FAX	20G1579FAX	20L0789FAX	21A0138FAX
		19L0408FAX	20C0499FAX	20D1067FAX	20H0073FAX	20L0830FAX	
		19L1288FAX	20C0625FAX	20D1332FAX	20H0074FAX	20M0286FAX	
352/5877TC	Adult-Pediatric Electrostatic Filter HME	19A1264FAX	19B1373FAX	19E1489FAX	19E1490FAX	19L0405FAX	19L0406FAX
		20A1187FAX	20D0385FAX	20F0723FAX	20J0462FAX	20J1163FAX	
		18J0351FAX	19F0751FAX	19M0624FAX	20E0160FAX	20G0717FAX	20J0743FAX
		18J0353FAX	19F0752FAX	19M0625FAX	20E0161FAX	20G0718FAX	20J0744FAX
		18J0487FAX	19F1226FAX	19M0933FAX	20E0162FAX	20G0719FAX	20J0745FAX
		19A1265FAX	19F1227FAX	19M1002FAX	20E0163FAX	20G0782FAX	20J0746FAX
		19B1363FAX	19F1228FAX	20B0116FAX	20E0164FAX	20G0783FAX	20J0747FAX
		19B1364FAX	19F1229FAX	20B0117FAX	20E0165FAX	20G0784FAX	20J0748FAX
		19B1365FAX	19G0117FAX	20B0118FAX	20E0390FAX	20G0785FAX	20J0749FAX
		19B1366FAX	19G0118FAX	20B0119FAX	20E0391FAX	20G0786FAX	20J1019FAX
		19B1367FAX	19G0119FAX	20B0120FAX	20E0392FAX	20G0787FAX	20J1020FAX
		19C0118FAX	19G0120FAX	20B0121FAX	20E0393FAX	20G1557FAX	20J1021FAX
		19C0119FAX	19G0121FAX	20B0122FAX	20E0585FAX	20G1558FAX	20J1022FAX
		19C0120FAX	19G0422FAX	20B0370FAX	20E0586FAX	20H0175FAX	20J1076FAX
		19C0121FAX	19G0423FAX	20B0371FAX	20E0587FAX	20H0178FAX	20J1077FAX
		19C0327FAX	19G1478FAX	20C0145FAX	20E0645FAX	20H0181FAX	20J1078FAX
		19D0424FAX	19G1479FAX	20C0146FAX	20E0646FAX	20H0182FAX	20J1326FAX
		19D0425FAX	19J0206FAX	20C0147FAX	20E0657FAX	20H0183FAX	20K0620FAX
		19D0447FAX	19J0207FAX	20C0148FAX	20E0871FAX	20H0184FAX	20K0621FAX
		19D0448FAX	19J0208FAX	20C0276FAX	20E0873FAX	20H0185FAX	20K0622FAX
		19D0449FAX	19J0209FAX	20C0369FAX	20E0874FAX	20H0331FAX	20K0623FAX
		19D0696FAX	19J0210FAX	20C0370FAX	20E0875FAX	20H0332FAX	20K0624FAX
		19D0697FAX	19J0211FAX	20C0462FAX	20E0935FAX	20H0333FAX	20K0625FAX
		19D0955FAX	19J0236FAX	20C0463FAX	20E0936FAX	20H0334FAX	20K0910FAX
		19D0958FAX	19J0237FAX	20C0464FAX	20E0997FAX	20H0335FAX	20K0911FAX
		19D0959FAX	19J0238FAX	20C0750FAX	20E0998FAX	20H0336FAX	20K0912FAX
		19D1251FAX	19J0239FAX	20C0751FAX	20E1072FAX	20H0337FAX	20K0924FAX
		19D1252FAX	19J0256FAX	20C0752FAX	20E1073FAX	20H0338FAX	20K0925FAX
		19D1253FAX	19J0257FAX	20C1016FAX	20F0876FAX	20H0339FAX	20K1036FAX

Item Code	Description	Affected Lot Numbers					
		19D1254FAX	19J0258FAX	20C1164FAX	20F0877FAX	20H0500FAX	20K1037FAX
		19E0180FAX	19K0186FAX	20C1165FAX	20F0881FAX	20H0501FAX	20K1038FAX
		19E0181FAX	19K0197FAX	20C1166FAX	20F0882FAX	20H0502FAX	20K1039FAX
		19E0182FAX	19K0198FAX	20C1235FAX	20F0883FAX	20H0503FAX	20K1040FAX
		19E0183FAX	19K0199FAX	20C1236FAX	20F0884FAX	20H0582FAX	20K1041FAX
		19E0184FAX	19K0368FAX	20D1187FAX	20F1132FAX	20H0583FAX	20K1042FAX
		19E0185FAX	19K0369FAX	20D1188FAX	20F1133FAX	20H0584FAX	20K1236FAX
		19E0354FAX	19K0370FAX	20D1189FAX	20F1134FAX	20H0585FAX	20K1241FAX
		19E0355FAX	19K1409FAX	20D1190FAX	20F1135FAX	20J0302FAX	20K1242FAX
		19E0356FAX	19K1410FAX	20D1191FAX	20F1155FAX	20J0303FAX	20K1243FAX
		19E0463FAX	19K1411FAX	20D1429FAX	20F1156FAX	20J0304FAX	20K1244FAX
		19E0464FAX	19K1412FAX	20D1430FAX	20G0006FAX	20J0305FAX	20K1245FAX
		19E0552FAX	19L0087FAX	20D1432FAX	20G0007FAX	20J0306FAX	20K1246FAX
		19E0553FAX	19L0088FAX	20D1433FAX	20G0224FAX	20J0307FAX	20K1247FAX
		19E1570FAX	19L0089FAX	20D1523FAX	20G0225FAX	20J0308FAX	20K1248FAX
		19E1571FAX	19L1289FAX	20D1524FAX	20G0226FAX	20J0309FAX	20K1249FAX
		19E1572FAX	19L1290FAX	20D1525FAX	20G0345FAX	20J0518FAX	20K1553FAX
		19E1573FAX	19L1291FAX	20D1526FAX	20G0346FAX	20J0519FAX	20K1554FAX
		19F0646FAX	19M0193FAX	20D1527FAX	20G0347FAX	20J0520FAX	
		19F0647FAX	19M0522FAX	20D1528FAX	20G0348FAX	20J0521FAX	
		19F0648FAX	19M0523FAX	20E0158FAX	20G0349FAX	20J0522FAX	
		19F0649FAX	19M0623FAX	20E0159FAX	20G0716FAX	20J0742FAX	
352/5881	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	20A1133FAX	21A0689FAX				
352/5893	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	20D1343FAX	21A0690FAX				
352/5978	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	20A0564FAX	20A1144FAX	20D1345FAX	20H0566FAX	20F0167FAX	20K0461FAX
		18K1370FAX	19B1024FAX	19G1287FAX	20A0677FAX	20F0168FAX	20K0462FAX
		18K1371FAX	19B1025FAX	19G1288FAX	20A0678FAX	20F0210FAX	20K0463FAX
		18K1372FAX	19B1327FAX	19G1289FAX	20A0679FAX	20F0211FAX	20K0464FAX
		18K1373FAX	19B1328FAX	19H0546FAX	20A0680FAX	20F0212FAX	20K0465FAX
		18K1374FAX	19B1329FAX	19H0547FAX	20A0772FAX	20F0214FAX	20K0466FAX
		18K1375FAX	19B1330FAX	19H0548FAX	20A0773FAX	20F0460FAX	20K0467FAX
		18K1376FAX	19B1331FAX	19H0549FAX	20A0997FAX	20F0462FAX	20K0668FAX
		18K1377FAX	19B1332FAX	19H0852FAX	20A0998FAX	20F0463FAX	20L0454FAX
		18L0051FAX	19B1488FAX	19H0853FAX	20A1008FAX	20F0542FAX	20L0455FAX
		18L0052FAX	19B1489FAX	19H0854FAX	20A1009FAX	20F0543FAX	20L0456FAX
		18L0053FAX	19B1490FAX	19H0893FAX	20A1010FAX	20F0544FAX	20L0571FAX
		18L0054FAX	19B1491FAX	19H0894FAX	20A1011FAX	20F0676FAX	20L0572FAX
		18L0055FAX	19C0148FAX	19H0955FAX	20A1012FAX	20F0677FAX	20L0573FAX
		18L0056FAX	19C0149FAX	19H0956FAX	20A1051FAX	20F0777FAX	20L0574FAX
		18L0057FAX	19C0150FAX	19J0644FAX	20A1052FAX	20F0779FAX	20L0717FAX
		18L0499FAX	19C0151FAX	19J0645FAX	20A1054FAX	20F0847FAX	20L0718FAX
		18L0500FAX	19C0152FAX	19J1106FAX	20A1055FAX	20F0848FAX	20L0719FAX
		18L0958FAX	19C0322FAX	19J1107FAX	20A1412FAX	20G0529FAX	20L0843FAX
		18L0959FAX	19C0323FAX	19J1108FAX	20A1413FAX	20G0530FAX	20L0844FAX
		18L0960FAX	19C0324FAX	19K0027FAX	20A1477FAX	20G0531FAX	20L0845FAX
		18L0961FAX	19C0325FAX	19K0028FAX	20A1478FAX	20G0895FAX	20L0846FAX
		18L0962FAX	19C0326FAX	19K0029FAX	20A1479FAX	20G0896FAX	20M0402FAX
		18L1275FAX	19C0437FAX	19K0075FAX	20A1480FAX	20G0897FAX	20M0403FAX
		18L1276FAX	19C0438FAX	19K0076FAX	20A1481FAX	20G0898FAX	20M0404FAX
		18L1277FAX	19C0439FAX	19K0077FAX	20B0557FAX	20G0899FAX	20M0563FAX
		18L1278FAX	19C0440FAX	19K0126FAX	20B0558FAX	20G0900FAX	20M0564FAX
		18L1279FAX	19C0441FAX	19K0127FAX	20B0642FAX	20G0901FAX	20M0565FAX
		18L1359FAX	19C0840FAX	19K0129FAX	20B0643FAX	20G0902FAX	20M0566FAX
		18L1360FAX	19D0121FAX	19K0130FAX	20B0644FAX	20G0903FAX	20M0636FAX
		18L1361FAX	19D0122FAX	19K0131FAX	20B0645FAX	20G0904FAX	20M0637FAX
		18L1362FAX	19D0123FAX	19K0976FAX	20B0698FAX	20G0905FAX	20M0638FAX
		18M0122FAX	19D0125FAX	19K0977FAX	20B0699FAX	20G1196FAX	20M0639FAX
		18M0123FAX	19D0351FAX	19K1062FAX	20B0878FAX	20G1197FAX	20M0640FAX
		18M0124FAX	19D0352FAX	19K1063FAX	20B0880FAX	20G1198FAX	20M0641FAX
		18M0125FAX	19E0269FAX	19K1064FAX	20B0988FAX	20G1199FAX	20M0642FAX
		18M0273FAX	19E0271FAX	19K1072FAX	20B0989FAX	20G1200FAX	20M0643FAX
		18M0274FAX	19E0272FAX	19K1073FAX	20B0990FAX	20G1201FAX	20M1339FAX
		18M1010FAX	19E0273FAX	19K1074FAX	20B0991FAX	20G1202FAX	20M1340FAX
		18M1011FAX	19E0702FAX	19K1075FAX	20B1097FAX	20H0787FAX	20M1341FAX
		18M1012FAX	19E0703FAX	19K1285FAX	20B1127FAX	20H0788FAX	20M1342FAX
		18M1013FAX	19E0801FAX	19K1286FAX	20B1128FAX	20H0789FAX	20M1343FAX
		18M1014FAX	19E0802FAX	19L0411FAX	20B1129FAX	20H0790FAX	20M1344FAX
		19A0149FAX	19E0803FAX	19L0412FAX	20B1292FAX	20H0791FAX	20M1345FAX
		19A0150FAX	19E0804FAX	19L0417FAX	20B1298FAX	20H0792FAX	20M1346FAX
		19A0151FAX	19E0966FAX	19L0814FAX	20B1405FAX	20H0793FAX	20M1347FAX
		19A0152FAX	19E0967FAX	19L0815FAX	20B1406FAX	20H1100FAX	21A0162FAX
		19A0153FAX	19E0968FAX	19L0816FAX	20C0073FAX	20H1101FAX	21A0163FAX
		19A0602FAX	19E0969FAX	19L0817FAX	20C0647FAX	20H1103FAX	21A0164FAX
		19A0603FAX	19E1271FAX	19L0818FAX	20C0648FAX	20H1104FAX	21A0165FAX
352/5996	Adult-Pediatric Electrostatic Filter HME						

Item Code	Description	Affected Lot Numbers					
		19A0604FAX	19E1272FAX	19L0819FAX	20C0649FAX	20H1105FAX	21A0166FAX
		19A1159FAX	19E1273FAX	19L0959FAX	20C0650FAX	20H1106FAX	21A0167FAX
		19A1160FAX	19E1274FAX	19L0960FAX	20C0803FAX	20H1107FAX	21A0237FAX
		19A1161FAX	19E1310FAX	19L0961FAX	20D0110FAX	20H1298FAX	21A0238FAX
		19A1162FAX	19E1311FAX	19L0962FAX	20D0111FAX	20H1299FAX	21A0240FAX
		19A1227FAX	19E1312FAX	19L1032FAX	20D0112FAX	20H1300FAX	21A0241FAX
		19A1228FAX	19E1320FAX	19L1033FAX	20D0154FAX	20H1301FAX	21A0242FAX
		19A1229FAX	19E1483FAX	19L1281FAX	20D0155FAX	20H1302FAX	21A0243FAX
		19A1327FAX	19E1484FAX	19L1282FAX	20D0156FAX	20H1303FAX	21A0244FAX
		19B0019FAX	19E1485FAX	19L1382FAX	20D0157FAX	20J0212FAX	21A0394FAX
		19B0020FAX	19E1486FAX	19M0056FAX	20D0158FAX	20J0213FAX	21A0395FAX
		19B0021FAX	19E1487FAX	19M0583FAX	20D0159FAX	20J0214FAX	21A0396FAX
		19B0022FAX	19E1590FAX	19M0585FAX	20D0287FAX	20J0215FAX	21A0397FAX
		19B0023FAX	19E1591FAX	19M0706FAX	20D0288FAX	20J0216FAX	21A0398FAX
		19B0103FAX	19F0030FAX	19M0936FAX	20D0289FAX	20J0217FAX	21A0399FAX
		19B0370FAX	19F0031FAX	20A0085FAX	20D0366FAX	20J0221FAX	21A0505FAX
		19B0371FAX	19F0189FAX	20A0086FAX	20D0527FAX	20J0222FAX	21A0506FAX
		19B0415FAX	19F0190FAX	20A0087FAX	20D0531FAX	20J1235FAX	21A0507FAX
		19B0416FAX	19F0191FAX	20A0088FAX	20E1106FAX	20J1236FAX	21A0508FAX
		19B0417FAX	19F0629FAX	20A0089FAX	20E1107FAX	20J1237FAX	21A0509FAX
		19B0419FAX	19G0135FAX	20A0091FAX	20E1108FAX	20J1413FAX	21A0510FAX
		19B0601FAX	19G0139FAX	20A0093FAX	20E1109FAX	20J1414FAX	21A0511FAX
		19B0602FAX	19G0140FAX	20A0224FAX	20E1260FAX	20J1415FAX	21A0999FAX
		19B0603FAX	19G0143FAX	20A0225FAX	20E1261FAX	20K0113FAX	21A1000FAX
		19B0610FAX	19G0146FAX	20A0293FAX	20E1262FAX	20K0114FAX	
		19B0611FAX	19G0549FAX	20A0294FAX	20E1263FAX	20K0257FAX	
		19B0612FAX	19G0550FAX	20A0468FAX	20E1264FAX	20K0258FAX	
		19B0642FAX	19G0551FAX	20A0471FAX	20E1265FAX	20K0259FAX	
		19B0643FAX	19G0552FAX	20A0570FAX	20E1266FAX	20K0260FAX	
		19B0644FAX	19G1285FAX	20A0630FAX	20E1267FAX	20K0261FAX	
		19B1023FAX	19G1286FAX	20A0631FAX	20E1268FAX	20K0262FAX	
352/5996TC	Adult-Pediatric Electrostatic Filter HME	18L0485FAX	19G1297FAX	20J0090FAX			
		19A0409FAX	19J1117FAX	20L0524FAX			
		19A1270FAX	19L1103FAX	20M0400FAX			
		19A1356FAX	20A1186FAX				
		19E1357FAX	20E1098FAX				
352/5996Z	Adult-Pediatric Electrostatic Filter HME	19A1285FAX	19J1268FAX	20F0372FAX	20J1234FAX		
		19A1286FAX	19J1269FAX	20F0373FAX	20K0268FAX		
		19C0143FAX	19L0413FAX	20H1297FAX	20L0453FAX		
		19E0270FAX	20B0877FAX	20J1232FAX	20M0401FAX		
		19G0223FAX	20E1097FAX	20J1233FAX	20M1336FAX		
352S19041Z	Adult-Pediatric Electrostatic Filter HME with Connector	19A0291FAX	19E1498FAX	20B0901FAX	20B1266FAX	20F1309FAX	
		19D0286FAX	19E1499FAX	20B0902FAX	20B1267FAX	20F1311FAX	
		19D0455FAX	19E1500FAX	20B0903FAX	20F0719FAX	20K0583FAX	
		19D0456FAX	19K0033FAX	20B0906FAX	20F0720FAX	20K0584FAX	
		19D0457FAX	19K1386FAX	20B1264FAX	20F0721FAX	20K0585FAX	
		19D0458FAX	19L1100FAX	20B1265FAX	20F0722FAX	20L0523FAX	
352S19043	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	21A0693FAX					
352S5855BR	Filtro e Umidificador Respiratório DAR com Cateter Mount	19B1147FAX	19L0124FAX	20C0263FAX	20H0795FAX	20J0540FAX	20L0290FAX
		19B1148FAX	19L0125FAX	20C0264FAX	20H0796FAX	20J0541FAX	20L0291FAX
		19B1149FAX	19M0035FAX	20C0504FAX	20H1154FAX	20J0794FAX	20L0322FAX
		19B1207FAX	19M0036FAX	20H0295FAX	20H1156FAX	20J0795FAX	20L0323FAX
		19B1208FAX	19M0050FAX	20H0296FAX	20H1157FAX	20J0796FAX	20L0641FAX
		19B1209FAX	19M0051FAX	20H0297FAX	20H1158FAX	20J0797FAX	20L0642FAX
		19B1210FAX	20A0113FAX	20H0341FAX	20H1159FAX	20J0798FAX	20L0643FAX
		19C0106FAX	20A0114FAX	20H0342FAX	20H1205FAX	20K0556FAX	20L0644FAX
		19C0107FAX	20B1229FAX	20H0343FAX	20J0370FAX	20K0557FAX	20L0645FAX
		19C0108FAX	20B1230FAX	20H0344FAX	20J0371FAX	20K0558FAX	20L0739FAX
		19C0109FAX	20B1418FAX	20H0345FAX	20J0372FAX	20K0559FAX	20L0740FAX
		19C0184FAX	20B1419FAX	20H0346FAX	20J0373FAX	20K0715FAX	20L0741FAX
		19C0185FAX	20B1420FAX	20H0586FAX	20J0374FAX	20K0859FAX	20L0742FAX
		19C0186FAX	20B1421FAX	20H0587FAX	20J0375FAX	20K0860FAX	20L0743FAX
		19C0187FAX	20B1422FAX	20H0588FAX	20J0376FAX	20L0048FAX	20L0799FAX
		19C0189FAX	20B1423FAX	20H0589FAX	20J0377FAX	20L0049FAX	20L0800FAX
		19C0190FAX	20B1424FAX	20H0691FAX	20J0378FAX	20L0050FAX	20L1013FAX
		19K1298FAX	20C0232FAX	20H0692FAX	20J0379FAX	20L0143FAX	20L1236FAX
		19K1299FAX	20C0233FAX	20H0693FAX	20J0537FAX	20L0144FAX	
		19K1300FAX	20C0234FAX	20H0694FAX	20J0538FAX	20L0288FAX	
		19L0123FAX	20C0262FAX	20H0794FAX	20J0539FAX	20L0289FAX	
		18J0048FAX	19B0525FAX	19F0461FAX	19L0272FAX	20C0795FAX	20J0615FAX
		18J0049FAX	19B0526FAX	19F0463FAX	19L0273FAX	20C0796FAX	20J0616FAX
		18J0144FAX	19B0527FAX	19F0464FAX	19L0274FAX	20C0934FAX	20J0617FAX
		18J0145FAX	19B0544FAX	19F0529FAX	19L0275FAX	20C0935FAX	20J0969FAX
		18J0364FAX	19B0545FAX	19F0653FAX	19L0555FAX	20C1123FAX	20J1086FAX
		18K0605FAX	19B0676FAX	19F0654FAX	19L0556FAX	20C1124FAX	20J1347FAX
		18K0606FAX	19B0677FAX	19F0655FAX	19L0557FAX	20C1125FAX	20J1348FAX

Item Code	Description	Affected Lot Numbers					
		18K0745FAX	19B0678FAX	19F0656FAX	19L0927FAX	20C1126FAX	20J1349FAX
		18K0746FAX	19B0679FAX	19F0739FAX	19L0928FAX	20D0101FAX	20K0102FAX
		18K0907FAX	19B0979FAX	19F0740FAX	19L0929FAX	20D0102FAX	20K0103FAX
		18K0908FAX	19B0980FAX	19F0741FAX	19L1268FAX	20D0103FAX	20K0104FAX
		18K0909FAX	19B0981FAX	19F0953FAX	19L1269FAX	20D0104FAX	20K0105FAX
		18K0910FAX	19B0982FAX	19F1200FAX	19L1271FAX	20D0105FAX	20K1052FAX
		18K0911FAX	19B1201FAX	19F1201FAX	19L1272FAX	20D0106FAX	20K1053FAX
		18K1191FAX	19B1202FAX	19F1266FAX	19L1274FAX	20D0653FAX	20K1054FAX
		18K1192FAX	19B1203FAX	19G0045FAX	19L1340FAX	20D0654FAX	20K1413FAX
		18K1193FAX	19B1204FAX	19G0224FAX	19M0195FAX	20D0759FAX	20K1414FAX
		18K1327FAX	19B1507FAX	19G0225FAX	19M0196FAX	20D0760FAX	20L0041FAX
		18K1328FAX	19B1508FAX	19G0226FAX	19M0197FAX	20D0761FAX	20L0042FAX
		18L0016FAX	19B1521FAX	19G0227FAX	19M0198FAX	20D0762FAX	20L0043FAX
		18L0017FAX	19B1522FAX	19G0228FAX	19M0410FAX	20D0763FAX	20L0363FAX
		18L0018FAX	19B1523FAX	19G0280FAX	19M0411FAX	20D0764FAX	20L0364FAX
		18L0264FAX	19C0126FAX	19G0281FAX	19M0517FAX	20D1076FAX	20L0366FAX
		18L0265FAX	19C0127FAX	19G0282FAX	19M0518FAX	20D1077FAX	20L0367FAX
		18L0266FAX	19C0128FAX	19G0479FAX	19M0519FAX	20D1078FAX	20L0368FAX
		18L0267FAX	19C0433FAX	19G0480FAX	19M0923FAX	20D1405FAX	20L0369FAX
		18L0560FAX	19C0434FAX	19G0481FAX	20A0083FAX	20D1406FAX	20L0546FAX
		18L0561FAX	19C0435FAX	19G0482FAX	20A0084FAX	20D1407FAX	20L0547FAX
		18L0562FAX	19C0436FAX	19G0483FAX	20A0105FAX	20D1592FAX	20L0548FAX
		18L0563FAX	19C0828FAX	19G0778FAX	20A0106FAX	20D1593FAX	20L0549FAX
		18L0590FAX	19C0837FAX	19G0840FAX	20A0107FAX	20D1594FAX	20L0550FAX
		18L0591FAX	19C0838FAX	19G1243FAX	20A0108FAX	20D1595FAX	20L0795FAX
		18L0592FAX	19C0839FAX	19G1244FAX	20A0374FAX	20E0083FAX	20L0796FAX
		18L0808FAX	19C0897FAX	19G1245FAX	20A0375FAX	20E0084FAX	20L0797FAX
		18L0809FAX	19C0898FAX	19G1246FAX	20A0376FAX	20E0085FAX	20L0798FAX
		18L0810FAX	19C0920FAX	19H0776FAX	20A0693FAX	20E0088FAX	20L1295FAX
		18L0811FAX	19C1262FAX	19H0777FAX	20A0694FAX	20E0089FAX	20L1296FAX
		18L0812FAX	19C1263FAX	19H0778FAX	20A0695FAX	20E0569FAX	20L1297FAX
		18L0813FAX	19C1265FAX	19H0895FAX	20A0717FAX	20E0570FAX	20L1298FAX
		18L1183FAX	19C1266FAX	19H0898FAX	20A1048FAX	20E0571FAX	20L1299FAX
		18L1184FAX	19C1364FAX	19H0899FAX	20A1049FAX	20E0643FAX	20L1334FAX
		18L1185FAX	19C1365FAX	19J0147FAX	20A1091FAX	20E0644FAX	20L1335FAX
		18L1186FAX	19C1366FAX	19J0148FAX	20A1258FAX	20E0724FAX	20L1336FAX
		18L1355FAX	19C1367FAX	19J0149FAX	20A1259FAX	20E0725FAX	20L1337FAX
353/19004	HME for Tracheostomized Patients	18L1356FAX	19C1368FAX	19J0150FAX	20A1260FAX	20E0726FAX	20L1338FAX
		18L1357FAX	19D0150FAX	19J0301FAX	20A1403FAX	20E0727FAX	20M0156FAX
		18L1358FAX	19D0151FAX	19J0302FAX	20A1404FAX	20E1274FAX	20M0609FAX
		18M0249FAX	19D0395FAX	19J0303FAX	20A1405FAX	20E1275FAX	20M0610FAX
		18M0250FAX	19D0397FAX	19J0304FAX	20A1505FAX	20F0036FAX	20M0611FAX
		18M0251FAX	19D0508FAX	19J0433FAX	20A1506FAX	20F0037FAX	20M0612FAX
		18M0252FAX	19D0509FAX	19J0434FAX	20B0239FAX	20F0303FAX	20M0613FAX
		18M0533FAX	19D0511FAX	19J0435FAX	20B0240FAX	20F0304FAX	20M0734FAX
		18M0534FAX	19D0512FAX	19J0466FAX	20B0241FAX	20F0305FAX	20M0735FAX
		18M0535FAX	19D0860FAX	19J0467FAX	20B0242FAX	20F0306FAX	20M0736FAX
		18M0536FAX	19D1052FAX	19J0469FAX	20B0243FAX	20F0452FAX	20M0737FAX
		18M0537FAX	19D1053FAX	19J0470FAX	20B0592FAX	20F0453FAX	20M0738FAX
		18M1054FAX	19D1146FAX	19J0729FAX	20B0593FAX	20F0454FAX	20M0739FAX
		18M1055FAX	19D1147FAX	19J0730FAX	20B0594FAX	20F0455FAX	20M0740FAX
		18M1056FAX	19D1148FAX	19J0731FAX	20B0633FAX	20F0909FAX	20M1041FAX
		18M1057FAX	19E0510FAX	19J1256FAX	20B0634FAX	20F0910FAX	20M1042FAX
		19A0491FAX	19E0511FAX	19J1257FAX	20B0635FAX	20F0912FAX	20M1043FAX
		19A0492FAX	19E0512FAX	19K0105FAX	20B0867FAX	20F0913FAX	20M1044FAX
		19A0493FAX	19E0513FAX	19K0106FAX	20B0869FAX	20F0914FAX	20M1103FAX
		19A0494FAX	19E0586FAX	19K0107FAX	20B0870FAX	20F1235FAX	20M1104FAX
		19A0613FAX	19E0587FAX	19K0477FAX	20B1331FAX	20F1236FAX	20M1105FAX
		19A0614FAX	19E0588FAX	19K0478FAX	20B1332FAX	20G0160FAX	20M1456FAX
		19A0618FAX	19E0589FAX	19K0479FAX	20B1333FAX	20G0161FAX	20M1457FAX
		19A0619FAX	19E0793FAX	19K0480FAX	20B1334FAX	20G0162FAX	20M1458FAX
		19A0620FAX	19E0794FAX	19K0481FAX	20B1336FAX	20G0211FAX	20M1461FAX
		19A0672FAX	19E0795FAX	19K0564FAX	20B1337FAX	20G0212FAX	20M1462FAX
		19A0673FAX	19E0796FAX	19K0565FAX	20C0040FAX	20G0213FAX	21A0434FAX
		19A0674FAX	19E1113FAX	19K0566FAX	20C0041FAX	20G0214FAX	21A0435FAX
		19A1153FAX	19E1114FAX	19K0567FAX	20C0042FAX	20G0304FAX	21A0436FAX
		19A1154FAX	19E1115FAX	19K0568FAX	20C0227FAX	20G0305FAX	21A0437FAX
		19A1155FAX	19E1275FAX	19K0801FAX	20C0228FAX	20G0306FAX	21A0438FAX
		19A1156FAX	19E1307FAX	19K0802FAX	20C0229FAX	20G0307FAX	21A0613FAX
		19A1271FAX	19E1308FAX	19K0895FAX	20C0420FAX	20G0313FAX	21A0618FAX
		19A1272FAX	19E1309FAX	19K0896FAX	20C0421FAX	20G0314FAX	21A0888FAX
		19A1273FAX	19E1446FAX	19K1050FAX	20C0422FAX	20G0315FAX	21A0889FAX
		19A1274FAX	19E1447FAX	19K1051FAX	20C0423FAX	20G0316FAX	21A0890FAX
		19A1275FAX	19E1574FAX	19K1128FAX	20C0424FAX	20H1148FAX	21A1119FAX
		19B0043FAX	19E1575FAX	19K1129FAX	20C0612FAX	20H1151FAX	21A1120FAX
		19B0044FAX	19E1576FAX	19K1247FAX	20C0614FAX	20H1152FAX	21A1121FAX
		19B0283FAX	19E1577FAX	19K1248FAX	20C0615FAX	20H1153FAX	21A1122FAX

Item Code	Description	Affected Lot Numbers					
		19B0284FAX	19F0158FAX	19K1249FAX	20C0791FAX	20J0243FAX	21B0161FAX
		19B0285FAX	19F0159FAX	19K1273FAX	20C0792FAX	20J0244FAX	
		19B0286FAX	19F0160FAX	19K1274FAX	20C0793FAX	20J0612FAX	
		19B0287FAX	19F0161FAX	19K1275FAX	20C0794FAX	20J0614FAX	
353/19004Z	HME for Tracheostomized Patients JAPAN	18M0063FAX	19B1206FAX	19D1034FAX	19J1188FAX	19M0929FAX	20J0454FAX
		18M0064FAX	19C1282FAX	19D1149FAX	19J1189FAX	20A0102FAX	20J0455FAX
		18M0416FAX	19C1283FAX	19D1347FAX	19J1190FAX	20A0103FAX	20J0456FAX
		18M0417FAX	19C1284FAX	19E0232FAX	19K1126FAX	20A0104FAX	20J0457FAX
		18M0578FAX	19D0153FAX	19E0233FAX	19K1127FAX	20D1408FAX	20K1560FAX
		18M0579FAX	19D0154FAX	19E0234FAX	19L0931FAX	20F0799FAX	20M0076FAX
		19A0675FAX	19D0513FAX	19E0585FAX	19L0932FAX	20F0800FAX	20M0078FAX
		19A1151FAX	19D0515FAX	19F0954FAX	19L1339FAX	20F0801FAX	21B0014FAX
		19A1152FAX	19D0516FAX	19F0955FAX	19M0516FAX	20G0218FAX	21B0015FAX
		19B1205FAX					
		17E1490FAX	19A0902FAX	19H0919FAX	19M0996FAX	20C0489FAX	20D1319FAX
		17E1491FAX	19B1154FAX	19H0920FAX	19M0997FAX	20D0596FAX	20E0304FAX
		17E1492FAX	19D1246FAX	19J1115FAX	19M0998FAX	20D0597FAX	20F1079FAX
		17E1493FAX	19D1247FAX	19L0053FAX	20B0996FAX	20D0598FAX	20J0443FAX
		18M1017FAX	19D1317FAX	19L0054FAX	20B1042FAX	20D0700FAX	20L1339FAX
		18M1018FAX	19D1318FAX	19L0055FAX	20C0488FAX	20D1318FAX	20M1335FAX
		18M1019FAX	19F1047FAX	19M0995FAX			
353S19007Z	Adult-Pediatric HME JAPAN	18L0379FAX	19B1155FAX	19E1393FAX	19K0779FAX	20B1224FAX	20J0444FAX
		18L1250FAX	19C0762FAX	19G1492FAX	19K1041FAX	20G1671FAX	20J0445FAX
		19B0605FAX	19D0092FAX				
353S19046	HME for Tracheostomized Patients	19B1473FAX	19F0870FAX	20B0890FAX	20G1669FAX	20G1670FAX	20H0377FAX
		19F0869FAX	19G0049FAX	20G0428FAX			
		18L0884FAX	19B1157FAX	19D0962FAX	19K0392FAX	20B0597FAX	20G0431FAX
		18M0502FAX	19B1544FAX	19D0963FAX	19K1250FAX	20B0598FAX	20J0752FAX
		18M0504FAX	19C0233FAX	19E0789FAX	19L0810FAX	20B1196FAX	20K0718FAX
		18M0739FAX	19C0234FAX	19F0327FAX	19L1011FAX	20B1198FAX	20K0720FAX
		18M0740FAX	19C0235FAX	19F1064FAX	19L1012FAX	20C0095FAX	20K0721FAX
		18M0939FAX	19C0236FAX	19G0813FAX	19L1341FAX	20D0375FAX	20K1530FAX
		18M0951FAX	19C0237FAX	19G0814FAX	19L1342FAX	20D0376FAX	20K1531FAX
		19A0608FAX	19C1136FAX	19G1038FAX	19M0617FAX	20D0377FAX	20L0010FAX
		19A1088FAX	19C1137FAX	19G1343FAX	19M0618FAX	20E0373FAX	20L0521FAX
		19A1090FAX	19C1555FAX	19G1344FAX	19M0619FAX	20E0374FAX	20L0522FAX
		19A1328FAX	19C1556FAX	19J0014FAX	20A0325FAX	20F0536FAX	20L1340FAX
		19A1329FAX	19D0093FAX	19J0815FAX	20A1015FAX	20G0429FAX	20M0535FAX
		19B0934FAX	19D0094FAX	19J1319FAX	20A1016FAX	20G0430FAX	20M1333FAX
		19B0935FAX	19D0774FAX	19J1337FAX	20A1491FAX		
354/5414	Adult-Pediatric Mechanical Filter HME with Catheter Mount	20D0499FAX					
354/5833	Adult-Pediatric Mechanical Filter HME with Catheter Mount	20D0500FAX					
354/5841	Adult-Pediatric Mechanical Filter HME with Catheter Mount	20D0501FAX					
		18D0307FAX	19B0717FAX	19H0394FAX	19M0456FAX	20D0657FAX	20K0711FAX
		18L0331FAX	19B1100FAX	19H0395FAX	19M0993FAX	20D0774FAX	20K0712FAX
		18L0332FAX	19B1101FAX	19H0877FAX	19M0994FAX	20D0775FAX	20K0851FAX
		18L0333FAX	19B1102FAX	19H0878FAX	20A0681FAX	20E0710FAX	20K0852FAX
		18L0334FAX	19B1103FAX	19H0879FAX	20A0683FAX	20E0711FAX	20K0984FAX
		18M0089FAX	19B1104FAX	19H0926FAX	20A0684FAX	20E0712FAX	20K0985FAX
		18M0090FAX	19C0486FAX	19H0927FAX	20A0787FAX	20E0953FAX	20K0986FAX
		18M0091FAX	19C0487FAX	19H0928FAX	20A0788FAX	20E1029FAX	20L0578FAX
		18M0092FAX	19C0488FAX	19J0151FAX	20A0790FAX	20E1030FAX	20L0579FAX
		18M0093FAX	19C0489FAX	19J0152FAX	20A0973FAX	20E1341FAX	20L0580FAX
		18M0094FAX	19C0490FAX	19J1261FAX	20A1493FAX	20E1342FAX	20L0650FAX
		18M0509FAX	19C0684FAX	19J1262FAX	20A1494FAX	20G0334FAX	20L0651FAX
		18M0961FAX	19D0689FAX	19K0034FAX	20B0251FAX	20G0335FAX	20L0670FAX
		18M0962FAX	19D0690FAX	19K0035FAX	20B0252FAX	20G0336FAX	20L0671FAX
		18M0963FAX	19D1135FAX	19K0207FAX	20B0253FAX	20G0337FAX	20M0672FAX
		18M0964FAX	19F0238FAX	19K0208FAX	20B0254FAX	20G0338FAX	20M0673FAX
		18M1015FAX	19F0239FAX	19K0666FAX	20B0255FAX	20G0339FAX	20M0674FAX
		18M1016FAX	19F0240FAX	19K0667FAX	20B0586FAX	20G0340FAX	20M0675FAX
		19A0314FAX	19F0241FAX	19K1180FAX	20B1313FAX	20G0341FAX	20M0676FAX
		19A0315FAX	19F0242FAX	19K1182FAX	20B1348FAX	20H0886FAX	20M0730FAX
		19A0475FAX	19F0243FAX	19K1339FAX	20B1349FAX	20H0887FAX	20M0731FAX
		19A0701FAX	19F0606FAX	19L1194FAX	20C0157FAX	20H0888FAX	20M0732FAX
		19A0702FAX	19F0607FAX	19L1195FAX	20C0440FAX	20H0889FAX	20M0733FAX
		19A0703FAX	19F0837FAX	19L1196FAX	20C0441FAX	20H1160FAX	20M1480FAX
		19A0704FAX	19G0692FAX	19L1197FAX	20D0113FAX	20H1161FAX	20M1481FAX
		19A0705FAX	19G0693FAX	19L1343FAX	20D0114FAX	20H1162FAX	20M1482FAX
		19A0706FAX	19G0694FAX	19L1344FAX	20D0115FAX	20H1163FAX	20M1483FAX
		19A1163FAX	19G0695FAX	19M0125FAX	20D0164FAX	20J0970FAX	
		19B0713FAX	19G0696FAX	19M0126FAX	20D0165FAX	20J1037FAX	
		19B0714FAX	19G0697FAX	19M0127FAX	20D0294FAX	20J1038FAX	
		19B0715FAX	19G0698FAX	19M0210FAX	20D0655FAX	20K0709FAX	

Item Code	Description	Affected Lot Numbers					
		19B0716FAX	19G0699FAX	19M0390FAX	20D0656FAX	20K0710FAX	
354/5876TC	Adult-Pediatric Mechanical Filter HME	18L0490FAX	19D0693FAX	20D0967FAX	20K0874FAX	20L0146FAX	21A0812FAX
		18M0520FAX	19D1300FAX	20G1580FAX	20K1091FAX	20M1045FAX	21A0813FAX
		19B0983FAX	19J1271FAX	20H1124FAX	20L0145FAX	20M1046FAX	21A1195FAX
		19C0585FAX					
354/5900	Adult-Pediatric Mechanical Filter HME with Connector	18L0336FAX	19F0247FAX	19K1269FAX	19M0957FAX	20B1155FAX	20D0502FAX
		20E0632FAX	20E0633FAX	20K0831FAX	20K0832FAX	20M0769FAX	
354/5900TC	Adult-Pediatric Mechanical Filter HME with Connector	20E0634FAX	20K0833FAX				
354/5902	Adult-Pediatric Mechanical Filter HME with Catheter Mount	20H0362FAX	20H0364FAX	20H0366FAX	20H1240FAX	20J0226FAX	20L0941FAX
		20H0363FAX	20H0365FAX	20H1239FAX			
354S19028	Adult-Pediatric Mechanical Filter HME	19B1172FAX	19C0528FAX	19C0590FAX	19C0591FAX	19C1432FAX	19G1184FAX
		19G1185FAX	19J0612FAX				
354S19028TC	Adult-Pediatric Mechanical Filter HME	19C0961FAX					
354S19031	Adult-Pediatric Mechanical Filter HME with Connector	19B1176FAX	19D0206FAX				
354S19031Z	Adult-Pediatric Mechanical Filter HME with Connector	19D0207FAX					
355/5427	Pediatric Neonatal Electrostatic Filter HME	18L0577FAX	19D0331FAX	19H0820FAX	20A1019FAX	20D1051FAX	20H1112FAX
		18M0055FAX	19E0974FAX	19H0821FAX	20A1093FAX	20D1240FAX	20H1113FAX
		19A0724FAX	19E0975FAX	19H0822FAX	20A1094FAX	20D1241FAX	20J0352FAX
		19A0725FAX	19E0976FAX	19H0824FAX	20A1095FAX	20D1242FAX	20J0353FAX
		19A0726FAX	19E0977FAX	19J0457FAX	20A1362FAX	20D1243FAX	20J0354FAX
		19A0727FAX	19E0978FAX	19K0042FAX	20A1363FAX	20E1100FAX	20J0355FAX
		19A0728FAX	19E0980FAX	19K0043FAX	20A1364FAX	20E1101FAX	20J0356FAX
		19A1411FAX	19E1436FAX	19K0044FAX	20C0200FAX	20E1102FAX	20J0357FAX
		19A1412FAX	19E1437FAX	19K0196FAX	20C0201FAX	20E1103FAX	20J0470FAX
		19B0016FAX	19E1566FAX	19K0246FAX	20C0454FAX	20E1318FAX	20J0471FAX
		19B0017FAX	19E1567FAX	19K0247FAX	20C0455FAX	20E1319FAX	20J0472FAX
		19B0018FAX	19E1568FAX	19K0248FAX	20D0504FAX	20F0157FAX	20K1126FAX
		19B0214FAX	19E1569FAX	19L0619FAX	20D0505FAX	20F0158FAX	20L0005FAX
		19B0215FAX	19F0041FAX	19L0620FAX	20D0506FAX	20F0159FAX	20L0006FAX
		19B0276FAX	19F0155FAX	19L0621FAX	20D0673FAX	20F0160FAX	20L0007FAX
		19B0277FAX	19F0156FAX	19L0622FAX	20D0674FAX	20F0161FAX	20M0799FAX
		19B0554FAX	19G0472FAX	19L0623FAX	20D0675FAX	20G0977FAX	20M0800FAX
		19B0967FAX	19G0473FAX	19L1003FAX	20D0812FAX	20G0978FAX	20M1469FAX
		19C1226FAX	19G0774FAX	19L1004FAX	20D0813FAX	20G0979FAX	20M1470FAX
		19C1227FAX	19G0775FAX	19L1005FAX	20D0814FAX	20G0980FAX	20M1471FAX
		19C1228FAX	19G0850FAX	19L1276FAX	20D0887FAX	20G0981FAX	20M1472FAX
		19C1229FAX	19G0851FAX	19L1277FAX	20D0888FAX	20G0982FAX	21A0328FAX
		19C1230FAX	19H0121FAX	20A0576FAX	20D0889FAX	20G0983FAX	21A0329FAX
		19C1231FAX	19H0122FAX	20A0577FAX	20D1046FAX	20G0984FAX	21A0392FAX
19D0087FAX	19H0123FAX	20A0578FAX	20D1047FAX	20H1108FAX	21A0393FAX		
19D0088FAX	19H0374FAX	20A0579FAX	20D1048FAX	20H1109FAX	21A0800FAX		
19D0089FAX	19H0375FAX	20A1017FAX	20D1049FAX	20H1110FAX			
19D0330FAX	19H0819FAX	20A1018FAX	20D1050FAX	20H1111FAX			
355/5427TC	Pediatric Neonatal Electrostatic Filter HME	19A0729FAX	19H1005FAX	19K0333FAX	19M0034FAX	20A0771FAX	20K1128FAX
		20L0964FAX	20M0839FAX				
355/5427Z	Pediatric-Neonatal Electrostatic Filter HME JAPAN	19C1232FAX	19G0852FAX	19L1006FAX	20E1099FAX	20J0611FAX	20K1127FAX
		19D0086FAX	19H1007FAX	19L1007FAX	20G0986FAX	20J0934FAX	20L0008FAX
		19D0328FAX	19K0040FAX	20C0202FAX	20G0987FAX	20J0935FAX	21A0799FAX
		19F0286FAX					
		18L0019FAX	19D0698FAX	19J0266FAX	20B0466FAX	20E0298FAX	20H0361FAX
		18L0020FAX	19D0699FAX	19J0267FAX	20B0467FAX	20E0299FAX	20J0821FAX
		18L0021FAX	19D0700FAX	19J0268FAX	20B1064FAX	20E0300FAX	20J0822FAX
		18M0511FAX	19D1150FAX	19J0269FAX	20B1065FAX	20E0301FAX	20J0823FAX
		18M0512FAX	19D1151FAX	19J0270FAX	20B1066FAX	20E0364FAX	20J0824FAX
		18M0513FAX	19D1152FAX	19J0669FAX	20B1067FAX	20E0366FAX	20J0825FAX
		18M0514FAX	19D1189FAX	19J0737FAX	20B1068FAX	20E0367FAX	20J1080FAX
		18M0539FAX	19D1279FAX	19J0738FAX	20B1069FAX	20E0368FAX	20J1081FAX
		18M0979FAX	19E0472FAX	19J0739FAX	20B1321FAX	20E0566FAX	20J1082FAX
		18M0980FAX	19E0473FAX	19J0740FAX	20B1323FAX	20E0664FAX	20J1083FAX
		18M0981FAX	19E0474FAX	19J0979FAX	20B1324FAX	20E0665FAX	20J1084FAX
		18M0982FAX	19E0570FAX	19J0980FAX	20B1325FAX	20E0666FAX	20J1085FAX
		19A0479FAX	19E0571FAX	19K0431FAX	20C0043FAX	20E0667FAX	20J1383FAX
		19A0480FAX	19E0572FAX	19K0432FAX	20C0044FAX	20E0668FAX	20J1384FAX
		19A0746FAX	19E1033FAX	19K0625FAX	20C0045FAX	20E1007FAX	20J1385FAX
		19A0747FAX	19E1034FAX	19K0626FAX	20C0046FAX	20F0513FAX	20J1386FAX
		19A0748FAX	19F0248FAX	19K0653FAX	20C0176FAX	20F0514FAX	20J1387FAX
		19A0749FAX	19F0249FAX	19K0803FAX	20C0177FAX	20F0515FAX	20K0097FAX
		19A0750FAX	19F0250FAX	19K0804FAX	20C0203FAX	20F0516FAX	20K0098FAX
		19B0505FAX	19F0251FAX	19K0894FAX	20C0230FAX	20F0673FAX	20K0099FAX
		19B0506FAX	19F0428FAX	19K1014FAX	20C0450FAX	20F0815FAX	20K0100FAX
		19B0507FAX	19F0429FAX	19K1015FAX	20C0451FAX	20F0816FAX	20K0101FAX
		19B0508FAX	19F0431FAX	19K1036FAX	20C0702FAX	20F0817FAX	20K0218FAX
		19B0718FAX	19F0774FAX	19K1037FAX	20C0703FAX	20F0818FAX	20K0219FAX

Item Code	Description	Affected Lot Numbers							
355/5430	Infant-Pediatric Electrostatic Filter HME	19B0719FAX	19F0775FAX	19K1038FAX	20C0704FAX	20F0819FAX	20K0220FAX		
		19B0720FAX	19F0776FAX	19K1254FAX	20C0705FAX	20F1142FAX	20K0532FAX		
		19B0797FAX	19F0970FAX	19K1255FAX	20C0706FAX	20F1143FAX	20K0533FAX		
		19B0798FAX	19F0971FAX	19K1359FAX	20C0868FAX	20F1144FAX	20K0534FAX		
		19B0799FAX	19F0972FAX	19K1360FAX	20C0869FAX	20F1145FAX	20K0535FAX		
		19B1072FAX	19F1197FAX	19M0112FAX	20C0870FAX	20F1146FAX	20L0339FAX		
		19B1073FAX	19F1198FAX	19M0113FAX	20C0946FAX	20F1147FAX	20L0340FAX		
		19B1075FAX	19F1199FAX	19M0114FAX	20C0947FAX	20F1148FAX	20L0396FAX		
		19B1077FAX	19G0149FAX	19M0392FAX	20C0948FAX	20G0010FAX	20L0397FAX		
		19B1078FAX	19G0150FAX	19M0393FAX	20C1142FAX	20G0011FAX	20L0430FAX		
		19B1524FAX	19G0151FAX	19M0442FAX	20C1143FAX	20G0012FAX	20L0431FAX		
		19B1525FAX	19G0152FAX	19M0909FAX	20C1144FAX	20G1372FAX	20L0432FAX		
		19B1527FAX	19G0153FAX	19M1003FAX	20C1237FAX	20G1373FAX	20M1072FAX		
		19B1528FAX	19G1015FAX	19M1004FAX	20C1239FAX	20G1374FAX	21A0478FAX		
		19B1529FAX	19G1016FAX	20A0284FAX	20C1240FAX	20G1375FAX	21A0479FAX		
		19B1530FAX	19G1018FAX	20A0285FAX	20C1241FAX	20G1376FAX	21A0480FAX		
		19C0134FAX	19G1019FAX	20A0377FAX	20D0123FAX	20G1377FAX	21A0481FAX		
		19C0135FAX	19G1350FAX	20A0378FAX	20D0124FAX	20H0048FAX	21A1114FAX		
		19C0136FAX	19G1351FAX	20A1500FAX	20D0125FAX	20H0049FAX	21A1115FAX		
		19C0376FAX	19G1352FAX	20B0244FAX	20D1621FAX	20H0050FAX	21A1116FAX		
		19C0377FAX	19J0056FAX	20B0245FAX	20D1622FAX	20H0051FAX			
		19C0378FAX	19J0057FAX	20B0246FAX	20D1623FAX	20H0052FAX			
		19C0485FAX	19J0058FAX	20B0364FAX	20D1624FAX	20H0358FAX			
		19D0630FAX	19J0059FAX	20B0365FAX	20D1625FAX	20H0359FAX			
		19D0631FAX	19J0060FAX	20B0366FAX	20E0297FAX	20H0360FAX			
		355/5430TC	Infant-Pediatric Electrostatic Filter HME	18L0357FAX	19A0507FAX	19A0508FAX	19F0753FAX	19J0053FAX	19J1118FAX
				19L0536FAX	20A0702FAX	20B0457FAX	20E1028FAX	20G0013FAX	
		355/5430TCZ	Infant-Pediatric Electrostatic Filter HME	18L0356FAX	19A0509FAX	19C0470FAX	19C0589FAX	19E0406FAX	19F0754FAX
				19L0537FAX	20G0014FAX	20G1584FAX	20H0037FAX	20J1125FAX	
		355/5430Z	Infant-Pediatric Electrostatic Filter HME JAPAN	18M0540FAX	19D0804FAX	19J0055FAX	20B0469FAX	20G0051FAX	20J0547FAX
19C0144FAX	19D0805FAX			19K0430FAX	20D1620FAX	20G0052FAX	20J0548FAX		
19D0632FAX	19D1153FAX			19K1344FAX	20E0320FAX	20H0357FAX	20J0757FAX		
19D0634FAX	19D1154FAX			19K1345FAX	20F0654FAX	20J0545FAX	20J0758FAX		
19D0701FAX	19F0253FAX			20B0061FAX	20F0655FAX	20J0546FAX	20K0347FAX		
19D0703FAX	19J0054FAX	20B0468FAX							
355/5814	Infant-Pediatric Electrostatic Filter HME with Catheter Mount	20H0618FAX	20H0619FAX	20H0620FAX					
355/5860	Pediatric-Neonatal Electrostatic Filter HME with Connector	19A1409FAX	19D0090FAX	19G0773FAX	19K0039FAX	19L1008FAX	20A1079FAX		
		20C0256FAX	20D0304FAX						
355/5884	Infant-Pediatric Electrostatic Filter HME with Connector	18M0541FAX	18M0543FAX	19B0511FAX	19C0455FAX	19G0157FAX	19J0050FAX		
		19M0924FAX	20A1502FAX	20D0063FAX	20D0251FAX	20D1649FAX	20L0338FAX		
		21A0879FAX							
355/5884Z	Infant-Pediatric Electrostatic Filter HME with Connector	18M0542FAX	19D1268FAX	19L0119FAX	20D1650FAX	20J0543FAX	20K0339FAX		
		19C0456FAX	19F0252FAX	19M0099FAX	20F0511FAX	20J0544FAX	20K0340FAX		
		19D1266FAX	19G0158FAX	19M0115FAX	20F0512FAX	20J0764FAX	20K0341FAX		
		19D1267FAX	19J0051FAX	20A1503FAX	20J0542FAX	20J0765FAX	20K0342FAX		
355/5916	Pediatric-Neonatal Electrostatic Filter HME with Connector	19A0723FAX	19D0091FAX	20A1085FAX	20C0452FAX	20D0305FAX	20F0313FAX		
		19A1410FAX	19L1009FAX						
444S02608	Pediatric-Neonatal Closed Suction System	20J0829FAX							
444SP01014	Closed Suction System 14CH FIXED	19C0257FAX							
444SP01016	Closed Suction System 16CH FIXED	20H1281FAX	20K1057FAX	20L0028FAX					
444SP01116	Closed Suction System 16CH FIXED CON MDI	20H0463FAX	20H0494FAX	20H1282FAX					
444SP01312	Closed Suction System 12CH FIXED TRACHEOSTOM	20G1527FAX							
444SP01314	Closed Suction System 14CH FIXED TRACHEOSTOM	20H0629FAX							
444SP01512	Closed Suction System 12CH FIXED CON MDI	20G1528FAX							
444SP02010	Closed Suction System, Replcble Catheter 10CH	20G1540FAX							
444SP02012	Closed Suction System, Replcble Catheter 12CH	18E1563FAX	19B0996FAX	20G1164FAX	20G1529FAX	20M0351FAX			
444SP02014	Closed Suction System, Replcble Catheter 14CH	20J0835FAX	20J0836FAX	20L0151FAX					
444SP02016	Closed Suction System, Replcble Catheter 16CH	20H0495FAX	20K1058FAX						
444SP02110	Closed Suction System, Replcble Cath. MDI 10CH	20G1541FAX							
444SP02112	Closed Suction System, Replcble Cath. MDI 12CH	20G1530FAX							
444SP02114	Closed Suction System, Replcble Cath. MDI 14CH	20J0837FAX	20L0155FAX						
444SP02116	Closed Suction System, Replcble Cath. MDI 16CH	19B1299FAX	20H0496FAX	20K1059FAX					
444SP02312	Closed Suction System, Replcble Cath. Trach. 12CH	20G1531FAX							
444SP02314	Closed Suction System, Replcble Cath. Trach. 14CH	20H0630FAX	20H0631FAX						
444SP02514	Closed Suction System, Replcble Cath. MDI 14CH	20J0842FAX							
444SP03014	Closed Suction System, Replcble Cath. 14CH X CPAP	20K0232FAX							
444SPY1314	Closed Suction System 14CH FIXED X SHILEY	20G0838FAX							
608/5092	Connector 22M-22M	20L0505FAX							
610/5257	Connector 22M/15F-15M	20H0557FAX							
616/5256	Adult Connector	20H0212FAX							
620/5173	Connector 'L'22M/15F-15M	19L1306FAX	20M0253FAX						