



No. DMD0001 : الرقم

Date : 27/07/2021 : التاريخ

FSN Letter

Damad-set Dressing pack; Oral Care Kit

Introduction:

In the interest of the National Company for Medical Products Ltd. (Damad) for the safety and security of its products and in compliance with the directives of the relevant regulatory authorities (the Food and Drug Authority), we inform you that we have received a notification that the Food and Drug Authority has recalled the product Mcare Cool mint Mouthwash from the Gulf Care Factory for health and cosmetic care products, as it contains: A high percentage of bacteria and fungi that exceeded the maximum allowed in addition to containing bacteria (Burkholdria cepacia complex).

Description:

Due to the presence of the aforementioned product as a component within our product oral care set that have been supplied to a number of healthcare facilities.

Thus as part of mutual collaboration & coordination with the NCMDR - SFDA, we decided to issue this FSN Letter for the below product & details:

- Product Name: Damad-set Dressing pack; Oral Care Kit).
- Manufacturer: The National Medical Products Co. Ltd.
- Reference number: 228022/1.
- Lot Number: all
- Distribution:
 - NUPCO Warehouses / Medical Supply - MOH
 - Public healthcare facilities

Potential Risk:

- It may cause respiratory infection.

بِسْمِ اللَّهِ الرَّحْمَنِ الرَّحِيمِ

THE NATIONAL MEDICAL PRODUCTS CO. LTD.

Limited Liability Company (**Damad**)
The Capital 30,000,000 S.R. Paid Full
Riyadh - Saudi Arabia
I. L. No. 78 Dated 1/4/1396H
C.C.R. No. 9468 Riyadh



الشركة الوطنية للمنتجات الطبية المحدودة
شركة ذات مسئولية محدودة (ضمان)
رأس المال ٣٠.٠٠٠.٠٠٠ ريال سعودي مدفوع بالكامل
الرياض - المملكة العربية السعودية
التصريح رقم ٧٨ بتاريخ ١٣٩٦/٤/١هـ
رقم الضريبة بالغرفة ٩٤٦٨ الرياض

No. DMD0001 : الرقم

Date : 27/07/2021 : التاريخ

Instructions:

We kindly, request all NUPCOMedical Supply Warehouses & Public healthcare facilities to:

- Stop using the Damad-set Oral Care Set.
- Coordinate to withdraw the remaining quantities of the product.
- Communicate in this regard through the contact details below.

Contact:

Contact information:

Mobile (0592219500) or e-mail (monem.alattar@damad.com).

Mustafa Ahmed
NCMDR Officer