

EDA Home<sup>3</sup> Medical Devices<sup>4</sup> Databases<sup>5</sup>
Class 2 Device Recall Stryker Sustainability Solutions Reprocessed BW Lasso 2515 NAV eco Variable Diagnostic EP Catheters,

6 510(k) |DeNovo<sup>8</sup>| Registration & | Adverse | |Recalls<sup>11</sup>|PMA<sup>12</sup>|HDE<sup>13</sup>|Classification<sup>14</sup>|Standards<sup>15</sup>

Listing<sup>9</sup> Events<sup>10</sup> SwoerSearch

CFR Title 21<sup>16</sup>|Radiation-Emitting Products<sup>17</sup>|X-Ray Assembler<sup>18</sup>|Medsun Reports<sup>19</sup>|CLIA<sup>20</sup>|TPLC<sup>21</sup>

New Search

Back to Search Results See Related

Class 2 Device Recall Stryker Sustainability Solutions Reprocessed BW Lasso 2515 NAV eco Variable Diagnostic EP

Catheters,

Date Initiated by Firm

January 02, 2018

**Create Date** 

March 06, 2018

Recall Status<sup>1</sup>

Open<sup>3</sup>, Classified

Recall Number

Z-0917-2018

Recall Event ID

7908723

510(K)Number

K112292<sup>24</sup>

**Product Classification** 

Catheter, recording, electrode, reprocessed<sup>25</sup> - Product Code NLH<sup>26</sup>

Product

Stryker Sustainability Solutions Reprocessed BW Lasso 2515 NAV eco Variable

Diagnostic EP Catheters, Ref D134302, Sterile, Rx.

The Reprocessed 2515 NAV eco Variable Electrophysiology (EP) Catheters are indicated for multiple electrode electrophysiological mapping of the cardiac structures of the heart, i.e. recording or stimulation only. They are designed to obtain electrograms in the atrial regions of the heart. The Reprocessed 2515 NAV eco Variable EP Catheters provide location information when used with compatible

CARTO EP Navigation Systems version 2.3 or higher.

**Code Information** 

Serial numbers 2468846, 2468849, 2509029, 2489074, 2491174, 2500529, 2465989, and

2500528

Recalling Firm/ Manufacturer

Stryker Sustainability Solutions

1810 W Drake Dr Tempe AZ 85283-4327

For Additional Information Contact Kara Madsen 888-888-3433

Manufacturer Reason for Recall

An EEPROM chip error code may occur when the catheters are used with CARTO EP

**Navigation Systems** 

**FDA Determined** Cause 2

Device Design

Action

Stryker sent an Urgent Medical Device Recall letter dated January 2, 2018, explains the reason for recall, risk to health, and instructs them to discontinue use of the affected products. The customer is to check their inventory for the affected serial numbers provided in the attachment and complete the recall effectiveness check form to indicate the amount

of product in inventory. Return instructions are provided. For further questions, please call (888) 763-8803.

Quantity in Commerce

8 devices

Distribution

USA (nationwide) Distribution was made to medical facilities in AZ, CA, CT, FL, GA, IL, KS,