



REF:131/12.21/GH.SA/FA

رقم المحفوظات:
بيروت في:

٢٤

٢١ كانون الأول ٢٠٢١

١٥٠٠
قرار رقم
يتعلق بتفاصيل توضيحية حول القرار رقم ١٤٠/٢١ تاريخ ١١/١١/٢٠٢١

لوائح باللوازم والمواد الطبية المدعومة لزوم القلب وغسيل الكلى

إن وزير الصحة العامة،
بناءً على المرسوم رقم ٨٣٧٦ تاريخ ٢٠٢١/٠٩/١٠ (تشكيل الحكومة)،
بناءً على المرسوم رقم ٨٣٧٧ تاريخ ١٩٩٢/١٢/٣١ وتعديلاته (تنظيم وزارة الصحة العامة)،
بناءً على القرار التنظيمي رقم ١/٤٥٥ تاريخ ٢٦ نيسان ٢٠١٣ (تنظيم أذخار واستعمال المواد الطبية المختلفة)
عطفًا على القرار رقم ١/١٣٤٠ تاريخ ١١/١١/٢٠٢١ المتعلق بتحديد اصناف المستلزمات المدعومة ونسب الدعم
المعتمدة،

بناءً على ضرورات المصلحة العامة،
بناءً على اقتراح مدير عام وزارة الصحة العامة،

يقرر ما يأتي:

المادة الأولى: تُحدد اللوازم والمواد الطبية المدعومة لزوم القلب وغسيل الكلى وفقاً للوائح المرفقة.

المادة الثانية: يبلغ هذا القرار حيث تدعو الحاجة ٥٠%.

وزير الصحة العامة
د. فراس الأبيض

نسخة إلى:

- مصرف لبنان
- نقابة مستوردي المستلزمات الطبية
- نقابة المختبرات الخاصة
- نقابة مستوردي الادوية
- نقابة الصيدلة
- المديرية العامة لوزارة الصحة
- مصلحة الصيدلة
- مصلحة الهندسة الصحية
- دائرة البرامج والمشاريع
- الموقع الالكتروني لوزارة الصحة
- المحفوظات.

Essential list MD for dialysis	
Description	GMDN Code
Dialysis fluid Filter	47072
Arterial & venous Fistula needles set	12741, 47572, 60879
Bloodlines, tubing systems	32142, 34999
Accessories, Catheter Adapters & Closure caps for peritoneal dialysis	35944, 35986
Acid Concentrate	35849
Bicarbonate cartridge	35849
Peritoneal Dialysis Solution	
Cholesterol Apheresis kit	

SA

GMDN_Code	GMDN_Term_Name
35039	Epicardial pacing lead
35223	Endocardial pacing lead
35852	Single-chamber implantable defibrillator
35853	Endocardial defibrillation lead
37265	Dual-chamber implantable defibrillator
46221	Pacing/defibrillation lead cap
47263	Cardiac resynchronization therapy implantable pacemaker
47265	Dual-chamber implantable pacemaker, rate-responsive
47267	Single-chamber implantable pacemaker, rate-responsive
47270	Cardiac resynchronization therapy implantable defibrillator
60910	Coronary venous pacing lead
45418	Cardiac occluder
45577	Mitral annuloplasty ring
45578	Mitral/tricuspid annuloplasty ring
53616	Bare-metal coronary artery stent
56284	Drug-eluting coronary artery stent, non-bioabsorbable-polymer-coated
58771	Drug-eluting coronary artery stent, bioabsorbable-polymer-coated
56280	Mitral valve tissue repair system
60240	Aortic bi-leaflet mechanical heart valve prosthesis
60241	Mitral bi-leaflet mechanical heart valve prosthesis

